

ELK-004 & 005


HYDRAULIC BATHING

TRANSFER





- A transfer stretcher with double hydraulic system. Functions operated by treading the foot pedals
- Steel framework, Bathing with environmental PVC cover and PE Foam
- **Dimension:**
Overall: 2045*770mm(L*W)
Height adjustment (Bed board to ground): 970-745mm
Trendelenburg & Rev. Trendelenburg angle: 0-13°
- Four pieces detachable protective railing
- Removing one of the long side rail, the bathing cloth able to be extended flat, easy for transfer patient to bed



ELK-005 Hydraulic Bathing Transfer



- A transfer stretcher with double hydraulic system
Functions operated by treading the foot pedals
- Bright arrow symbols direct the user how to operate clearly
- Central-lock castors, diameter 150mm. Matched with direction castor under the chassis help user easy control the stretcher direction
- **Dimension:**
Overall dimension: 1900*740mm(L*W)
Height adjustment(Bed board to ground): 830-595mm
Size of bed board: 1900*610mm(L*W)
Back-rest lifting angle: 0-55° ±15°
Trendelenburg & Rev. Trendelenburg angle: 0-13°
- Drop-down protective railing, PP bed board produced by blowing molding
- Net Weight: 82kg
- Standard accessories: Mattress, I.V. pole

