

# ***ELK-124 & 125***



## ELK-124 3 Cranks Manual Bed



- Steel framework, top class ABS & PP plastic
- Overall size: 2160mm(L) 960mm(W)\*575mm-750mm  
Bed board dimension: 1910mm(L)\*830mm(W)
- Three functions controlled by cranks  
Back-rest: 0-65°±5°  
Leg-rest: 0-45° ±5°  
Hi-Low: 750mm-575mm
- PP head & foot board and bed board, detachable and easy to clean
- 6 columns Collapsible Aluminum alloy protective railing
- 4pcs universal castor, Φ125mm, individual brake
- Net Weight: 102kg
- Standard accessory: I.V.Pole





## ELK-125 2 Cranks Manual Bed



- Steel framework, top class ABS & PP plastic
- Overall size: 2140mm(L)\*1020mm(W)\*510mm  
Bed board dimension: 1940mm(L) 830mm(W)
- Two functions controlled by cranks  
Back-rest lifting angle: 0-65° ±5°  
Knee-rest lifting angle: 0-30° ±5°  
Height(Bed board to ground): 510mm
- Detachable PP head & foot board
- 4pcs drop-down PP protective railing, controlled by gas spring
- Central brake castors, diameter 125mm. Brake pedal locate below the foot board
- Standard accessory: I.V.Pole

