


# ***ELK-132***

## ***MANUAL BED***



## ELK-132 Manual Bed

LS-MA 185 2 Crank Manual Bed  
( $\phi$  125mm castor individual brake)



LS-MA 185B 2 Crank Manual Bed  
( $\phi$  150mm castor, central brake)



LS-MA 185C 3 Crank Manual Bed  
( $\phi$  125mm castor, individual brake)



LS-MA 185D 3 Crank Manual Bed  
( $\phi$  150mm castor, central brake)



# ELK-132 Manual Bed

Item	Manual bed			
Code	LS-MA 185	LS-MA 185B	LS-MA 185C	LS-MA 185D
Info	2/3 functions manual bed, functions controlled by cranks			
Main material	Steel framework, top class ABS & PP plastic			
Bed board	Powder coated cold-rolled punching panel, with ventilation holes			
Head & foot board	Top class PP & ABS plastic, with comer bumper			
Protective railing	Collapsible 6 Columns Aluminum alloy protective railing			
Overall dimension	With head & foot board; 2150mm(L) *985mm (W) Without head & foot board: 2000mm(L) *900mm (W)			
Bed board dimension	1875mm(L)*770mm(W)			
Height (bed board to ground)	500mm	480mm	490mm-700mm	500mm-715mm
Castor	φ125mm Individual brake	φ150mm Central brake	φ125mm Individual brake	φ150mm Central brake
Crank	2 cranks	2 cranks	3 cranks	3 cranks
Net Weight	83kg	93Kg	106Kg	115Kg
Functions	Back-rest: 0-70°±5° Leg-rest: 0-40° ±5° Sitting part sinking: 0-10°	Back-rest: 0-70° ±5° Leg-rest: 0-40° ±5° Sitting part sinking: 0-10°	ack-rest: 0-70° ±5° Leg-rest: 0-40° ±5° Sitting part sinking:0-10° Hi-low: 490mm-700mm	Back-rest: 0-70° ±5° Leg-rest: 0-40° ±5° Sitting part sinking: 0-10° Hi-low: 500mm-715mm
Load capacity	250kg			