

# ***ELK-342*** ***ELECTRIC BED***



## ELK-342 Electric Bed

LS-EA185B 2 function electric bed  
( $\Phi$ 150mm castor, central brake)



LS-EA185C 3 function electric bed  
( $\Phi$ 125mm castor, central brake)



LS-EA185D 3 function electric bed  
( $\Phi$ 150mm castor, central brake)



## ELK-342 Electric Bed

|                                 |   |   |   |
|---------------------------------|---|---|---|
| Item                            | Manual bed  |   |   |
| Code                            | LS-EA185B   | LS-EA185C   | LS-EA185D   |
| Info                            | 2/3 functions electric bed .Functions controlled by remote control                              |   |   |
| Main material                   | Steel framework, top class ABS & PP plastic   |   |   |
| Bed board                       | Powder coated cold-rolled punching panel, with ventilation holes                                |   |   |
| Head & foot board               | Top class PP & ABS plastic, with corner bumper  |   |   |
| Protective railing              | Collapsible 6 Columns Aluminum alloy protective railing   |   |   |
| Overall dimension               | With head & foot board: 2150mm(L) *985mm (W)<br>Without head & foot board: 2000mm(L) *900mm (W) |   |   |
| Bed board dimension             | 1875mm(L)*770mm(W)  |   |   |
| Height<br>(bed board to ground) | 480mm   | 485mm-695mm   | 500mm-710mm   |
| Castor                          | φ150mm<br>Central brake 150mm   | φ150mm<br>Central brake 150mm   | φ150mm<br>Central brake 150mm   |
| Net Weight                      | 105kg   | 114Kg   | 117Kg   |
| Functions                       | Back-rest: 0-70° ±5°<br>Leg-rest: 0-40° ±5°<br>Sitting part sinking: 0-10°                      | Back-rest: 0-70° ±5°<br>Leg-rest: 0-40° ±5°<br>Sitting part sinking: 0-10°<br>Hi-low: 485mm-695mm | Back-rest: 0-70° ±5°<br>Leg-rest: 0-40° ±5°<br>Sitting part sinking: 0-10°<br>Hi-low: 500mm-710mm |
| Load capacity                   | 250kg   |   |   |