

ELK-394G

SLIT LAMP MICROSCOPE



ELK-394G Slit Lamp Microscope

SPECIFICATIONS

- All the optical lenses are moisture proof, mildew proof and anti-reflected treated
- Large spot diameter, slit length reach to 14mm, more comprehensive view of the eye condition
- Multi-coated optical glass technology, imaging is clear and bright, the view field is wide
- Osram light bulbs, small heat, high-brightness, spot uniformity
- Handrail design, user-friendly and the convenience of patients
- All the joints can be connected by accessories with international standards

Microscope	Microscope Type	Galilean stereoscope microscope
	Magnification Selection	3 steps by drum rotation
	Eyepiece	12.5X
	Magnification Ratio (Field of View)	10X(22.5mm), 16X(14mm), 25X(8.8mm)
	Pupillary Adjustment	50mm-82mm
	Diopter Adjustment	+7D ~ -7D
Base	Side Shift	110mm
	Depth Shift	90mm
	Height Shift	30mm
Chin rest	Height Shift	80mm
	Fixation Target	Red LED
Illumination system	Slit Width	0-14mm Continues (at 14mm, slit becomes a circle)
	Slit Height	0-14mm Continues
	Slit Angle	0°~180° adjustable
	Slit Inclination	5°, 10°, 15°, 20°
	Diameter of Light Spot	Φ14mm, Φ10mm, Φ45mm, Φ3mm, Φ2mm, Φ1mm, Φ0.2mm and 1~14mm continuous
	Filter	Heat Absorption, Grey, Red-Free, Cobalt Blue
	Illumination Bulb	12V 50W halogen lamp
Power	Input Voltage	110V/220V AC
	Frequency	50/60Hz
	Power Consumption	68VA
Weight and Size	Gross Weight	18Kg
	Dimension	690mm(L)X440mm(W)X429mm(H)