

ELK-546 & 547

3 FUNCTIONS ELECTRIC BED



ELK-546 3 Functions Electric Bed



- Steel framework, top class ABS & PP plastic
- Overall dimension: 2160mm(L)*1120mm(W)*440mm-655mm
- Bed board dimension: 1935mm(L)*900mm(W)
- 3 functions controlled by remote control
Back-rest : 0-70° ±5°
Foot-rest : 0-30° ±5°
Hi-Low: 655mm-440mm
- Net Weight: 139KGS
- Load capacity: <250KGS
- Castor with central brake system, wheel diameter: 150mm
- Standard accessory: I.V. Pole



ELK-547

3 Functions Electric Bed



- Steel framework, top class ABS & PP plastic Overall
- Dimension: 2160mm(L) 970mm(W)*440mm-655mm
- Bed board dimension: 1935mm(L)*900mm(W)

3 functions controlled by remote control

Back-rest :0-70° ±5°

Foot-rest :0-30° ±5°

Hi-Low:655mm-440mm

- Net Weight: 112KGS
- Load capacity: <250KGS
- Castor with central brake system, wheel diameter: 150mm(Castor(125mm) with individual brake for option)
- Standard accessory: I.V. Pole

