

ELK-856

5-FUNCTIONS ELECTRIC HOSPITAL BED



ELK-856

5-functions electric hospital bed

Function

Back-rest lifting angle: 0° - 80°

Knee-rest lifting angle: 0° - 40°

Height lifting range: 500-750mm

Trendelenburg tilting angle: 0° - 15°

Reverse trendelenburg tilting angle: 0° - 15°

Size

L2120×W1100×H500-750mm

Feature

Cold steel plate wholly molded surface, 4 sections metal bed surface, connected with 2 soft links

Maximum load capacity: 250kg

Bed functions control by control remote

Standard Accessories

2 pcs PP head and foot board, detachable

4 pcs engineering plastic PP side rails

4 pcs central control brake castors, diameter 125mm

1 set advanced motor(include 4 linear actuators, 1 control box, 1 remote control),

1 I.V pole

4 I.V pole holders, 4 Drainage bag hooks