



**Elektro**  
GENESIS



**Dental  
Instrument  
Catalogue**



[www.electrogenesis.com](http://www.electrogenesis.com)







### Prosthodontic

89-99



### Accessories

128-129



### Orthodontic

100-115



### Index

130-145



### Instrument Management

116-127







# Diagnostics

Handle Options

Mirror Handles

Explorers

Micro Surgical Mirrors

Tweezers



HOW TO USE THIS SECTION

Instrument Name & Pattern: 23 Shepherd's Hook  
Part Code of pictured instrument: EXS6XL  
Available handle designs: Handle options : #4, #3, #2, #1

HANDLE OPTION

Most popular instruments available in all handle options.

DOUBLE END



SINGLE END



**NOTE :**  
This is not a mouth mirror handle.  
It's just only for the selection of instrument handle which you want to use.




Diagnostics

MIRROR HANDLE






MIRROR HANDLE SIGNATURE SERIES # GRIPS



**DID YOU KNOW ?**  
Colour coded instruments rings are silicon rubber rings for identifying and organising instruments. Useful for maintaining proper instruments order in cassettes and identifying instruments by procedure.

### PROBES

A double-ended instrument that combines two critical instruments required for every effective diagnostic set-up-with an explorer on one end and a probe on the other, an Probe is excellent for dental examinations and charting.










Diagnostics

MICRO SURGICAL MIRRORS


Double-sided-mirrors are manufactured from highly polished stainless steel and come in round-end or rectangular shapes.




	<p><b>Plain</b></p> <p>MMTP3 MMTP4 MMTP5</p> <p>— Set of 12 pcs</p>
	<p><b>Magnifying</b></p> <p>MMTM4 MMTM5</p> <p>— Set of 12 pcs</p>
	<p><b>Front surface</b></p> <p>MMTFS4 MMTFS5</p> <p>— Set of 12 pcs</p>
	<p><b>Front surface on both sides</b></p> <p>MMTDS5 — Set of 1 pc. MMTDC/6 — Set of 6 pcs</p>
	<p><b>Front surface on both sides</b> <b>CONE SOCKET MIRROR</b></p> <p>MMT5DSCS — Set of 1 pc. MMT5DSCS6 — Set of 6 pcs</p>

	<p>2.5mm Small Round</p> <p>MMRS</p>	<p>4.5mm Large Round</p> <p>MMR</p>	<p>4 x 12mm Large Rectangular</p> <p>MM3</p>
	<p>2 x 7mm Small Rectangular</p> <p>MM1</p>	<p>3 x 9mm Medium Rectangular</p> <p>MM2</p>	
	<p><b>Front surface on both sides with Handle</b> <b>CONE SOCKET</b> MMH5DSCS</p>		

**DID YOU KNOW ?**

 **Back surface**  
These mirrors have their reflective coatings on their rear surfaces.

 **Front surface**  
These mirrors have their reflective coating on their front surface, yielding undistorted imaging with no disturbing double images.


TWEEZER

Used to grasp or transfer materials in and out of the oral cavity. All dressing pliers are made of heavy-gauge, heat-treated stainless steel to prevent tip bending and misalignment.

London College  
DPU17



Mariam  
DP318



London College Premium  
DP2



Mariam Premium  
DP3



Locking Premium  
DP18L  
Locking Serrated  
DP18LS



GP Holding Premium  
EPL1  
GP Holding  
EPL1S



Brush Holding  
BHT5



Perry  
DP1



Curved Tweezer  
DP90





# Oral Surgery

Extraction Forceps

Root Elevators

Scissors & Needle Holders

Retractors

Periosteal Elevators







**S**  
Standard

**P**  
Premium


**E**  
Ergonomic

### EXTRACTION FORCEPS

Most popular instruments available in all handle options.




**1# Upper Anteriors**  
FX1S  
FX1P  
FX1E




**7# Upper Premolars**  
FX7S  
FX7P  
FX7E




**17# Upper Molars Right**  
FX17S  
FX17P  
FX17E




**18# Upper Molars Left**  
FX18S  
FX18P  
FX18E




**18A# Upper Molars**  
FX18AS




**29# Upper Roots**  
FX29S  
FX29P




**29N# Upper Roots Narrow**  
FX29NS  
FX29NP




**30# Upper Roots**  
FX30S  
FX30P



**55# Separating Upper Molars**  
FX55S



**56# Separating Lower Molars**  
FX56S




EXTRACTION FORCEPS



41# Upper Roots  
FX41S  
FX41P



44# Upper Roots  
FX44S  
FX44P



51# Upper Roots  
FX51S



51A# Upper Roots  
FX51AS  
FX51AP  
FX51AE



67A# Upper Third Molars  
FX67AS  
FX67AP  
FX67AE



76N# Upper Roots  
FX76NS



89# Upper Molars Right  
FX89S  
FX89P  
FX89E



90# Upper Molars Left  
FX90S  
FX90P  
FX90E



13# Lower Premolars  
FX13S  
FX13P  
FX13E



22# Lower Molars  
FX22S  
FX22P



EXTRACTION FORCEPS



31# Lower Roots  
FX31S



33# Lower Roots  
FX33S  
FX33P



33A# Lower Roots  
FX33AS  
FX33AP



33L# Lower Roots  
FX33LS  
FX33LP



73# Lower Molars  
FX73S  
FX73P  
FX73E




74# Lower Anteriors  
and Roots  
FX74S



74N# Lower Roots  
FX74NS  
FX74NP  
FX74NE



75# Lower Premolars  
FX75S



79# Lower Third Molars  
FX79S  
FX79P  
FX79E



45# Lower Roots  
FX45S  
FX45P



86# Lower Molars  
FX86S  
FX86P  
FX86E



86A# Lower Molars  
FX86AS  
FX86AP



86B# Lower Molars  
FX86BS  
FX86BP



86C# Lower Molars  
FX86CS  
FX86CP



87# Lower Molars  
FX87S





**EXTRACTION FORCEPS**  
AMERICAN PATTERN



American Pattern

Suggested Pair



**53R# Molars**  
First & Second Upper  
Molar Teeth Right  
F53R



**53L# Molars**  
First & Second Upper  
Molar Teeth Left  
F53L



Suggested Pair



**88R# Nevius**  
First & Second Upper  
Molar Teeth Right  
F88R



**88L# Nevius**  
First & Second Upper  
Molar Teeth Left  
F88L



**23# Cowhorn**  
First & Second  
Lower Molars  
F23



Peedo Universal



**150# Cryer Universal**  
Upper Incisors, Canines,  
Premolars & Roots  
F150



**151# Cryer Universal**  
Lower Incisors, Canines,  
Premolars & Roots  
F151



**69# Tomes**  
Upper & Lower  
Fragments & Roots  
F69



**150S# Peedo Universal**  
Upper Primary Teeth  
& Roots  
F150S



**151S# Peedo Universal**  
Lower Primary Teeth  
& Roots  
F151S



**DID YOU KNOW ?**

All instruments should be cleaned and sterilised in an open position.

Universal Forceps



221# Universal for Upper Roots  
FX221



223# Universal for Lower Roots  
FX223



EXTRACTION FORCEPS # SECURE

Secure-Line Extraction Forceps

These forceps will allow you to complete any extraction and work more efficiently even if complications arise. If teeth break during extraction, these are the tools you need to rapidly, safely solve your problems. These forceps have special, allpurpose, jaws designed to cover a broad range of applications and have sharp pointed tips that allow separating the tips of roots from surrounding tissues by lightly twisting them and then securely trapping the extracted roots in the gaps between their jaws.



845.00# Lower Roots  
SFX845.00



945.00# Lower Roots  
SFX945.00




959.01# Lower Roots  
SFX959.01



974.00# Lower Roots  
SFX974.00



849.00# Upper Roots  
SFX849.00



944.01# Upper Roots  
SFX944.01



944.00# Upper Roots  
SFX944.00



149.11# Upper Roots  
SFX149.11



897.00# Upper Roots  
SFX897.00



197.11# Upper Roots  
SFX197.11



951.00# Upper Roots  
SFX951.00



Oral Surgery

EXTRACTION FORCEPS # ATRAUMATIC

Forceps "Blade Beaks"

- "Blade Beaks" smoothly penetrate in contact with the root.
- Large contact surface for a secure grip.
- Take care of the tissues and assure tooth fracture risk reduction.
- Better control over the root during the extraction operations.



1# Upper Centrals & Canines  
AFX1



7# Upper Premolars  
AFX7



13# Lower Premolars  
AFX13



51A# Upper Roots  
AFX51A



74N# Lower Roots  
AFX74N



79# Lower Third Molars  
AFX79



Albrecht  
Extraction Forceps



1# Albrecht Premolar  
FXA1



2# Albrecht Molar  
FXA2

Special Wisdom Routerier  
Extraction Forceps



22 1/2L#  
Lower Molars & Wisdom Left  
FX22 1/2LS  
FX22 1/2LP



22 1/2R#  
Lower Molars & Wisdom Right  
FX22 1/2RS  
FX22 1/2RP



Their offset design will allow you to reach around obstacles and protect the corners of patient's mouth whenever their mouth openings are restricted or their jaws have been firmly clamped in position. Their slim compact design allows readily accessing molar areas.



**EXTRACTION FORCEPS  
BLACK LINE**

The surface shows excellent sliding features to reliably counteract contamination. An incrustation of proteins etc. is not possible. The matt-finished black surface not only prevents reflection of the OT light, but also guarantees an optimal contrast in the OR environment during extraction.



1# Upper Anterior  
FX1B



7# Upper Premolar  
FX7B




17# Upper Molar Right  
FX17B



18# Upper Molar Left  
FX18B



51A# Upper Root  
FX51AB



67A# Upper Third Molar  
FX67AB



89# Upper Molar Right  
FX89B



90# Upper Molar Left  
FX90B



13# Lower Premolar  
FX13B



73# Lower Molar  
FX73B



74N# Lower Root  
FX74NB



86# Lower Molar  
FX86B



EXTRACTION FORCEPS # PEEDO

STANDARD



1# Lower Incisors  
FX1CS



2# Lower Primary Molars  
FX2CS



3# Lower Roots  
FX3CS



5# Upper Incisors & Canines  
FX5CS



4# Upper Primary Molars  
FX4CS



6# Upper Molars  
FX6CS



7# Upper Roots  
FX7CS



PREMIUM



1# Lower Incisors  
FX1CP



2# Lower Primary Molars  
FX2CP



3# Lower Roots  
FX3CP



7# Upper Roots  
FX7CP



5# Upper Incisors & Canines  
FX5CP



4# Upper Primary Molars  
FX4CP



6# Upper Molars  
FX6CP















**PROXITIP - LUXATING HYBRIDS**

Designed with thin, precisely angled tips and a large diameter handle, the Luxating Hybrids provide improved apical and interproximal access for enhanced oral adaptation, aiding in the preservation of the buccal plate by severing the periodontal ligaments during atraumatic extractions.




**Straight Spade**  
EPTSSP  
For initial loosening  
the desmodontal  
fibres

		 <b>Small Straight</b> EPTSMS For front teeth & upper jaw	 <b>Small Curved</b> EPTSMC For lower jaw	
	<p>Suggested Pair</p>  	<p>Suggested Pair</p>  		
	<b>Small Distal</b> EPTSD For distal luxating	<b>Small Mesial</b> EPTSM For mesial luxating	<b>Large Distal</b> EPTLD For large distal luxating	<b>Large Mesial</b> EPTLM For large mesial luxating
	<p>Suggested Pair</p>  	<p>Suggested Pair</p>  		
	<b>Warwick James Left</b> EPTWJL	<b>Warwick James Right</b> EPTWJR	<b>Bernard Right</b> EPTBER6 For lingual / buccal luxating	<b>Bernard Left</b> EPTBER7 For lingual / buccal luxating



ROOT ELEVATORS



HANDLE OPTIONS



Standard



Premium

ROOT ELEVATORS

Bein



Cryer



Coupland



ROOT ELEVATORS



Root Tip Pick Heidbrink

Suggested Trio



EHB1 EHB2 EHB3 Heidbrink EHB13/14 #6 Howard EHB2/3 #6 London Hospital E20H E20HP

Root Elevators Lindo Levian



LLL1 4x8mm LLM2 3x8mm LLS3 2.5x5mm LLLC4 4x8mm LLMC5 3x8mm LLSC6 2.5x5mm

Also available in pouch set of 6 pcs. LLSP6

Warwick James Suggested Trio



EWJS EWJL EWJR

Flohr



E1F E1FP E2F E2FP E3F E3FP



HANDLE OPTIONS



**S**  
Standard



**P**  
Premium

ROOT ELEVATORS




Lacluse  
EWL

EWC12L  
EWC12LP

EWC12R  
EWC12RP

Root Elevators Premium with Cassette




- Flohr #E1FP
- Flohr Right #E3FP
- Flohr Left #E2FP
- Cryer Right #E32P
- Cryer Left #E31P
- London Hospital #E20HP
- Coupland #EA1CP
- Winter Cryer Right #EWC12RP
- Winter Cryer Left #EWC12LP
- Cassette #IMSCDRE

REPWC9

Oral Surgery

**BONE RONGEURS**

Manufactured from the finest grade materials for optimal performance and longevity. GDC features and extensive portfolio of rongeurs, including full round end cutting, end cutting and side cutting designs, engineered for trimming and recontouring alveolar bone and gross tissue removal.



**Micro Bone Rongeurs**

Small beaks for use in limited access areas and for trimming fibrous tissue or small bony spicules.





### SCISSORS

GDC scissors are fabricated from special high-tempered surgical grade stainless steel providing longer lasting cutting edges and corrosion resistance. Designed to exacting specifications for balanced and smooth control. GDC offers an extensive selection including microsurgical, suture, crown and gold, bandage and general surgery scissors.



Iris Straight 11.5cm  
S17



Iris Curved 11.5cm  
S18



Micro Iris Straight 11.5cm  
S26



Micro Iris Curved 11.5cm  
S26C



Iris Side Curved 11.5cm  
S27



Iris Angular 11.5cm  
S7



Quinby Straight 12.5cm  
S8S



Quinby Curved 12.5cm  
S8



Goldman Fox Straight 13.0cm  
S16



Goldman Fox Curved 13.0cm  
S16C



Goldman Fox Double Curved 13.0cm  
S10



Legrange Double Curved 12.0cm  
S14



Oral Surgery

SCISSORS



Kelly Straight 16.0cm  
S2



Kelly Curved 16.0cm  
S1



Locklin Straight Handle 16.0cm  
S12



Locklin Curved Handle 16.0cm  
S11



Micro Heath for Suture Cutting 10.0cm  
S25S



Heath for Suture Cutting 15.5cm  
S25



Metzenbaum Straight 14.5cm  
S28



Metzenbaum Curved 14.5cm  
S29



Mayo Straight 14.5cm  
S4



Mayo Curved 14.5cm  
S3



Crown & Band Straight 12.0cm  
SCGS



Crown & Band Curved 12.0cm  
SCGC

For Suture Cutting



Spencer 9.0cm  
S13S



Spencer Side Curved 9.0cm  
S13SA



Spencer 13.0 cm  
S13



Spencer Side Curved 13.0cm  
S13A

Marilyn Scissors



Straight 10.0cm  
S3547



Curved 10.0cm  
S3548



Straight 9.5cm  
S3541



Curved 9.5cm  
S3542



Dean Angular 16.5cm  
S9



Noyes Angular 11.0cm  
S42



Westcott Curved 11.0cm  
S35



Oral Surgery

SCISSORS



Castroviejo Straight 12.0cm  
S32



Castroviejo Curved 12.0cm  
S33



Curved

Straight



Castroviejo Straight 14.0cm  
S39



Castroviejo Curved 14.0cm  
S38



Micro Castroviejo Curved 18.0cm  
SPV



Noyes Straight 11.0cm  
S30



Noyes Curved 11.0cm  
S31



SCISSORS

Featuring tungsten carbide inserts bonded to surgical-grade stainless steel on both cutting edges for smooth, controlled cutting and maximum durability.



Legrange TC  
Double Curved 12.0cm  
S5084

Goldman Fox TC  
Double Curved 13.0 cm  
S5079



Spencer TC Angular 10.0cm  
S3481



Spencer TC Angular 13.0cm  
S3482



Spencer TC Straight 10.0 cm  
S3483



Spencer TC Straight 13.0 cm  
S3484



Iris TC Straight 11.5cm  
S5082



Iris TC Curved 11.5cm  
S5083



Goldman Fox TC Straight 13.0 cm  
S5080



Goldman Fox TC Curved 13.0 cm  
S5081



Metzenbaum TC Straight 12.0cm  
S5054



Metzenbaum TC Cruved 12.0cm  
S5055



Crown & Band TC Straight 12.0cm  
S5038



Crown & Band TC Curved 12.0cm  
S5039

Oral Surgery

HEMOSTATS

Each hemostat is precision engineered from surgical grade stainless steel to exacting specifications featuring a secure ratchet and various working and geometries. Designed to optimise clinical efficiencies as a multipurpose instrument, the GDC hemostats are used for clamping off blood vessels, removing small root tips and grasping loose objects.



Kelly Straight 14.0cm  
H1



Kelly Curved 14.0cm  
H2



Crile Straight 16.0cm  
HRPS



Crile Curved 16.0cm  
HRPC



Mosquito Straight 12.0cm  
H4



Mosquito Curved 12.0cm  
H3



Mosquito Straight 1x2 12.5cm  
H5



Mosquito Curved 1x2 12.5cm  
H6



Baby Mosquito Straight 10.0cm  
HHS



Baby Mosquito Curved 10.0cm  
HHC

Allis Forcep



Allis Baby 12.5cm  
TFB



Allis 15.0cm  
TFS

NEEDLE HOLDER

Used to grasp and guide the needle when suturing



Mayo-Hegar Straight 16.0cm  
NHMH



Mayo-Hegar Curved 16.0cm  
NHMHC



Crile-Wood 15.0cm  
NHCW



Adson Straight 17.0cm  
NHA



Needle Holder 18.0cm  
NHMH18

(Also available in 20.0cm #NHMH20)



Olsen-Hegar Needle Holder  
with Scissors 17.0cm  
NHOH



Derf Straight 12.5cm  
NHD



Derf Curved 12.5cm  
NHDC

Castroviejo  
Needle Holder



14.0cm  
Curved  
NHCVC



14.0cm  
Straight  
NHCVS

Backhaus Towel Clamp



Backhaus 9.0cm  
TC3



Backhaus 11.0cm  
TC4



Backhaus 13.0cm  
TC5



Oral Surgery

NEEDLE HOLDER TUNGSTEN CARBIDE

Needle holder feature tungsten carbide insert with a ground pyramidal profile to assure firm, positive needle grip. Rotation of needle and /or slipping of suture material is virtually eliminated



Mayo-Hegar 16.0cm  
NH5042



Crile-Wood 15.0cm  
NH5038



Halsey 13.0cm  
NH5037



Derf 11.5cm  
NH5032

Sterilising Forceps



Cheatele Forcep 27.0cm  
CF

Micro Castroviejo Needle Holder Tungsten Carbide Insert



Castroviejo TC Straight 14.0cm  
NH5020



Castroviejo TC Curved 14.0cm  
NH5021



Castroviejo TC Straight 18.0cm  
NH5024R

Dressing Forceps



Gross Maier Straight 20.0cm  
SDFS



Gross Maier Curved 20.0cm  
SDFC



Foerster Straight 20.0cm  
SDFP



Rogge 20.0cm / 30.0cm  
RF20  
RF30

TISSUE FORCEPS

Used to grasp and stabilize soft tissue flaps during suturing.



Adson Serrated 15.0cm  
TP45



Adson 1x2 15.0cm  
TP46



Adson Serrated 12.0cm  
TP41



Adson 1x2 12.0cm  
TP42



Adson Curved 1x2 12.0cm  
TP44



Semken Serrated  
Straight 12.5cm  
TP31



Semken Serrated  
Curved 12.5cm  
TP32



Semken Straight  
1x2 12.5cm  
TP33



Semken Curved  
1x2 12.5cm  
TP34

Tissue  
Forcep Fine



12.5cm  
TP35



11.5cm  
TP36



11.5cm  
TP37

Oral Surgery

MICRO TISSUE FORCEPS



DENTAL SYRINGE

GDC syringes are available in two styles, the cooke-waite and with large diameter finger rings for enhanced control, large window for enhanced carpule visibility and sharp harpoons for easy cartridge engagement.

Dental Syringe with colour coded ring



SYRACC1.8R  
(Also available in different colours of ring)



Self Aspirating Syringe



SYRSA1.8

Intraligamental Syringe - Inject



SYRI1

"Suitable for upper and lowers, just turning the front barrel"



Not Aspirating

SYR1.8  
SYR2.2

Aspirating

SYRA1.8  
SYRA2.2

Aspirating

SYRF1.8

(Also available in European & American thread needle type)

ADVANTAGES OF INTRALIGAMENTAL ANESTHESIA :

- Automatic limitation of the power pressure.
- Anesthetizing of a single tooth.
- No numbness of the tongue, lips or cheek.
- Simultaneous treatment of teeth in different quadrants.
- Supports regional block anesthesia.
- Eliminates the risk of haematoma.
- Lower emotional impact on patient.

WHY INTRALIGAMENTAL IS BETTER


Because it's chrome free.  
Because you can get a better sterilization.  
Because you can use for long time  
Because has anti-reflection surface.

OPERATING FEATURES:

- Fully autoclavable.
- Removable headpiece with angled threaded nozzle for the injection needle.
- Viewing chamber with safety sleeve for carpule cartridge.
- Handle servo dosing lever with optimal support for index finger for a gradual advance of the dosing plunger resetting key.
- No waiting time between injections and beginning of treatment.
- Excellent depth of anesthesia combined with shorter duration.
- Much less pain than standard injections (infiltration anesthesia in hard palate, intrapulpal or intraseptal injections or for papillectomy).
- A smaller quantity of given anaesthetic (~0.06 ml for single injection/click.)
- Needle suggested: 30G.



SUCTION TUBES & ADAPTOR



BONE CHISEL

Used to scoop bone and make trough




MALLETS

Used with a chisel to split teeth and remove bone



BONE FILE

Bone File used for final smoothing of bone



SCALPEL HANDLES

GDC scalpel handles feature balanced designs and are manufactured from the finest grade surgical stainless steel for optimal control and enhanced longevity.



Scalpel Handle Straight 14.5cm  
10-130-5EM



Scalpel Handle Angular 14.5cm  
10-130-5AE



Scalpel Handle Straight 14.5cm  
10-130-5EM #4



Scalpel Handle Angular 14.5cm  
10-130-5AE #4



Bard Parker No. 3 12.5cm  
10-130-03



Bard Parker Style with  
Metric Ruler No. 3 12.5cm  
10-130-03E



No. 4 13.5cm  
10-100-04E



Orientable Scalpel Handle  
10-130-180



No. 7 16.5cm  
10-130-07E



Moon's Probe  
P13M



Double Side Blade  
Scalpel Handle  
10-130-5DE

Oral Surgery

CHEEK RETRACTORS / TISSUE RETRACTORS

Used to hold mucoperiosteal flaps, cheeks, lips, and tongue away from the surgical area.

Austin Cheek Retractors



Mouth Gags



Heister 13.0cm  
HMG



Large 14.0cm  
MGA



Small 11.0cm  
MGC

Hartmann Retractor



15mm 72x24mm  
CRH

Weider Cheek & Tongue Depressor



Large 15.0cm  
CRW2

Small 14.0cm  
CRW1

Orringer Retractors



60x115mm Large  
CRO3

Also Available in CRO2 Medium & CRO1 Small

CHEEK RETRACTORS / TISSUE RETRACTORS



Bishop 13.5cm  
CRB



Vestibulum 14.0cm  
CRV



Martin Simplex 10.5cm  
CRMS



MSY 11.5cm  
MSY



Columbia 14.0cm  
CRC



Sternberg 16.0cm  
CRCS



Bruenings 19.0cm  
TRB



University of Minnesota  
14.0cm  
CRM

Cheek Retractor with Suction



Fork  
TRA2S



Half Moon  
TRA3S



With Suction Large 15.0cm  
CRW2S




University of Minnesota  
with Suction 14.0cm  
CRMS



CHEEK RETRACTORS / TISSUE RETRACTORS



Langenbeck 21.5cm



Middledorf 19.0cm



Ramus 21.5cm



Kocher 21.5cm



Kocher 21.0cm

CHEEK RETRACTORS / TISSUE RETRACTORS



Oral Surgery

PERIOSTEAL ELEVATOR

For reflecting and retracting the mucoperiosteum after incisions of the gingival tissue.



Hirschfeld  
P20  
#6



Glickman  
P24G  
#6



Kramer-Nevins  
PKN1  
#6



Prichard  
PPR3  
#6



Benque  
PBQ  
#6



Buser  
PPBUSER  
#6



Seldin  
P23



Goldman Fox  
P14  
#6



P149  
#6



P8D  
#6



PH2  
#6



Howarth  
PH

PERIOSTEAL ELEVATOR



PPLABANCA2  
#6

P856  
#6

Molt  
P883  
#6

PHP3  
#6

Freer  
P884  
#6



Molt-9 Large  
P9  
#3

Molt-9 Small  
P9S  
#3

Freer  
P16  
#3

PPLABANCA2  
#3

BEAL1  
#3

Kramer Naviens  
PKN  
#3





# Oral Maxillofacial

Retractors

Trauma Kit Instruments

Cleft Lip & Palate Instruments

Impaction Kit Instruments



BONE HOLDING, RETRACTORS & OTHER INSTRUMENTS

Asch Nasal Septum Forcep



Large ANS

Walsham Nasal Septum Forcep



Left WNSL

Right WNSR

Rowe Maxillary Disimpaction Forcep



Left RMDL

Right RMDR

Nasal Chisel with Guard



Left NCGL

Right NCGR

Ferguson Mouth Gag



Large FMGA

Small FMGC

Hayton William Forcep



Left HWL

Right HWR

Oral Maxillofacial

BONE HOLDING, RETRACTORS  
& OTHER INSTRUMENTS



Bone Holding Forcep (Mandible)  
BHM



Channel Retractor Regular  
CHR



Kilner Retractor C Shape  
KRC



Chin Retractor  
CR



Perigoidous Chisel  
PC



Nasal Septum  
Osteotome with  
Guard U Shape  
NSOU



Bristow Elevator  
BE



Bone Elevator  
Hohman Type  
BEH



Micro Miller Bone  
File Straight  
BFMS



Barsky Pharyngeal  
Flap Elevator  
BPF

BASIC BONE PLATING INSTRUMENTS



Modelling Lever  
ML



Plier  
BPP



Bone Cutting Forcep Single Action  
RBC



Zygomatic Bone Awl  
ZBA



Ordinary Screw Driver  
OSD1.5  
OSD2.0  
OSD2.5



Self Holding Screw Driver  
SSH1.5  
SSH2.0  
SSH2.5




Wire Twister  
WTL



Bone Holding Forcep Pointed  
BHP



Trauma Kit  
TRAUMA



SS Drill Bite for Micromotor  
DSS1.5  
DSS2.0  
DSS2.5



2 Holes  
P2H1.5SS  
P2H2.0SS  
P2H2.5SS



4 Holes  
P4H1.5SS  
P4H2.0SS  
P4H2.5SS



Screw  
S1.5SS  
S2.0SS  
S2.5SS



L SHAPE



Right 4 Holes  
P4HLR1.5SS  
P4HLR2.0SS



16 Holes  
P16H1.5SS  
P16H2.0SS  
P16H2.5SS

Orbital 4 Holes  
P4HO1.5SS  
P4HO2.0SS

Left 4 Holes  
P4HLL1.5SS  
P4HLL2.0SS



CLEFT LIP & PALATE INSTRUMENTS



Palate Dissector Right  
CPDR



Palate Dissector Left  
CPDL



Cleft Palate  
Raspatory Up  
CPRU



Cleft Palate  
Raspatory Down  
CPRD



Senn - Miller Retractor (Kilner)  
TRSM1



McIndoe Cleft Palate Raspatory  
MCPR



Cleft Palate Elevator  
CPE2



Spoon Shape Elevator  
SSE



Cleft Palate Elevator  
CPE1



Miller Bone File Curved  
BF52

CLEFT LIP & PALATE INSTRUMENTS



Dingman Mouth Gag with 4 Baldes  
DMG



Cleft Palate  
Down Hook  
CPHD

Cleft Knife  
Triangle  
CPKT

Cronin Cleft  
Palate Elevator  
CCP

Cleft Palate Raspatory  
French Type  
CPRFL  
CPRFR

Howarth  
PH

Michel Trimmer  
TM

Cleft Palate  
Raspatory Barsky  
CPRB

IMPACTION KIT



Mouth Props with Chain (Set of 3)  
MP3



IMPACTION



London Hospital  
E20H



Coupland  
EA1C



Chisel 7.5mm  
OSS6520S



Cryer  
E31  
E32



Root Tip Pick Heidbrink  
EHB2  
EHB3



Warwick James  
EWJS



Blumenthal  
RBL



Heath Suture  
Cutting 15.5cm  
S25



Crile Straight 16.0cm  
HRPS



Mayo-Hegar  
Straight 16.0cm  
NHMH



Backhaus 13.0cm  
TC5



Adson 1x2 15.0cm  
TP46



Fork  
TRA2



Scalpel Handle  
Straight 14.5cm  
10-130-5EM



Curette Lucas  
CL87  
#3



Miller 52  
BF52



Molt-9 Large  
P9  
#3



Frazier Suction  
ASPFR14



Howarth  
PH

# LUXATIP



### Also available in dual edge blades

- L-2S Straight Blade- Green-2 mm- Apical and/or Interproximal
  - L-3C Curved Blade- Black- 3 mm General use/Starter
  - L-3S Straight Blade- Grey- 3 mm Interproximal
  - L-3CA Contra Angle- Blue- 3 mm Lingual and/or distal-Molars
  - L-3IC Inverted Curved- Light Blue- 3 mm-General, Lingual and/or Distal
  - L5C Curved Blade- Brown- 5 mm- Large molar roots-General use
  - L5S Straight Blade-Yellow-5mm- Interproximal
  - DE3 dual edge/straight grey 3 mm/1.5 mm
  - DE5 dual edge/straight brown 5 mm/3 mm
- Deep fractured/Decayed Roots





# Implantology

Osteotomes

Sinus Lift Instruments

Caliper & Chisel

Bone Injector & Collector

Bone Mill & Crusher




Osteotomes Bone Pushers Straight

Allows for the transport of bone graft chips from the implant cavity or for condensing them apically.



Osteotomes Bone Pushers Angulated

Designed to condense bone material and/or autogenous bone under the sinus. Minimise the risk of perforation of the mucous membrane in sinus elevation procedures.




OSTEOTOMES CONVEX & CONCAVE

Osteotomes Tapered Convex Tip Angulated

Specially designed for closed sinus elevation operations in the lateral upper jaw region or in difficult to access maxillary areas.



Osteotomes Tapered Concave Tip Angulated





BONE MILL & SCRAPER

**Bone Mill**

**For gaining autogenic bone material.**

In order to correct bone defects, an intraoral harvesting of bone is often necessary. Pieces of bone may be gained with trephines or chisels. It is the aim of the Bone Mill Forceps to cut bone material into small pieces so that a maximum volume of bone chips can be achieved. The "grainy" consistency of bone chips created will ensure that it is adapted most favourably to osseous embedding.



**Bone Crusher**  
IMPBMS



**Bone Cracker for Large Bone Pieces**  
IMPBMR



IMPBSB2 IMPBSB1

**Bone Scraper Blade**



IMPBS IMPBC

**Bone Scraper**

**Bone Scraper**

**Collecting autologous bone material.**

In oral surgery bone replacements and bone structures are often combined with autologous bone pieces during augmentations in order to use their osteoinductive effect. The Bone Scraper is the perfect instrument for an easy and quick collection of autologous bone structures. It enables you to scrape, collect and transplant the patient's own bone. Everything can be done without the usage of bone filter, trephine bur, saw or bone mill.

The blade makes a collection of cortical and spongiöse bone possible. The bone is collected in a chamber during the scraping and at the same time the bone is being mixed with blood. The bone material can be implanted directly out of the Bone Scraper which has been sterilized before.

**DID YOU KNOW ?**  
Bone Scraper Blades can be purchased seperately.



IMPLANTOLOGY CALIPER & BONE COLLECTOR

Bone Caliper



BC15  
15.5cm



BC35  
15.5cm

Bone Collector



7.0mm  
BC1331/1F

BC1331/1  
18.5cm



12.0mm  
BC1331/2F

BC1331/2  
18.5cm

Castroviejo Caliper



CLC40L  
18.0cm



CLC40C  
18.0cm



CLC20C  
8.5cm



CLC20L  
8.5cm



Trephine Size (in mm) :  
3 | 4 | 5 | 6 | 7 | 8 | 9 | 10

Trephine Kit  
IMPTS8



Drill Size (in mm) :  
2 | 2.5 | 2.8 | 3.2 | 3.65 | 4.2 | 4.8 | 5.2

Drills Kit  
IMPDS8

TISSUE PUNCH



TPC4  
10.0cm

TPC5  
10.0cm

TPC6  
10.0cm

TPS4  
10.0cm

TPS5  
10.0cm

TPS6  
10.0cm



Instrument to measure the Depth of bony tunnel

PCP1301  
Tip Size : 1.2mm

Bone Graft Condenser  
BC14614

Bone Graft Carrier  
BC14615

Bone Compactor & Carrier  
BCC14616

Membrane Placement  
PMPIS



PRF BOX  
Palatelet Rich Fibrin  
PRF

OSTEOTOMES BONE CHISEL




OCHSENBEIN CHISEL



BONE PUNCH



BONE INJECTORS



MALLETS

Used with a chisel to split teeth and/remove bone



### SINUS LIFT INSTRUMENTS

Featuring a wide range of sinus lift instrument designs, GDC offers a variety of size, working end and angulation options to provide optimal access to separate and reflect the schneiderian membrane from the maxillary bone and to elevate the membrane during lateral at sinus lift procedures.



IMPSL1

IMPSL2

IMPSL3

IMPSL4

IMPSL5



IMPSL6

IMPSL7

IMPSL8

IMPSL9



Micro Surgery

MICRO SURGERY INSTRUMENTS



Castroviejo Needle Holder 18.0cm NH5024R



Castroviejo Scissor 18.0cm SPV



Micro Surgery Instruments With Cassette MSIWC6



Micro Tissue Forcep Straight TPSSTMBH



Micro Scalpel Handle 12.0cm 1-015



Papila Elevator PH26M



Glickman P24G



# Endodontic

Composite Instruments

Spoon Excavator

Plugger / Condenser / Burnisher

Root Canal Plugger & Spreader

Rubber Dam Clamp, Forcep & Punch





COMPOSITE INSTRUMENTS

Plasma Black



CIPB1  
#8

CIPB2  
#8

CIPB3  
#8

CIPB4  
#8

CIPB5  
#8

CIPB6  
#8

Titanium Blue



CIB1  
#8

CIB2  
#8

CIB3  
#8


CIB4  
#8

CIB5  
#8

CIB6  
#8

COMPOSITE INSTRUMENTS

Titanium Gold



CIG1  
#6

CIG2  
#6


CIG3  
#6

CIG4  
#6

CIG5  
#6

CIG6  
#6

Titanium Blue



CIB1  
#6

CIB2  
#6

CIB3  
#6

CIB4  
#6


CIB5  
#6

CIB6  
#6



COMPOSITE INSTRUMENTS

Titanium Blue # Anterior



TNCIGFT3  
#6

TNCIPCL  
#6


TNCIPCM  
#6

TNCIPCS  
#6

TNCIGFTMI1  
#6

Also available in Plasma Black

Titanium Blue # Posterior



TNCFIS/M  
#6

TNCFIM/L  
#6

TNPLG5A  
#6

TNPLGOT  
#6

TNPLGH3  
#6

Also available in Plasma Black

COMPOSITE INSTRUMENTS

Titanium Goldstein Flexi Thin



CIGFT1  
#6

CIGFT2  
#6

CIGFT3  
#6

CIGFT4  
#6

CIGFTMINI3  
#6

CIGFTMINI1  
#6

Anterior Kit

**TNCIGFT3**  
Flexible, reversed, flared paddle design for shaping and placement of Class III and IV restorations.

**TNCIPCL**  
For final placement in Class I and II restorations. The larger, round ball end is used for condensing and shaping in Class I and II restorations and on lingual surfaces of anterior teeth.

**TNCIPCM**  
For small pits and fissures, as well as placement and condensing with limited access.

**TNCIPCS**  
Micro-Mini for extremely small pits and fissures.

**TNCIGFTMI1**  
Mini versions of the TNCIGFTM1 for small pits and fissures, tunnel preparations or minor tooth defects on lower anteriors.

Posterior Kit

**TNCFIS/M**  
Rounded cone-shaped paired instrument designed to provide improved contact forming for small/medium Class II restorations.

**TNCFIM/L**  
Rounded cone-shaped instrument to provide improved contact forming for medium/large Class II restorations.

**TNPLG5A**  
Small, round, inverted-cone plugger for use with condensable composite material in posterior restorations.

**TNPLGOT**  
Rhomboid-shaped plugger for use with condensable composite material in posterior restorations.

**TNPLGH3**  
Rectangular-shaped plugger for use with condensable composite material in posterior restorations.

Gold Flexi Thin

**CIGFT1**  
Small universal style with rounded plugger tip and a narrow paddle for initial placement and contouring of Class I, II and III restorations.

**CIGFT2**  
Large universal style for final placement and contouring of Class I, II and III restorations.

**CIGFT3**  
Flexible, reversed, flared paddle design for shaping and placement of Class III and IV restorations.

**CIGFT4**  
Flexible, paired, offset, paddle-shaped bladed for placing and shaping material on posterior, mesial and distal surfaces. Reverse angle is also useful for placing and shaping anterior bonded restorations.


**CIGFTMINI3**  
Mini version of the CIGFT3. Can also be used for packing gingival retraction cord.

**CIGFTMINI1**  
Mini version of the CIGFT1 for small pits and fissures, tunnel preparations or minor tooth defects on lower anteriors.


Endodontic

SPOON EXCAVATOR

GDC endodontic excavators feature extended length shanks and are meticulously hand crafted from the highest grade stainless steel. They help the clinician access the pulp chamber for easy and efficient removal of dentin. The result is a superior clinical outcome and improved efficiency.



Oval Spoon



SPOON EXCAVATOR



EXC31L  
#3 #6

EXC32L  
#3 #6

EXC33L  
#3 #6

EXC31W  
#6

EXC816  
#3

EXC817  
#3

EXC819  
#3

Darby-Perry



EXC820  
#3

EXC210/211  
#3 #6

EXC212/213  
#3

EXC153/154  
#3 #6

EXC38/39  
#6

EXCGL2  
#6



Endodontic

PLUGGER / CONDENSERS &  
PLASTIC FILLING INSTRUMENTS

Used to pack filling materials into cavity preparations.



PLG1  
#3 #4

PLG2  
#3 #4

PLG3  
#3 #4

PLG4  
#3 #4

PLG1L  
#3 #4

PLG2L  
#3 #4

PLG3L  
#3 #4

PLG4L  
#3 #4

Plastic Filling Instruments



PFGFT1  
#3

PFGFT2  
#3

PF2  
#3 #4

PF21  
#3 #4

PF11  
#3

PF12  
#3

PF13  
#3

PLGH3  
#3 #4

**BURNISHER & CARVERS**

Designed to condense, smooth, carve and polish amalgam



BB21B  
#3 #4

BB21  
#3 #4

BBL3  
#3

PKT3  
#3 #4

BB22/23  
#3 #4

BB18  
#3 #4

BB27/29  
#3 #4

TB2  
#3 #4

**Carvers**

Used to carve anatomical features and trim excess materials.



Discoid-cleoid  
CD89/92  
#3

Discoid-cleoid  
CD4/5  
#3 #4

Interproximal  
Off-angle  
CVIPCOA  
#3

Hollenback  
CVHL3S  
#3 #4

Wall W-Shaped  
CVHL3W  
#3

Ward's  
CVWR1S  
#3 #4

Diamond Shape  
Frahm  
CVFR2/3H  
#3 #4

Diamond Shape  
Frahm  
CVFR2/3  
#3

**MARGINAL TRIMMER**

Used to produce proper bevel on enamel margins.




MT26H  
#6

MT27H  
#6

MT28H  
#3 #4 #6

MT29H  
#3 #4 #6

**ENDO GUNS FOR RETROGRADS APPLICATIONS SYRINGES**




MTA4865  
13.0cm

MTA4866  
13.0cm

MTA4869  
13.0cm

**PLACEMENT INSTRUMENTS**

Used to place base or liner within cavity preparations.



Calcium hydroxide or glass ionomer base/liner placement instrument. Also useful as a small burnisher.

Calcium hydroxide or glass ionomer base/liner placement instrument combined with a mini-spatula for efficient mixing.

Calcium hydroxide or glass ionomer base/liner placement instrument.

PICH6  
#3

SP60616  
#3

PICH7  
#3

PICH8

CHISELS AND HOE

Used to refine and cavity preparations. Forming lines angles on anterior preparations.

Enamel Chisel



Wedelstaedt



Hatchet



Hoe



CP7 #3 #4

CP7/8 #3

CP1/2 #3

CP3/4 #3

CP5/6 #3

CP53/54 #3 #4 #6

CP24 #3 #4

MATRIX RETAINER



Ivory MRI



Tofflemire Universal MRTU



Squivland Narrow MR4972S



Squivland Wide MR4972L

AMALGAM CARRIER

To carry and dispense amalgam filling materials.



AC504



Large/Regular AC5202



Jumbo/Regular AC5203



Mini/Regular AC5201



Amalgam / Bone well WA

Stainless steel, with removable non-slip base ring. Designed for easy amalgam/bone graft handling.





Endodontic

ISO\*SIZED PLUGGER COLOUR CODED

When performing vertical condensation, choose from ISO-sized anterior or posterior. These pluggers are handcrafted to exact specifications using the finest materials.

**\* DID YOU KNOW ?**  
 ISO Pluggers have an ISO colour-coded ring on the handle for ease of identification during usage.  
 #30 Blue #40 Green #50 Yellow #60 Red



### SPREADERS, HEAT CARRIERS

ISO sized spreaders with color coded handles for easy identification, heat carriers, and a variety of other popular spreaders. Our spreaders help the clinician gain access into the most difficult canals for optimal compaction of filling material during lateral condensation.



### NERVE-CANAL PLIERS

Nerve-Canal Pliers for grasping fractured root-canal instruments or silver pins. The very fine, concave cutout on the working tip allows to remove even deep-lying fractured files.



RFES

### ROOT FRAGMENT STIEGLITZ

Post and Silver Point Removal forceps



RFSTR

RF45

RF90



RFEC

Endodontic

RUBBER DAM

Rubber Dam Clamps # Gold Line



2A  
RDCG2A



3  
RDCG3



4  
RDCG4



12A  
RDCG12A



13A  
RDCG13A



14A  
RDCG14A



201  
RDCG201



202  
RDCG202



205  
RDCG205



210  
RDCG210



211  
RDCG211

Rubber Dam Clamps # Black Line



2A  
RDCB2A



3  
RDCB3



4  
RDCB4



12A  
RDCB12A



13A  
RDCB13A



14A  
RDCB14A



201  
RDCB201



202  
RDCB202



205  
RDCB205



210  
RDCB210



211  
RDCB211

POINTS OF PERFORMANCE

- Which resists corrosion while aiding in flexibility and memory
- A dark coating allowing for improved intraoral contrast and visual acuity
- Matt finish for decreased light reflection and glare



Rubber Dam Clamps # Matt Finish



2A  
RDC2A



3  
RDC3



4  
RDC4



12A  
RDC12A



13A  
RDC13A



14A  
RDC14A



201  
RDC201



202  
RDC202



205  
RDC205



210  
RDC210



211  
RDC211

Rubber Dam Clamps Set



Gold  
RDCOGB11



Black  
RDCOBB11



Matt Finish  
RDCOB11

Rubber Dam Clamps Peedo # Matt Finish



A1  
RDCA1



A2  
RDCA2



A3  
RDCA3



A4  
RDCA4



A5  
RDCA5



A6  
RDCA6

DESCRIPTIONS

2A  
Flat jaws for larger premolars.

3  
Flat narrow jaws for small lower molars.

4  
Flat narrow jaws for small upper molars.

12A  
Serrated offset jaws for lower right molars. Paired with 13A

13A  
Serrated offset jaws for lower left molars. Paired with 12A

14A  
Deeply festooned jaws for partially erupted for irregularly shaped molars.

201  
Slightly festooned jaws for upper molars.

202  
Universal molars.

205  
Slightly festooned jaws for large upper molars.

210  
Universal double bowed anterior clamps.

211  
Universal double bowed anterior clamps.



Winged Clamps



Wingless Clamps



RUBBER DAM CLAMP, FORCEP & PUNCH



Ivory 17.0cm  
RDF



Stoke 17.5cm  
RDF4E



Ainsworth 17.0cm  
RDPA



Rubber Dam Punch 16.0cm  
RDP



Matrix Ring Holding Forcep 17.0cm  
MRDF



Dental Dam Kit Adult DDK  
Dental Dam Kit Peedo DDKP


RUBBER DAM

Plastic Frame



RDAF6P

Stainless Steel Frame



Large 6"  
RDAF6

Medium 5"  
RDCF5

Small 4"  
RDCF4



Tooth Separator Elliot  
TSE



Tooth Separator Ivory  
TSI

Rubber Dam Sheet  
6" x 6" - RDLA  
5" x 5" - RDLC  
Thickness - Fine / Medium / Heavy  
Colour - Green / Blue / Black

Template 6" x 6"  
RDT  
Adult template shown.  
Pedodontic template  
also included.

ART Instruments



ART1

ART2

CP53/54  
#3 #6

EXC125/126  
#3 #6

EXC133/134  
#3 #6

The State of the art  
- ART Instruments

ART restoration step-by-step using a high-viscosity glass-ionomer (Ketac Molar EasyMix). Dentine carious cavity in tooth LLE.

(a) Notice the discoloration around the cavity opening. It indicates that the carious lesion has extended under the enamel. This unsupported enamel is demineralized and will break off easily under slight pressure.

(b) Opening of cavity further for access with the Enamel Access Cutter (ART instrument).

(c) Carious tissue removal using a small and sharp excavator.

(d) The conditioner is applied to the cleaned cavity, pits and fissures with a cotton wool pellet.

(e) Carefully dried cavity before placing the high-viscosity glass-ionomer.

(f) Firm finger pressure has been applied over the occlusal surface and excess filling material has been removed with a carver instrument. The ART restoration presented after the bite has been adjusted showing the sealant restoration: the cavity is restored, fissures are sealed.



# Periodontal

Curettes


Scalers

Probes

Knives & Chisels

Periotomes





STANDARD

**Gracey Curette Standard**  
Originally designed to be finishing curettes for fine scaling and root planing. The combination of this unique blade and 9 different shank designs for specific tooth surfaces provide improved adaptation and deposit removal in scaling procedures.



MINI FIVE

**Gracey Curette Mini Five**

- Terminal shank is 3mm longer than Standard Gracey
- Compared to the standard Gracey curette, blade is 50% shorter and 10% thinner
- Good for scaling in deep, narrow pockets

HANDLE OPTION




#3

#6

#7

GRACEY CURETTES # STANDARD



SG1/2  
#3 #6

SG3/4  
#3 #6

SG5/6  
#3 #6

SG7/8  
#3 #6

SG9/10  
#3 #6

SG11/12  
#3 #6

SG13/14  
#3 #6


SG15/16  
#3 #6

SG17/18  
#3 #6




CURETTES

GRACEY CURETTES COLOUR CODED # STANDARD



Also available Gracey Currettes Colour Coded # RIGID

GRACEY CURETTES COLOUR CODED # MINI FIVE




Mc CALL CURETTES



SM11/12 #6    SM11/12S #6    SM11/12A #6    SM13/14 #6    SM13/14S #6    SM17/18 #6    SM19/20 #6


LANGER CURETTES



SL1/2 #6    SL3/4 #6    SL5/6 #6

Langer 1/2 for mandibular posterior.  
Langer 3/4 for maxillary posterior.  
Langer 5/6 for maxillary and mandibular anterior.

COLUMBIA CURETTES



SC2R/2L #3 #6    SC4R/4L #3 #6    SC2R/2L #7    SC4R/4L #7

Columbia 2R/2L anterior application  
Columbia 4R/4L posterior application

Derby Perry  
Chisel

MANIPAL SCALERS



Posterior  
S129  
#3 #6

Hoe Orban  
SO7/9  
#3

Bifid  
SB4/8  
#3 #4

Derby  
SDP5/6K  
#6

SM137


SM2/3

SM5/6

TUNNELLING KNIVES

FURCATION CURETTES

Designed for scaling root concavities and furcation.



Tunneling Knife 1

Tunneling Knife 2

Quetin Buccal/Lingual

Quetin Mesial/Distal

TKN1  
#6

TKN2  
#6

SQBL2  
#6

SQMD2  
#6


The shank of the TKN1 has a smaller angle so it is suitable for tunnelling work in less complex areas.

The small radius shank of the TKN2 allows incision in convex areas with relatively large undercuts.

SUB GINGIVAL SCALERS



SUPRA GINGIVAL SCALER





PERIO PROBES

Used to assess periodontal pocket depths, attachment levels, anatomy configurations and gingival bleeding.



Thin William  
POW6  
#3 #4 #6


Goldman Fox  
PGF  
#6

University of North Carolina  
PCPUNC15  
#3 #6


Marquis Colour Coded  
PCP12  
#6

Screening Probe  
PCP11.5B  
#3 #6

PERIO PROBES



EXP23/QW6  
#6



EXP23/UNC15  
#6



EXP23/11.5  
#6



University of North  
Carolina/William  
PCPUNC15/W6  
#3 #6



Naber's  
Colour Coded  
PQ2N  
#3 #6 #7



Thin William/  
Screening Probe  
PCPQOW/11.5B  
#3 #6



Goldman Fox/  
William Probe  
PCPGF/W6  
#6



UNC15/  
Screening Probe  
PCP15/11.5B  
#6



CP12/  
Thin William  
PCP12/QOW6  
#6



CP12/  
Screening Probe  
PCP12/11.5B  
#6

### HIRSCHFELD FILE

Used to crush and remove heavy calculus deposits.



**Buccal/Lingual**  
FH3/7  
#6

**Mesial/Distal**  
FH5/11  
#6

**Anterior**  
FH9/10  
#6

### PERIODONTAL FILES

Used interproximally. Files surfaces on both sides allow for push or pull application.



**Sugarman**  
FS1/2S  
#6

**Sugarman**  
FS3/4S  
#6

**Schluger**  
FS9/10S  
#6

For removal of granulation tissue in supra-bony and infra-bony pockets.

### OCHSENBEIN CHISEL

Useful for reshaping bone. The semi-circular notches are beveled on each side. Can also be used to reflect flaps or remove secondary palatal flaps.



CO1

CO2



CO3  
#6

CO4  
#6

### RHODES BACK ACTION

Designed for recontouring bone. Particularly valuable for removal of 'widows' peaks on difficult to reach distal surfaces of posterior teeth.



C36/37  
#6

Periodontal

PERIODONTAL KNIVES

Goldman-fox



KGF8  
#6

KGF9  
#6

KGF11  
#6

Buck



KB3/4  
#6

KB5/6  
#6

Crane Kaplan



KCK3  
#6

Corn Suture Plier



SP20SLEM  
18.0cm

SPM20  
15.0cm

SP20

Orban



KO1/2  
#6

KO12KP03A  
#6  
30% thinner

Kirkland



KK15/16  
#6



Tissue Nipper



NIPS

Paquette  
Blade Handle



KPO


Goldman Fox  
Pocket Marker



Left  
PMGF1

Right  
PMGF2

Surgical Curette Lucas



CL84 #3 #4 #6    CL85 #3 #4 #6    CL86 #3 #4 #6    CL87 #3 #4 #6    CL88 #3 #4 #6

Bone Curette Hemingway



CH876 #3 #4    CH877 #3 #4    CH878 #3 #4    CH879 #3 #4

Periodontal

PERIOTOMES # FLEXI


This product is used for minimizing damage to the alveolar bone as it helps cut the periodontal ligament to detach teeth from the periodontal ligament for implant immediately after tooth extraction or general tooth extraction.

CAUTIONS

PERIOTOME FLEXI is used for cutting off only the periodontal ligament. Do not lift it up or twist with excessive force such as ELEVATOR.




PERIOTOMES NON FLEXIBLE



### VISTA TUNNELING PROCEDURE

The first incision is made at the base of the sulcuso, then sever the attachment from the osseous crest raise the mucoperiosteal pouch with a sub periosteal dissection and follow with a curette to separate the Vista from the osseous crest.



VISTA1

VISTA2

VISTA3

VISTA4

VISTA5

VISTA6

- Corn Suture Plier #SP20
- Scaler #S204S
- Surgical Curette Lucas #CL86
- Castroviejo Scissor Curved #S33
- Scalpel Handle Angular #10-130-5AE
- Perio Probe Thin William #POW6
- Suction Tubes Coupland #ASPCT3
- Periosteal Elevator Benque #PBQ
- Micro Tissue Forceps Gerald Straight Serrated #TPG1
- Micro Tissue Forceps Gerald Curved Serrated #TPG2
- Needle Holder Mayo-Hegar Straight #NHMH
- Tissue Forcep Adson 1x2 #TP42
- Instrument Management System #IMSCA20E







# Prosthodontic

Wax & Modelling Carvers

Cement Spatulas

Wax Knives & Spatulas

Impression Trays

Crown Removers





WAX AND MODELLING CARVERS



Zahle  
WCZ1

Zahle  
WCZ2

Zahle  
WCZ3

Hyllin  
WCH1

Hyllin  
WCH2


Hyllin  
WCH3

Gritman  
WCG


Lecron  
WCL1

Lecron  
WCL4

PKT CARVERS COLOUR CODED



PKT CARVERS



WAX KNIVES & SPATULAS



Fahnenstock  
WKFS



Lessmann  
WKLS



Fahnenstock  
WKFL



Lessmann  
WKLL



Knife  
OK5A



Gritman  
OKG



Flexible  
PS6



Spatula  
PS3



Spatula  
PS4

HEIDEMANN SPATULA

CEMENT SPATULA



PFIHS1 #1 #6  
PFIHS2 #1 #6  
PFIHS3 #1 #6



CS6 #3  
CS5 #3  
CS4 #3 #4  
CS3 #3

GINGIVAL CORD PACKER



Serrated GCP113 #3 #4  
Non-Serrated GCP113N #3 #4  
Serrated GCPAS113 #3 #4  
Serrated GCPR7 #3 #4  
Non-Serrated GCPR7N #3 #4

IMPRESSION TRAYS DENTULOUS PERFORATED



DENTULOUS NON-PERFORATED





IMPRESSION TRAYS EDENTULOUS PERFORATED

UPPER



ITRLEDPU1

ITRLEDPU2

ITRLEDPU3

ITRLEDPU4

LOWER



ITRLEDPL1

ITRLEDPL2

ITRLEDPL3

ITRLEDPL4

EDENTULOUS NON-PERFORATED

UPPER



ITRLEDNPU1

ITRLEDNPU2

ITRLEDNPU3

ITRLEDNPU4

LOWER



ITRLEDNPL1

ITRLEDNPL2

ITRLEDNPL3

ITRLEDNPL4

IMPRESSION TRAY DENTULOUS PERFORATED PEEDO



IMPRESSION TRAYS HALF



ANTERIOR AND POSTERIOR



ITRLDPUA

ITRLDPLA

ADJUSTABLE



ITRLDPA

ITRLDNPA

CALIPER IWANSON



Metal  
CLP1

Wax  
CLP2

ARTICULATING PAPER HOLDER



Miller Straight  
APF2E

Miller Curved  
APF2C

CROWN REMOVER

TURKEY PATTERN



MANUAL PATTERN





CROWN REMOVER

AUTOMATIC



CRS

CRP

Crown Remover Tips with CRS & CRP



CR1

CR2

CR3

CR4

CR5

CR6

CR7  
(For CRS)

CR8  
(For CRP)

CR9

CR10

**Basic Kit Includes :**

- Crown and Bridge Remover # CRS/CRP
- Single crown remover tip # CR1
- Crown remover tip Bifid # CR2
- Crown & bridge removing tip Sickle # CR3
- Single crown remover Contra angle tip # CR4
- Single crown remover tip Inverted curved # CR5
- Single Crown remover tip straight # CR6
- Bridge remover Flexi-wire set # CR7 (For CRS)
- Bridge remover tip with lock # CR8 (For CRP)
- Tips adjustment Key # CR9
- Open-end adjustment wrench # CR10

**Terms of Use**

Removal of :

- Permanently cemented crowns and/or bridges
- Temporarily cemented provisional or final restorations
- Maryland bridges
- Orthodontic brackets/bands
- Removable appliances "stuck" at attachments
- Multi-unit restoration "stuck" on working lab model

**Precautions :**

1. To remove crowns placed over pins or screw type posts, use lowest "strike intensity" to prevent fracture.
2. This instrument could cause damage to the margins of soft gold crowns. If re-use is anticipated, cut a small groove parallel to the lingual margin.
3. Highest strike intensity could cause root fracture; be careful in long span restorations.
4. Only low intensity strike force should be used with teeth having high mobility and little bone support.
5. For highest strike intensity settings, let the individual situation and common sense prevail.

ARTICULATOR



A5222

CROWN HOLDING FORCEP



CHFRS  
12.5cm

CHFRC  
12.5cm

PUTTY CUTTER



Right  
Handress  
PCR

Left  
Handress  
PCL

PLASTER SHEARS



PCS  
21cm

CROWN OPENER



CR4566

TEMPORARY  
CROWN REMOVING PLIER



CR4562



**Elektro**  
GENESIS


# Orthodontic

Pliers

Cutters

Black Line Cutters

Accessories





Orthodontic

ORTHODONTICS PLIERS

Every instrument we create becomes, in essence, a critical extension of each practitioner's skill. To achieve this, we handcraft each instrument individually, assuring instruments of the highest quality, durability and performance available to the profession.



**MATT FINISH**  
Dramatically reduces light reflection into the clinician's eye because the surface absorbs light versus high polish which reflects light.

**ORBIT JOINT**

- Consistent smooth functionality — glides smoothly with no wobble.
- True alignment of cutter tips — rattling at the joint means tips do not come together properly.
- Single part construction — no shims to wear out, seize up or corrode.

**STAINLESS STEEL**

- A premium blend of surgical stainless steel which provides corrosion resistance, optimal strength and long product life.
- 100% stainless steel from tip to tip — no dissimilar metals reaction.
- No chrome plating to scratch off, exposing lower quality metals that rust.

**LASER ENGRAVING OF APPLICABLE WIRE SIZES**  
All GDC Ortho cutters and pliers are laser engraved with the applicable wire minimums and maximums.



Weingart

How



Short  
3000/31



Slim  
3000/32



2mm  
Straight  
3000/84



2mm  
Curved  
3000/72

Also available in Wide Tip  
4mm Straight # 3000/85  
4mm Curved # 3000/73

Bracket Remover

Posterior Band Remover



Curved  
3000/82



Straight  
3000/83



Short  
3000/49



Long  
3000/50

UTILITY PLIERS



Anterior Band Remover  
3000/48



Jhonson Contouring  
3000/59



Ball & Socket  
3000/4



Band / Crown Crimping  
3000/225



Canine Contouring  
3000/310



Separating Module  
3000/391



Separating Module Angular  
3000/392

WIRE BENDING PLIERS

Bird Beak



De La Rosa

Ribbon/Tweed Arch



WIRE BENDING PLIERS



**Torquing Plier Male/Female**  
Suggested Pair

**Torquing Plier with Key**  
Suggested Pair





Band Forming



Right  
3000/42

Left  
3000/41

Band Pinching



Left  
3000/5

Right  
3000/6

LABORATORY PLIERS



Universal  
3000/55

Adam's  
3000/33

Optical  
3000/53

Ligature Tying  
3000/210

PLIERS WITH TUNGSTEN CARBIDE INSERTS



Weingart TC



Ribbon/Tweed Arch



Bird Beak



PLIERS WITH TUNGSTEN CARBIDE INSERTS



Orthodontic

LIGATURE CUTTER TUNGSTEN CARBIDE INSERTS



TUNGSTEN CARBIDE INSERTS



Micro Ligature Cutter  
3000/110



Ligature Cutter  
3000/708



Ligature Cutter  
3000/106



Ligature Cutter  
3000/107



Lingual Ligature Cutter 45°  
3000/111



Ligature Cutter 15°  
3000/109



DISTAL END CUTTER TUNGSTEN CARBIDE INSERTS



Long Handle  
3000/67



Short Handle  
3000/67S



With Hole Short Handle  
3000/67H



Flush Cut Long Handle  
3000/68



Micro Angular  
3000/100



Micro  
3000/101

HEAVY WIRE CUTTER



Hard Wire Cutter  
3000/7



Hard Wire Cutter  
3000/8



Hard Wire Cutter  
3000/9

Orthodontic

BLACK LINE CUTTER

BLACK LINE

- A performance engineered coating.

REDUCED LIGHT REFLECTION

- The black surface not only prevents reflection of the OT light, but also guarantees optimal contrast in the environment during treatment

TC Material

- Tungsten Carbide material insert in the tip to cut the wire easily.

Distal End Cutter  
3200/67B

RIVET JOINT

- Consistent smooth functionality — glides smoothly with no wobble.
- True alignment of cutter tips — rattling at the joint means tips do not come together properly.
- Single part construction — no shims to wear out, seize up or corrode.

Laser Engraving of Applicable Wire Size

- All GDC Ortho cutters & pliers are laser engraved with the applicable wire minimums and maximum

arch wire upto .021"x.025" (.55mm)



Pin Ligature Cutter  
3200/106



Pin Ligature Cutter 15°  
3200/109



Pin Ligature Cutter  
3200/708



Hard Wire Cutter  
3200/7

MATHIEU PLIERS



Smaha  
NHMS

Also available in mathieu smaha  
with ring # NHMSR



Smaha Curved  
NHMC



Fine  
NHMF



With Ring  
NHMR

BRACKET HOLDING & PLACING TWEEZER



1006



1008



1132

BUCCAL  
TUBE HOLDER



1134



Boley Gauge  
4982



Band Pusher  
4870



Wire Twister TC 15.0cm  
WTC



Wire Twister WT  
WT

LIGATURE INSTRUMENTS



Ligature Tucker  
672/2  
#3 #4



Ligature Directional Instrument  
672/1  
#3



Utility Instrument  
673  
#3



Elastic Applicator  
674  
#3



HEIGHT GAUGE WICK TYPE



018  
SS 4876/1



022  
SS 4876/2



Adjustable  
4875



Band Seater  
4871



Specifically designed for posterior positioning of any direct bond attachment excellent view afforded by the small head. Serrated tips give you a firm grip. The hard anodized aluminium allows you to work in comfort

Light Weight Positioner (Aluminium)  
2817



Dontrix Gauge  
2816



Boon's Gauge SS  
2812

For different heights for bracket positioning. 3.5mm, 4mm, 4.5mm & 5mm. Star points sit inside the brackets face and measure to the incisal edge.



Boon's Gauge with Lid  
2811

PLIER STAND  
PS



Distal Bending Instrument  
These instruments for the distal bending of the arches



2815





**Elektro**  
GENESIS  
**Instruments**  
MANAGEMENT

Stainless steel trays with locking instruments in silicon, resistant to corrosion.

GDC trays are simple to use and improve the protection and the storage of your valuable instruments.

Maximize the efficiency and performance of washing and sterilization.




CASSETTES FOR INSTRUMENTS

CASSETTES



For 4 Instruments  
IMSCA4E



For 5 Instruments  
IMSCA5E  
IMSCA5PPL



For 6 Instruments  
IMSCA6E



For 7 Instruments  
IMSCA7E



For 10 Instruments  
IMSCA10E  
IMSCA10PPL



For 6 Instruments Utility  
IMSCA6U



CASSETTES FOR INSTRUMENTS



For Extraction Forceps  
IMSCDEF



For Luxatip  
IMSCDL



For 20 Instruments  
IMSCA20E  
IMSCA20PPL




For Surgical Instruments  
IMSCDS



For Root Elevator  
IMSCDRE

INSTRUMENTS KIT

instruments kit




Extraction Forceps #Premium  
EFPWC12



Root Elevators #Premium  
REPWC9




Surgical Instruments  
SIWC10



Perio Surgical  
PSIWC12



Endo Surgical  
ESIWC13



Conservative Instruments #Premium  
CIWC19  
#4



Basic Implant Surgery Kit  
IMPBSWC19



Micro Surgery Kit  
MSIWC6



Luxatip New  
LWC7



Luxatip  
LGWC7



INSTRUMENTS KIT

Composite Instruments Titanium Coated



PMTWC4  
#4



Blue  
CIBWC6  
#4 #6



Gold  
CIGWC6  
#4 #6

Composite Instruments Titanium Coated



Filling Instruments  
FIWC10  
#1 #3 #4



Posterior  
CIPWC5  
#6



Anterior  
CIAWC5  
#6

Gracey Curettes



GCWC9  
#6



Colour Coded  
GCWC7  
#7



Colour Coded Mini Five  
GCMFWC7  
#7





Osteotomes Flat Angular  
OSTMPUAWC6

Osteotomes Flat Straight  
OSTMPUWC6

Also available in concave and convex with cassettes


Osteotomes Concave Angulated  
OSTMSHAWC6

Osteotomes Convex Angulated  
OSTMSPAWC5



Impaction Kit  
IMPACTION

Cleft & Palate Instrument Kit  
CLEFT



Trauma Kit  
TRAUMA

INSTRUMENTS KIT



Extraction Forceps #Standard  
EFSP12



Extraction Forceps #Premium  
EFPP12  
EFEP12 (Ergonomic forceps)



Extraction Forceps Peedo #Standard  
EFSP7



Extraction Forceps Peedo #Premium  
EFPP7



Root Elevator #Standard  
RESP9



Root Elevator #Premium  
REPP9



Root Elevator with Warwick James  
RESP10



Surgical Instruments  
SISP10

OFFER PACKAGE STANDARD  
OPS



OFFER PACKAGE PREMIUM  
OPP





INSTRUMENTS KIT



Composite Inst. Gold Flexi Thin  
CIGFTP6  
#6




Composite Inst. Gold  
CIGP6 #6  
CIBP6 #6  
(Also available in composite  
instruments titanium blue)



Composite Inst. Anterior Blue  
CIAP5  
#6



Composite Inst. Posterior Blue  
CIPP5  
#6



Conservative Instruments  
CIP19  
#1 #3 #4



Gracey Curette Premium  
GCPP9  
#6



Perio Scaler  
PSP8  
#3 #4 #6



Supra Gingival Scaler  
SGSP8





Sub Gingival Scaler  
SGSP7



Sinus Lift Instruments  
IMPSLP8  
#6



PKT Carver Colour Coded  
PKTCCP5



Impression Trays Dentulous Perforated  
ITRLDPP8



Impression Trays Edentulous Perforated  
ITRLEDPP8



Impression Trays Dentulous Non-Perforated  
ITRLDNPP8



Impression Trays Edentulous Non-Perforated  
ITRLEDNPP8

INSTRUMENTS KIT



PMT  
PMTS4  
#3



PMT Premium  
PMTP4  
#4



Filling Instruments  
FIS12  
#3 #4



RCT Pluggers  
RCPS6  
#3



RCT Spreaders  
RCSS6  
#3



Periotomes #5mm  
PTS3



Manipal Scalers  
MS3  
#3



Dental Dam Kit  
DDK (ADULT)  
DDKP (PEEDO)



Gracey Curettes  
GCS9  
#3



PKT Carvers  
PKTS5  
#3

INSTRUMENTS KIT



Plier Stand  
PS



Extraction Forcep Stand  
EFS



Endo Box  
AEB



Endo Box Premium  
AEBP



Cotton Holder  
ACH



Spirit Lamp  
ASL



Kidney Tray Autoclavable  
KT



Endo Box Round  
AEBR



Mouth Props with Chain  
MP3

X-RAY FILM HANGERS



AFH1



AFH6



AFH10



AFH14

NAPKIN HOLDERS



ANH2



ANH1

NOVELTIES



KCEF



KCRE



KCMMH



KCC



KCIMPT



KCSYR



ITEM NAME	HANDLE	PAGE #	ITEM NAME	HANDLE	PAGE #
<b>A</b>					
Allis Forcep Large Straight 15.0cm 5x6 # TFS		27	Bone Plating S S 4 Holes L-Shape Right 1.5 # P4HLR1.5SS		44
Allis Forcep Small Curved 12.5cm 4x5 # TFB		27	Bone Plating S S 4 Holes L-Shape Right 2.0 # P4HLR2.0SS		44
Amalgam Bone/Wall # WA		68	Bone Plating S S 4 Holes Orbital Type 1.5 # P4HO1.5SS		44
Amalgam Carrier D/E Jumbo/Regular # ACS203		68	Bone Plating S S 4 Holes Orbital Type 2.0 # P4HO2.0SS		44
Amalgam Carrier D/E Regular/Mini # ACS201		68	Bone Plating S S Drill Bite For Micromotor 1.5 # DSS1.5		44
Amalgam Carrier D/E Large/Regular # ACS202		68	Bone Plating S S Drill Bite For Micromotor 2.0 # DSS2.0		44
Amalgam Carrier S/E # AC504		68	Bone Plating S S Drill Bite For Micromotor 2.5 # DSS2.5		44
Art Instruments # Art1	#6	74	Bone Plating S S Screw 1.5 # S1.5SS		44
Art Instruments # Art2	#6	74	Bone Plating S S Screw 2.0 # S2.0SS		44
Art Instruments # CPS3/54	#6	74	Bone Plating S S Screw 2.5 # S2.5SS		44
Articulator Anatomical # AS222		99	Bone Plating Self Hoding Screw Driver 1.5 # SSH1.5		44
Articulating Paper Holder Straight # APF2E		96	Bone Plating Self Hoding Screw Driver 2.0 # SSH2.0		44
Articulating Paper Holder Curved # APF2C		96	Bone Plating Self Hoding Screw Driver 2.5 # SSH2.5		44
Aspirator Adaptor With Silicon Tubing # PTAA		33	Bone Plating Wire Twister # WTL		44
<b>B</b>			Bone Punch 4mm # BP6333		55
Band Pusher # 4870		113	Bone Punch 7mm # BP6373		55
Band Sealer # 4871		114	Bone Ronger Blumenthal # RBL		21
Basic Bone Plating Instruments Kit # Trauma	44, 122		Bone Ronger Blumenthal 45 # RBL45		21
Basic Implant Surgry Kit With Cassette # IMPBWC19		120	Bone Ronger Blumenthal 90 # RBL90		21
Boley Gauge # 4982		113	Bone Ronger Double Action # RBOE		21
Bone Caliper 15.5cm # BC15		53	Bone Ronger Friedman # RF		21
Bone Caliper 15.5cm # BC35		53	Bone Ronger Mead # R1A		21
Bone Chisel Chandler 5mm # CC7		33	Bone Ronger Micro Friedman 30 # RMFM30		21
Bone Chisel Chandler 7mm # CC9		33	Bone Ronger Micro Friedman 45 # RMFM45		21
Bone Chisel Gardner 3mm # CG1		33	Bone Ronger Micro Friedman 90 # RMFM90		21
Bone Chisel Gardner 4mm # CG2		33	Bone Ronger Micro Friedman With La Grange Bending 30 # RMFML		21
Bone Collector Filter # BC1331/1F		53	Bone Ronger Mini Friedman # RMF		21
Bone Collector Filter # BC1331/2F		53	Bone Scraper Blade Curved # IMPBSB2		52
Bone Collector With Filter 12.0cm # BC1331/2		53	Bone Scraper Blade Straight # IMPBSB1		52
Bone Collector With Filter 7.0cm # BC1331/1		53	Bone Scraper Cvd # IMPBC		52
Bone Compressor & Carrier # BCC14616	#6	54	Bone Scraper Str # IMPBS		52
Bone Cracker For Large Bone Pieces # IMPBMR		52	Bong Holding Forcep # BHP		44
Bone Crusher # IMPBMS		52	Boots Gauge With Lid # 2811		114
Bone Curette Hemingway 876 # CH876	#3,#4	86	Boots Gauge Without Lid # 2812		114
Bone Curette Hemingway 877 # CH877	#3,#4	86	Bracket Holding Tweezer Fine Tip # 1008		112
Bone Curette Hemingway 878 # CH878	#3,#4	86	Bracket Holding Tweezer Square Tip # 1006		112
Bone Curette Hemingway 879 # CH879	#3,#4	86	Bracket Placing Tweezer # 1132		112
Bone Cutting Forcep Single Action # RBC		44	Buccal Tube Holder # 1134		112
Bone File Howard # BF12	#3	33	Burnisher # BB21	#3,#4	66
Bone File Miller # BF64	#3	33	Burnisher # BB21B	#3,#4	66
Bone File Miller Large # BF52	#3	33	Burnisher # BBL3	#3	66
Bone File Miller Micro Straight # BFMS		43	Burnisher Ball # BB1B	#3,#4	66
Bone File Miller Small # BF52S	#3	33	Burnisher Ball # BB22/23	#3,#4	66
Bone Graft Carrier # BC14615	#6	54	Burnisher Cone # PKT3	#3,#4	66
Bone Graft Condenser # BC14614	#6	54	Burnisher Egg # BB27/29	#3,#4	66
Bone Injector With Cone 2.5mm/15.0cm # BIC4882		55	Burnisher T # TB2	#3,#4	66
Bone Injector With Cone 3.5mm/15.0cm # BIC4887		55	<b>C</b>		
Bone Injector With Cone 4.4mm/15.0cm # BIC4883		55	Caliper Iwanson Metal # CLP1		96
Bone Mill Forcep # IMPBM		52	Caliper Iwanson Wax # CLP2		96
Bone Modelling Lever # ML		44	Carver Cleoid-Discoid # CD 4/5	#3,#4	66
Bone Plating Ordinary Screw Driver 1.5 # OSD1.5		44	Carver Cleoid-Discoid # CD99/92	#3	66
Bone Plating Ordinary Screw Driver 2.0 # OSD2.0		44	Carver Diamond Shape Large # CVFR2/3	#3	66
Bone Plating Ordinary Screw Driver 2.5 # OSD2.5		44	Carver Diamond Shape Small # CVFR2/3H	#3,#4	66
Bone Plating Plier # BPP		44	Carver Hollen Back # CVHL3S	#3,#4	66
Bone Plating S S 16 Holes 1.5 # P16H1.5SS		44	Carver Interproximal Carver off Angle # CVPCOA	#3	66
Bone Plating S S 16 Holes 2.0 # P16H2.0SS		44	Carver Wall's # CVHL3W	#3	66
Bone Plating S S 16 Holes 2.5 # P16H2.5SS		44	Carver Ward's 1S # CVWR1S	#3,#4	66
Bone Plating S S 2 Holes 1.5 # P2H1.5SS		44	Castroviejo Caliper Large Curved 18.0cm # CL40C		53
Bone Plating S S 2 Holes 2.0 # P2H2.0SS		44	Castroviejo Caliper Large Straight 18.0cm # CLC40L		53
Bone Plating S S 2 Holes 2.5 # P2H2.5SS		44	Castroviejo Caliper Small Curved 8.5cm # CLC20C		53
Bone Plating S S 4 Holes 1.5 # P4H1.5SS		44	Castroviejo Caliper Small Straight 8.5cm # CLC20L		53
Bone Plating S S 4 Holes 2.0 # P4H2.0SS		44	Cement Spatula D/E # CS3	#3	92
Bone Plating S S 4 Holes 2.5 # P4H2.5SS		44	Cement Spatula D/E # CS4	#3,#4	92
Bone Plating S S 4 Holes L-Shape Left 1.5 # P4HL1.5SS		44	Cement Spatula S/E # CS5	#3	92
Bone Plating S S 4 Holes L-Shape Left 2.0 # P4HL2.0SS		44	Cement Spatula S/E # CS6	#3	92

INDEX

ITEM NAME	HANDLE	PAGE #
Chisel Dreyb Perry # SOPS/6K	#6	80
Chisel Wedelstaedt # CP1/2	#3	88
Chisel Wedelstaedt # CP3/4	#3	88
Chisel Wedelstaedt # CP5/6	#3	88
Cleft Lip & Palate Instruments Kit # CLEFT		46,122
Cleft Palate Dissector Left # CPDL		45
Cleft Palate Dissector Right # CPDR		45
Cleft Palate Down Hook # CPHD		46
Cleft Palate Elevator Cronin # CCP		46
Cleft Palate Elevator Right & Left # CPE1		45
Cleft Palate Elevator Side & Flate # CPE2		45
Cleft Palate Michel Trimmer # TM		46
Cleft Palate Resporator McIndoe # MCPR		45
Cleft Palate Respiratory Barsky # CPRB		46
Cleft Palate Respiratory Down # CPRD		45
Cleft Palate Respiratory French Type Left # CPRFL		46
Cleft Palate Respiratory French Type Right # CPRFR		46
Cleft Palate Respiratory Up # CPRU		45
Cleft Palate Triangle # CPKT		46
Columbia Curette 2R/2L # SC2R/2L	#3,#6,#7	78
Columbia Curette 4R/4L # SC4R/4L	#3,#6,#7	78
Composite Instrument Titanium Goldstein Flex thin # CIGT1	#6	62
Composite Instrument Titanium Goldstein Flex thin # CIGT2	#6	62
Composite Instrument Titanium Goldstein Flex thin # CIGT3	#6	62
Composite Instrument Titanium Goldstein Flex thin # CIGT4	#6	62
Composite Instrument Titanium Goldstein Flex thin # CIGTMM#1	#6	62
Composite Instrument Titanium Goldstein Flex thin # CIGTMM#3	#6	62
Composite Instruments Anterior Blue in Pouch # CIAF5	#6	125
Composite Instruments Anterior Blue With Cassette # CIAWC5	#6	121
Composite Instruments Blue With Cassette # CIBWC6	#6	121
Composite Instruments Gold Flex Thin In Pouch # CIGFTP6	#6	125
Composite Instruments Gold in Pouch # CIGP6	#6	125
Composite Instruments Gold With Cassette # CIGWC6	#6	121
Composite Instruments Plasma Black Silicon # CIPB1	#8	59
Composite Instruments Plasma Black Silicon # CIPB2	#8	59
Composite Instruments Plasma Black Silicon # CIPB3	#8	59
Composite Instruments Plasma Black Silicon # CIPB4	#8	59
Composite Instruments Plasma Black Silicon # CIPB5	#8	59
Composite Instruments Plasma Black Silicon # CIPB6	#8	59
Composite Instruments Posterior Blue in Pouch # CIPP5	#6	125
Composite Instruments Posterior Blue With Cassette # CIPWC5	#6	121
Composite Instruments Titanium Blue # CIB1	#6	60
Composite Instruments Titanium Blue # CIB2	#6	60
Composite Instruments Titanium Blue # CIB3	#6	60
Composite Instruments Titanium Blue # CIB4	#6	60
Composite Instruments Titanium Blue # CIB5	#6	60
Composite Instruments Titanium Blue # CIB6	#6	60
Composite Instruments Titanium Blue Anterior # TNCIGFT3	#6	61
Composite Instruments Titanium Blue Anterior # TNCIGFTM#1	#6	61
Composite Instruments Titanium Blue Anterior # TNCIPCL	#6	61
Composite Instruments Titanium Blue Anterior # TNCIPCM	#6	61
Composite Instruments Titanium Blue Anterior # TNCIPCS	#6	61
Composite Instruments Titanium Blue Posterior # TNCRML	#6	61
Composite Instruments Titanium Blue Posterior # TNCRSM	#6	61
Composite Instruments Titanium Blue Posterior # TNPLGSA	#6	61
Composite Instruments Titanium Blue Posterior # TNPLGH3	#6	61
Composite Instruments Titanium Blue Posterior # TNPLGOT	#6	61
Composite Instruments Titanium Blue Silicon # CIB1	#8	59
Composite Instruments Titanium Blue Silicon # CIB2	#8	59
Composite Instruments Titanium Blue Silicon # CIB3	#8	59
Composite Instruments Titanium Blue Silicon # CIB4	#8	59
Composite Instruments Titanium Blue Silicon # CIB5	#8	59
Composite Instruments Titanium Blue Silicon # CIB6	#8	59
Composite Instruments Titanium Gold # CIG1	#6	60
Composite Instruments Titanium Gold # CIG2	#6	60
Composite Instruments Titanium Gold # CIG3	#6	60

ITEM NAME	HANDLE	PAGE #
Composite Instruments Titanium Gold # CIG4	#6	60
Composite Instruments Titanium Gold # CIG5	#6	60
Composite Instruments Titanium Gold # CIG6	#6	60
Cone Socket Mirror Front Surface on Both Side # MMT5DSCS (Single)		5
Cone Socket Mirror Front Surface on Both Side # MMT5DSCS6 (Box of 6pcs)		5
Cone Socket Mirror Front Surface on Both Side With Handle # MMH5DSCS	#4	5
Conservative Instruments In Pouch # CIP19	#3	125
Conservative Instruments Set With Cassette # CINC19	#4	120
Corn Suture Plier 15.0cm # SP20		85
Corn Suture Plier 15.0cm # SPM20		85
Corn Suture Plier 18.0cm # SP20SLEM		85
Cotton Holder # ACH		128
Crown Holding Forcep With Ring Curved 12.5cm # CHFR		99
Crown Holding Forcep With Ring Straight 12.5cm # CHFRS		99
Crown Opener # CR4566		99
Crown Remover Fully Automatic Premium # CRP		98
Crown Remover Fully Automatic Standard # CRS		98
Crown Remover Manual Pattern # CRMP		97
Crown Remover Manual Pattern Tip # CRMP1		97
Crown Remover Manual Pattern Tip # CRMP2		97
Crown Remover Manual Pattern Tip # CRMP3		97
Crown Remover Tip # CR1		98
Crown Remover Tip # CR10		98
Crown Remover Tip # CR2		98
Crown Remover Tip # CR3		98
Crown Remover Tip # CR4		98
Crown Remover Tip # CR5		98
Crown Remover Tip # CR6		98
Crown Remover Tip # CR7 (for Standard Crown Remover Tip)		98
Crown Remover Tip # CR8 (for Premium Crown Remover Tip)		98
Crown Remover Tip # CR9		98
Crown Remover Turkey Pattern # CRT		97
Crown Remover Turkey Pattern Tip # CRT1		97
Crown Remover Turkey Pattern Tip # CRT2		97
Crown Remover Turkey Pattern Tip # CRT3		97
Crown Remover Turkey Pattern Tip # CRT4		97
Crown Remover Turkey Pattern Tip # CRT5		97
Crown Remover Turkey Pattern Tip # CRT6		97
Crown Remover Turkey Pattern Tip # CR8		97
Crown Remover Turkey Pattern Tip # CR9		97
Cutter Distal End Black Line TC # 3200/67B		111
Cutter Distal End Flush Cut TC # 3000/68		110
Cutter Distal End Micro Angular TC # 3000/100		110
Cutter Distal End Micro TC # 3000/101		110
Cutter Distal End Short Handle TC # 3000/67S		110
Cutter Distal End TC # 3000/67		110
Cutter Distal End With Hole TC # 3000/67H		110
Cutter Hard Wire Black Line TC # 3200/7		111
Cutter Hard Wire TC # 3000/7		110
Cutter Hard Wire TC # 3000/8		110
Cutter Hard Wire TC # 3000/9		110
Cutter Ligature 15" TC # 3000/109		109
Cutter Ligature Black Line 15" TC # 3200/109		111
Cutter Ligature Black Line TC # 3200/106		111
Cutter Ligature Long Handle Black Line TC # 3200/708		111
Cutter Ligature Long Handle TC # 3000/708		109
Cutter Ligature TC # 3000/106		109
Cutter Ligature TC # 3000/107		109
Cutter Micro Ligature TC # 3000/110		109
Cutter Pin Ligature Lingual 45° TC # 3000/111		109

D		
Dental Dam Kit Adult # DDK		73,127
Dental Dam Kit Peedo # DDKP		73,127
Dental Syringe Aspirating 1.8 # SYRA1.8		32
Dental Syringe Aspirating 1.8 # SYRF1.8		32



ITEM NAME	HANDLE	PAGE #	ITEM NAME	PAGE #
Dental Syringe Aspirating 2.2 # SYRA2.2		32	Extraction Forcep American Pattern Peedo Universal Lower 151S # F151S	11
Dental Syringe Intraalveolar # SYR1		32	Extraction Forcep American Pattern Peedo Universal Upper 150S # F150S	11
Dental Syringe Not Aspirating 1.8 # SYR1.8		32	Extraction Forcep American Pattern Tomes Upper & Lower Fragments & Roots 89 # F89	11
Dental Syringe Not Aspirating 2.2 # SYR2.2		32	Extraction Forcep Atraumatic Lower Premolars 13 # AFX13	13
Dental Syringe Self Aspirating With Colour Coded Ring # SYRSA1.8		32	Extraction Forcep Atraumatic Lower Roots 74N # AFX74N	13
Dental Syringe With Colour Coded Ring 1.8 # SYRACC1.8R		32	Extraction Forcep Atraumatic Lower Third Molars 79 # AFX79	13
Dingman Mouth Gag With 4 Blades # DMG		46	Extraction Forcep Atraumatic Upper Centrals & Canines 1 # AFX1	13
Distal Wire Bending Instruments 55.0mm/75.0mm # 2815		115	Extraction Forcep Atraumatic Upper Premolars 7 # AFX7	13
Donrix Gauge # 2816		114	Extraction Forcep Atraumatic Upper Roots 51A # AFX51A	13
Dressing Forcep Foerster Straight 20.0cm # SDF		29	Extraction Forcep Black Line Lower Molars 73 # FX73B	14
Dressing Forcep Gross Maier Curved 20.0cm # SDFC		29	Extraction Forcep Black Line Lower Molars 86 # FX86B	14
Dressing Forcep Gross Maier Straight 20cm # SDFS		29	Extraction Forcep Black Line Lower Premolars 13 # FX13B	14
Drill Kit # IMFDS8		53	Extraction Forcep Black Line Lower Roots 74N # FX74NB	14
Dycal Applicator D/E # PICH6	#3	67	Extraction Forcep Black Line Upper Anterior 1 # FX1B	14
Dycal Applicator D/E # PICH7	#3	67	Extraction Forcep Black Line Upper Molars Left 18 # FX18B	14
Dycal Applicator S/E # PICH8		67	Extraction Forcep Black Line Upper Molars Right 17 # FX17B	14
Dycal Applicator With Spatula D/E # SP60616	#3	67	Extraction Forcep Black Line Upper Molars 89 # FX89B (Right)	14
			Extraction Forcep Black Line Upper Molars 90 # FX90B (Left)	14
			Extraction Forcep Black Line Upper Premolars 7 # FX7B	14
			Extraction Forcep Black Line Upper Roots 51A # FX51AB	14
			Extraction Forcep Black Line Upper Third Molars 67A # FX67AB	14
			Extraction Forcep Ergonomic in Pouch # EFEP12	123
			Extraction Forcep Lower Anterior Roots 74 # FX74S (Standard)	10
			Extraction Forcep Lower Molars & Wisdom Left 22 1/2L # FX22 1/2LP (Premium)	13
			Extraction Forcep Lower Molars & Wisdom Left 22 1/2L # FX22 1/2LS (Standard)	13
			Extraction Forcep Lower Molars & Wisdom Right 22 1/2R # FX22 1/2RP (Premium)	13
			Extraction Forcep Lower Molars & Wisdom Right 22 1/2R # FX22 1/2RS (Standard)	13
			Extraction Forcep Lower Molars 22 # FX22P (Premium)	9
			Extraction Forcep Lower Molars 22 # FX22S (Standard)	9
			Extraction Forcep Lower Molars 73 # FX73E (Ergonomic)	10
			Extraction Forcep Lower Molars 73 # FX73P (Premium)	10
			Extraction Forcep Lower Molars 73 # FX73S (Standard)	10
			Extraction Forcep Lower Molars 86 # FX86E (Ergonomic)	10
			Extraction Forcep Lower Molars 86 # FX86P (Premium)	10
			Extraction Forcep Lower Molars 86 # FX86S (Standard)	10
			Extraction Forcep Lower Molars 86A # FX86AP (Premium)	10
			Extraction Forcep Lower Molars 86A # FX86AS (Standard)	10
			Extraction Forcep Lower Molars 86B # FX86BP (Premium)	10
			Extraction Forcep Lower Molars 86B # FX86BS (Standard)	10
			Extraction Forcep Lower Molars 86C # FX86CP (Premium)	10
			Extraction Forcep Lower Molars 86C # FX86CS (Standard)	10
			Extraction Forcep Lower Molars 87 # FX87S (Standard)	10
			Extraction Forcep Lower Premolars 13 # FX13E (Ergonomic)	9
			Extraction Forcep Lower Premolars 13 # FX13P (Premium)	9
			Extraction Forcep Lower Premolars 13 # FX13S (Standard)	9
			Extraction Forcep Lower Premolars 75 # FX75S (Standard)	10
			Extraction Forcep Lower Roots 31 # FX31S (Standard)	10
			Extraction Forcep Lower Roots 33 # FX33P (Premium)	10
			Extraction Forcep Lower Roots 33 # FX33S (Standard)	10
			Extraction Forcep Lower Roots 33A # FX33AP (Premium)	10
			Extraction Forcep Lower Roots 33A # FX33AS (Standard)	10
			Extraction Forcep Lower Roots 33L # FX33LP (Premium)	10
			Extraction Forcep Lower Roots 33L # FX33LS (Standard)	10
			Extraction Forcep Lower Roots 45 # FX45P (Premium)	10
			Extraction Forcep Lower Roots 45 # FX45S (Standard)	10
			Extraction Forcep Lower Roots 74N # FX74NE (Ergonomic)	10
			Extraction Forcep Lower Roots 74N # FX74NP (Premium)	10
			Extraction Forcep Lower Roots 74N # FX74NS (Standard)	10
			Extraction Forcep Lower Third Molars 79 # FX79E (Ergonomic)	10
			Extraction Forcep Lower Third Molars 79 # FX79P (Premium)	10
			Extraction Forcep Lower Third Molars 79 # FX79S (Standard)	10
			Extraction Forcep Peedo in Pouch # EFPSP7	123
			Extraction Forcep Peedo Lower Incisors 1 # FX1CP (Premium)	15
			Extraction Forcep Peedo Lower Incisors 1 # FX1CS (Standard)	15
			Extraction Forcep Peedo Lower Primary Molars 2 # FX2CP (Premium)	15
			Extraction Forcep Peedo Lower Primary Molars 2 # FX2CS (Standard)	15
			Extraction Forcep Peedo Lower Roots 3 # FX3CP (Premium)	15
<b>E</b>				
Elastic Placing Instrument # 674	#3	113		
Enamel Chisel D/E # CP7/8	#3	68		
Enamel Chisel S/E 1mm # CP7	#3,#4	68		
Enamel Chisel S/E 2mm # CP8	#3	68		
Enamel Chisel S/E 3mm # CP9	#3	68		
Enamel Chisel S/E 4mm # CP10	#3	68		
Enamel Hatchet # CP53/54	#3,#4	68		
Enamel Hoe S/E # CP24	#3,#4	68		
Endo Box Premium # AEBP		128		
Endo Box Round # AEBR		128		
Endo Box Standard # AEB		128		
Endo Explorer # EXDG16	#3,#6,#7	4		
Endo Surgical Instruments Set With Cassette # ESIWC13		119		
Excavator # EXC125/126	#3,#6	63,74		
Excavator # EXC127/128	#3,#6	63		
Excavator # EXC129/130	#3,#6	63		
Excavator # EXC131/132	#3,#6	63		
Excavator # EXC133/134	#3,#6	63,74		
Excavator # EXC14	#3,#6	63		
Excavator # EXC153/154	#3,#6	64		
Excavator # EXC17H	#3,#6	63		
Excavator # EXC17WH	#3,#6	63		
Excavator # EXC18H	#3,#6	63		
Excavator # EXC18WH	#3,#6	63		
Excavator # EXC19H	#3,#6	63		
Excavator # EXC19WH	#3,#6	63		
Excavator # EXC21G/211	#3,#6	64		
Excavator # EXC212/213	#3	64		
Excavator # EXC31L	#3,#6	64		
Excavator # EXC31W	#3,#6	64		
Excavator # EXC32L	#3,#6	64		
Excavator # EXC33L	#3,#6	64		
Excavator # EXC33L	#3,#6	64		
Excavator # EXC38/39	#6	64		
Excavator # EXC816	#3	64		
Excavator # EXC817	#3	64		
Excavator # EXC819	#3	64		
Excavator # EXC820	#3	64		
Excavator # EXCGL2	#6	64		
Extracting Forcep Albrecht Molars 2 # FXA2		13		
Extracting Forcep Albrecht Premolars 1 # FXA1		13		
Extraction Forcep American Pattern Cowhorn 1st & 2nd Upper Molars 23 # F23		11		
Extraction Forcep American Pattern Cryer Universal Lower 151 # F151		11		
Extraction Forcep American Pattern Cryer Universal Upper 150 # F150		11		
Extraction Forcep American Pattern Molars 1st & 2nd Upper Molars 53L # F53L (Left)		11		
Extraction Forcep American Pattern Molars 1st & 2nd Upper Molars 53R # F53R (Right)		11		
Extraction Forcep American Pattern Nevius 1st & 2nd Upper Molars 88L # F88L (Left)		11		
Extraction Forcep American Pattern Nevius 1st & 2nd Upper Molars 88R # F88R (Right)		11		

INDEX

ITEM NAME	PAGE #
Extraction Forcep Peedo Lower Roots 3 # FX3CS (Standard)	15
Extraction Forcep Peedo Premium in Pouch # EFPFP7	123
Extraction Forcep Peedo Upper Incisors & Canines 5 # FX5CP (Premium)	15
Extraction Forcep Peedo Upper Incisors & Canines 5 # FX5CS (Standard)	15
Extraction Forcep Peedo Upper Molars 6 # FX6CP (Premium)	15
Extraction Forcep Peedo Upper Molars 6 # FX6CS (Standard)	15
Extraction Forcep Peedo Upper Primary Molars 4 # FX4CP (Premium)	15
Extraction Forcep Peedo Upper Primary Molars 4 # FX4CS (Standard)	15
Extraction Forcep Peedo Upper Roots 7 # FX7CP (Premium)	15
Extraction Forcep Peedo Upper Roots 7 # FX7CS (Standard)	15
Extraction Forcep Premium in Pouch # EFPFP12	123
Extraction Forcep Secure Lower Root 845.00 # SFX845.00	12
Extraction Forcep Secure Lower Root 945.00 # SFX945.00	12
Extraction Forcep Secure Lower Root 959.01 # SFX959.01	12
Extraction Forcep Secure Lower Root 974.00 # SFX974.00	12
Extraction Forcep Secure Upper Root 149.11 # SFX149.11	12
Extraction Forcep Secure Upper Root 197.11 # SFX197.11	12
Extraction Forcep Secure Upper Root 849.00 # SFX849.00	12
Extraction Forcep Secure Upper Root 897.00 # SFX897.00	12
Extraction Forcep Secure Upper Root 944.00 # SFX944.00	12
Extraction Forcep Secure Upper Root 944.01 # SFX944.01	12
Extraction Forcep Secure Upper Root 951.00 # SFX951.00	12
Extraction Forcep Separating Lower Molars 56 # FX56S (Standard)	8
Extraction Forcep Separating Upper Molars 55 # FX55S (Standard)	8
Extraction Forcep Stand Acrylic # EFS	128
Extraction Forcep Standard in Pouch # EFPSP12	123
Extraction Forcep Universal For Lower Roots # FX223	12
Extraction Forcep Universal For Upper Roots # FX221	12
Extraction Forcep Upper Anterior 1 # FX1E (Ergonomic)	8
Extraction Forcep Upper Anterior 1 # FX1P (Premium)	8
Extraction Forcep Upper Anterior 1 # FX1S (Standard)	8
Extraction Forcep Upper Molars 18S # FX18AS (Standard)	8
Extraction Forcep Upper Molars Left 18 # FX18E (Ergonomic)	8
Extraction Forcep Upper Molars Left 18 # FX18P (Premium)	8
Extraction Forcep Upper Molars Left 18 # FX18S (Standard)	8
Extraction Forcep Upper Molars Right 17 # FX17E (Ergonomic)	8
Extraction Forcep Upper Molars Right 17 # FX17P (Premium)	8
Extraction Forcep Upper Molars Right 17 # FX17S (Standard)	8
Extraction Forcep Upper Molars Right 8B # FX8BE (Ergonomic)	9
Extraction Forcep Upper Molars Right 8B # FX8BP (Premium)	9
Extraction Forcep Upper Molars Right 8B # FX8BS (Standard)	9
Extraction Forcep Upper Molars Left 90 # FX90E (Ergonomic)	9
Extraction Forcep Upper Molars Left 90 # FX90P (Premium)	9
Extraction Forcep Upper Molars Left 90 # FX90S (Standard)	9
Extraction Forcep Upper Premolars 7 # FX7E (Ergonomic)	8
Extraction Forcep Upper Premolars 7 # FX7P (Premium)	8
Extraction Forcep Upper Premolars 7 # FX7S (Standard)	8
Extraction Forcep Upper Roots 29 # FX29P (Premium)	8
Extraction Forcep Upper Roots 29 # FX29S (Standard)	8
Extraction Forcep Upper Roots 30 # FX30P (Premium)	8
Extraction Forcep Upper Roots 30 # FX30S (Standard)	8
Extraction Forcep Upper Roots 41 # FX41P (Premium)	9
Extraction Forcep Upper Roots 41 # FX41S (Standard)	9
Extraction Forcep Upper Roots 44 # FX44P (Premium)	9
Extraction Forcep Upper Roots 44 # FX44S (Standard)	9
Extraction Forcep Upper Roots 51 # FX51S (Standard)	9
Extraction Forcep Upper Roots 51A # FX51AE (Ergonomic)	9
Extraction Forcep Upper Roots 51A # FX51AP (Premium)	9
Extraction Forcep Upper Roots 51A # FX51AS (Standard)	9
Extraction Forcep Upper Roots 76N # FX76NS (Standard)	9
Extraction Forcep Upper Roots Narrow 29N # FX29NP (Premium)	8
Extraction Forcep Upper Roots Narrow 29N # FX29NS (Standard)	8
Extraction Forcep Upper Third Molars 67A # FX67AE (Ergonomic)	9
Extraction Forcep Upper Third Molars 67A # FX67AP (Premium)	9
Extraction Forcep Upper Third Molars 67A # FX67AS (Standard)	9
Extraction Forcep With Cassette # EFPWC12	119

ITEM NAME	HANDLE	PAGE #
<b>F</b>		
Filing Instruments Set # FS12	#3	127
Filing Instruments Set With Cassette # FWC10	#4	121
Furcation Curette Buccal/Lingual # SQBL2	#6	80
Furcation Curette Mesial/Distal # SQMD2	#6	80
<b>G</b>		
Gingival Cord Packer # GCP113	#3,#4	92
Gingival Cord Packer # GCP113N	#3,#4	92
Gingival Cord Packer # GCPAS113	#3,#4	92
Gingival Cord Packer # GCPFR7	#3,#4	92
Gingival Cord Packer # GCPFR7N	#3,#4	92
Goldman Fox Pocket Marker Left # PMGF1		86
Goldman Fox Pocket Marker Right # PMGF2		86
Gracey Curette Colour Coded # SG1/2	#7	77
Gracey Curette Colour Coded # SG11/12	#7	77
Gracey Curette Colour Coded # SG13/14	#7	77
Gracey Curette Colour Coded # SG3/4	#7	77
Gracey Curette Colour Coded # SG5/6	#7	77
Gracey Curette Colour Coded # SG7/8	#7	77
Gracey Curette Colour Coded # SG9/10	#7	77
Gracey Curette Colour Coded Instruments Mini Five With Cassette # GCMFWC7	#7	121
Gracey Curette Colour Coded Instruments With Cassette # GCWC7	#7	121
Gracey Curette Colour Coded Mini Five # SAS1/2	#7	77
Gracey Curette Colour Coded Mini Five # SAS11/12	#7	77
Gracey Curette Colour Coded Mini Five # SAS13/14	#7	77
Gracey Curette Colour Coded Mini Five # SAS15/16	#7	77
Gracey Curette Colour Coded Mini Five # SAS3/4	#7	77
Gracey Curette Colour Coded Mini Five # SAS5/6	#7	77
Gracey Curette Colour Coded Mini Five # SAS7/8	#7	77
Gracey Curette Instruments With Cassette # GCWC9	#6	121
Gracey Curette Premium in Pouch # GCPFP9	#6	125
Gracey Curette Set # GCS9	#3	127
Gracey Curette Standard # SG1/2	#6	76
Gracey Curette Standard # SG11/12	#6	76
Gracey Curette Standard # SG13/14	#6	76
Gracey Curette Standard # SG15/16	#6	76
Gracey Curette Standard # SG17/18	#6	76
Gracey Curette Standard # SG3/4	#6	76
Gracey Curette Standard # SG5/6	#6	76
Gracey Curette Standard # SG7/8	#6	76
Gracey Curette Standard # SG9/10	#6	76
Gracey Curette Standard Mini Five # SAS1/2	#6	77
Gracey Curette Standard Mini Five # SAS11/12	#6	77
Gracey Curette Standard Mini Five # SAS13/14	#6	77
Gracey Curette Standard Mini Five # SAS15/16	#6	77
Gracey Curette Standard Mini Five # SAS3/4	#6	77
Gracey Curette Standard Mini Five # SAS5/6	#6	77
Gracey Curette Standard Mini Five # SAS7/8	#6	77
<b>H</b>		
Heidemann Spatula # PFRHS1	#1,#6	92
Heidemann Spatula # PFRHS2	#1,#6	92
Heidemann Spatula # PFRHS3	#1,#6	92
Height Gauge Wick S S Adjustable 018/022 # 4875		114
Height Gauge Wick S S 018 # SS 4876/1		114
Height Gauge Wick S S 022 # SS 4876/2		114
Hemostat Crile Curved 16.0cm # HRPC		27
Hemostat Crile Straight 16.0cm # HRPS		27
Hemostat Mosquito Forcep Curved 12.5cm # H3		27
Hemostat Mosquito Forcep Curved 1X2 12.5cm # H6		27
Hemostat Mosquito Forcep Straight 12.5cm # H4		27
Hemostat Mosquito Forcep Straight 1X2 12.5cm # H5		27
Hemostate Forcep Baby Curved 10.0cm # HHC		27
Hemostate Forcep Baby Straight 10.0cm # HHS		27
Hemostate Forcep Kelly Curved 14.0cm # H2		27
Hemostate Forcep Kelly Straight 14.0cm # H1		27



ITEM NAME	HANDLE	PAGE #
Hirschfeld File Anterior # FH3/10	#6	84
Hirschfeld File Buccal/Lingual # FH3/7	#6	84
Hirschfeld File Mesial/Distal # FH5/11	#6	84
<b>I</b>		
Impaction Kit With Box # IMPACTION		47,122
Impression Tray Dentulous Non Perforated Adjustable # ITRLDNPA		96
Impression Tray Dentulous Non Perforated in Pouch # ITRLDNPPB		126
Impression Tray Dentulous Non Perforated Lower # ITRLDNPL1		93
Impression Tray Dentulous Non Perforated Lower # ITRLDNPL2		93
Impression Tray Dentulous Non Perforated Lower # ITRLDNPL3		93
Impression Tray Dentulous Non Perforated Lower # ITRLDNPL4		93
Impression Tray Dentulous Non Perforated Upper # ITRLDNPU1		93
Impression Tray Dentulous Non Perforated Upper # ITRLDNPU2		93
Impression Tray Dentulous Non Perforated Upper # ITRLDNPU3		93
Impression Tray Dentulous Non Perforated Upper # ITRLDNPU4		93
Impression Tray Dentulous Perforated Adjustable # ITRLDPA		96
Impression Tray Dentulous Perforated Half Left # ITRLDPHL1		95
Impression Tray Dentulous Perforated Half Left # ITRLDPHL2		95
Impression Tray Dentulous Perforated Half Left # ITRLDPHL3		95
Impression Tray Dentulous Perforated Half Left # ITRLDPHL4		95
Impression Tray Dentulous Perforated Half Right # ITRLDPHR1		95
Impression Tray Dentulous Perforated Half Right # ITRLDPHR2		95
Impression Tray Dentulous Perforated Half Right # ITRLDPHR3		95
Impression Tray Dentulous Perforated Half Right # ITRLDPHR4		95
Impression Tray Dentulous Perforated in Pouch # ITRLDPPB		126
Impression Tray Dentulous Perforated Lower # ITRLDPL1		93
Impression Tray Dentulous Perforated Lower # ITRLDPL2		93
Impression Tray Dentulous Perforated Lower # ITRLDPL3		93
Impression Tray Dentulous Perforated Lower # ITRLDPL4		93
Impression Tray Dentulous Perforated Lower Anterior # ITRLDPLA		96
Impression Tray Dentulous Perforated Peedo Lower # ITRLDPLL1		95
Impression Tray Dentulous Perforated Peedo Lower # ITRLDPLL2		95
Impression Tray Dentulous Perforated Peedo Lower # ITRLDPLL3		95
Impression Tray Dentulous Perforated Peedo Lower # ITRLDPLL4		95
Impression Tray Dentulous Perforated Peedo Upper # ITRLDPLU1		95
Impression Tray Dentulous Perforated Peedo Upper # ITRLDPLU2		95
Impression Tray Dentulous Perforated Peedo Upper # ITRLDPLU3		95
Impression Tray Dentulous Perforated Peedo Upper # ITRLDPLU4		95
Impression Tray Dentulous Perforated Upper # ITRLDPU1		93
Impression Tray Dentulous Perforated Upper # ITRLDPU2		93
Impression Tray Dentulous Perforated Upper # ITRLDPU3		93
Impression Tray Dentulous Perforated Upper # ITRLDPU4		93
Impression Tray Dentulous Perforated Upper Anterior # ITRLDPUA		96
Impression Tray Edentulous Non Perforated in Pouch # ITRLEDNPPB		126
Impression Tray Edentulous Non Perforated Lower # ITRLEDNPL1		94
Impression Tray Edentulous Non Perforated Lower # ITRLEDNPL2		94
Impression Tray Edentulous Non Perforated Lower # ITRLEDNPL3		94
Impression Tray Edentulous Non Perforated Lower # ITRLEDNPL4		94
Impression Tray Edentulous Non Perforated Upper # ITRLEDNPU1		94
Impression Tray Edentulous Non Perforated Upper # ITRLEDNPU2		94
Impression Tray Edentulous Non Perforated Upper # ITRLEDNPU3		94
Impression Tray Edentulous Non Perforated Upper # ITRLEDNPU4		94
Impression Tray Edentulous Perforated in Pouch # ITRLEDPB		126
Impression Tray Edentulous Perforated Lower # ITRLEDPL1		94
Impression Tray Edentulous Perforated Lower # ITRLEDPL2		94
Impression Tray Edentulous Perforated Lower # ITRLEDPL3		94
Impression Tray Edentulous Perforated Lower # ITRLEDPL4		94
Impression Tray Edentulous Perforated Upper # ITRLEDPU1		94
Impression Tray Edentulous Perforated Upper # ITRLEDPU2		94
Impression Tray Edentulous Perforated Upper # ITRLEDPU3		94
Impression Tray Edentulous Perforated Upper # ITRLEDPU4		94
Instrument Management System For 10 Instrument # IMSCA10E		117
Instrument Management System For 10 Instrument # IMSCA10PPL		117
Instrument Management System For 20 Instrument # IMSCA20E		118
Instrument Management System For 20 Instrument # IMSCA20PPL		118
Instrument Management System For 4 Instrument # IMSCA4E		117

ITEM NAME	HANDLE	PAGE #
Instrument Management System For 5 Instrument # IMSCA5E		117
Instrument Management System For 5 Instrument # IMSCA5PPL		117
Instrument Management System For 6 Instrument # IMSCA6E		117
Instrument Management System For 7 Instrument # IMSCA7E		117
Instrument Management System For Extraction Forcep # IMSCDEF		118
Instrument Management System For Luxatip # IMSCDL		118
Instrument Management System For Root Elevator # IMSCDRE		118
Instrument Management System For Surgical Instruments # IMSCDS		118
Instrument Management System Utility For 6 Instrument # IMSCAGU		117
Instruments Measure Depth of Bony Tunnel # PCP1301		54

**K**

Key Chain Cutter # KCC		129
Key Chain Dental Cartridge Syringe # KCSYR		129
Key Chain Extraction Forcep # KCEF		129
Key Chain Impression Tray # KCMPT		129
Key Chain Mouth Mirror Handle # KCMWH		129
Key Chain Root Elevator # KDRE		129
Kidney Tray # KT		128
Knife Buck's # 3/4 # KB3/4	#6	85
Knife Buck's # 5/6 # KB5/6	#6	85
Knife Crane Kaplan # KCK3	#6	85
Knife Goldman Fox # KGF11	#6	85
Knife Goldman Fox # KGF8	#6	85
Knife Goldman Fox # KGF9	#6	85
Knife Kirkland # KK15/16	#6	85
Knife Orban's # KD1/2	#6	85
Knife Orban's # KD12KPO3A	#6	85

**L**

Larger Curette # SL1/2	#6	78
Larger Curette # SL3/4	#6	78
Larger Curette # SL5/6	#6	78
Ligature Instrument # 672/1	#3	113
Ligature Instrument Utility # 673	#3	113
Ligature Tucker # 672/2	#3, #4	113
Light Weight Positioner (Aluminium) # 2817		114
Lingual Tying Plier # 3000/210		106
Luxatip Contra Angle 3mm # L3CA		48
Luxatip Curved Blade 3mm # L3C		48
Luxatip Curved Blade 5mm # L5C		48
Luxatip Dual Edge Straight 3mm # DE3		48
Luxatip Dual Edge Straight 5mm # DE5		48
Luxatip Grey With Cassette # LGWC7		120
Luxatip Inverted Curved 3mm # L3IC		48
Luxatip New With Cassette # LWC7		120
Luxatip Straight Blade 2mm # L2S		48
Luxatip Straight Blade 3mm # L3S		48
Luxatip Straight Blade 5mm # L5S		48

**M**

Mallet Mead 1 19.0cm # MAL1		33,55
Mallet Mead 2 17.0cm # MAL2		33,55
Manipal Scaler Set # MS3		127
Margin Trimmer # MT26H	#6	67
Margin Trimmer # MT27H	#6	67
Margin Trimmer # MT28H	#3, #4, #6	67
Margin Trimmer # MT29H	#3, #4, #6	67
Mathiew Plier Fine # NHMF		112
Mathiew Plier Smaha Curved # NHMC		112
Mathiew Plier Smaha Straight # NHMS		112
Mathiew Plier With Ring # NHMR		112
Matrix Ring Holder # MRDF		73
Matrix Retainer Ivory # MRI		68

INDEX

ITEM NAME	HANDLE	PAGE #
Matrix Retainer Squivland Narrow # MR4972S		68
Matrix Retainer Squivland Wide # MR4972L		68
Matrix Retainer Toffelmire Universal # MRTU		68
Mc Call Curette # SM11/12	#6	78
Mc Call Curette # SM11/12A	#6	78
Mc Call Curette # SM11/12S	#6	78
Mc Call Curette # SM13/14	#6	78
Mc Call Curette # SM13/14S	#6	78
Mc Call Curette # SM17/18	#6	78
Mc Call Curette # SM19/20	#6	78
Membrane Placement # PMPIS	#6	54
Micro Castroviejo Needle Holder Curved 14.0cm # NHCVC		28
Micro Castroviejo Needle Holder Curved 14.0cm T C # NH5021		29
Micro Castroviejo Needle Holder Straight 14.0cm # NHCVS		28
Micro Castroviejo Needle Holder Straight 14.0cm T C # NH5020		29
Micro Castroviejo Needle Holder Straight 18.0cm T C # NHMS024R		29,57
Micro Castroviejo Scissor Curved 18.0cm # SPV		25, 57
Micro Mirror Large Rectangular 4x12mm # MM3	#6	5
Micro Mirror Large Round 4.5mm # MMR	#6	5
Micro Mirror Medium Rectangular 3x9mm # MM2	#6	5
Micro Mirror Small Rectangular 2x7mm # MM1	#6	5
Micro Mirror Small Round 2.5mm # MMRS	#6	5
Micro Scalpel Handle 12.0cm # 1-015		57
Micro Surgery Instruments with Cassette # MSINC6		120
Micro Tissue Forcep 1x2 Straight 18.0cm # TPSSSTMBH		31,57
Micro Tissue Forcep Curved 18.0cm # TPSLCOCMN		31
Micro Tissue Forcep Gerald 1x2 Curved 18.0cm # TPG4		31
Micro Tissue Forcep Gerald 1x2 Straight 18.0cm # TPG3		31
Micro Tissue Forcep Gerald Curved 18.0cm # TPG2		31
Micro Tissue Forcep Gerald Straight 18.0cm # TPG1		31
Micro Tissue Forcep Straight 18.0cm # TPSCOSM		31
Mirror Handle # MH1	#1	3
Mirror Handle # MH2	#2	3
Mirror Handle # MH3	#3	3
Mirror Handle # MH4	#4	3
Mirror Handle Signature Series #Grips # MHG1	#7	3
Mirror Handle Signature Series #Grips # MHG10	#7	3
Mirror Handle Signature Series #Grips # MHG2	#7	3
Mirror Handle Signature Series #Grips # MHG3	#7	3
Mirror Handle Signature Series #Grips # MHG4	#7	3
Mirror Handle Signature Series #Grips # MHG5	#7	3
Mirror Handle Signature Series #Grips # MHG6	#7	3
Mirror Handle Signature Series #Grips # MHG7	#7	3
Mirror Handle Signature Series #Grips # MHG8	#7	3
Mirror Handle Signature Series #Grips # MHG9	#7	3
Mirror Handle with Millimeter Scale # MHS		3
Moon's Probe # P13M	#3	34
Mouth Mirror Top Front Surface on Both Side # MMTDC/6 (Box of 6pcs)		5
Mouth Mirror Top Front Surface on Both Side # MMTDSS (Single)		5
Mouth Mirror Top Front Surface Rhodium 4 # MMTFS4 (Box of 12pcs)		5
Mouth Mirror Top Front Surface Rhodium 5 # MMTFS5 (Box of 12pcs)		5
Mouth Mirror Top Magnified 4 # MMTM4 (Box of 12pcs)		5
Mouth Mirror Top Magnified 5 # MMTM5 (Box of 12pcs)		5
Mouth Mirror Top Plain 3 # MMTTP3 (Box of 12pcs)		5
Mouth Mirror Top Plain 4 # MMTTP4 (Box of 12pcs)		5
Mouth Mirror Top Plain 5 # MMTTP5 (Box of 12pcs)		5
Mouth Prop With Chain # MP3		128
MTA Carrier Curved 0.8mm/13.0cm # MTA4869		67
MTA Carrier Curved 1.2mm/13.0cm # MTA4866		67
MTA Carrier Straight 1.2mm/13.0cm # MTA4865		67

**N**

Napkin Holder # ANH1		129
Napkin Holder # ANH2		129
Nasal Chisel with Guard Left # NCGL		42
Nasal Chisel with Guard Right # NCGR		42
Needle Holder Adson Straight 17.0cm # NHA		28

ITEM NAME	HANDLE	PAGE #
Needle Holder Crile Wood Straight 15.0cm T C # NH5038		29
Needle Holder Crile-Wood Straight 15.0cm # NHCW		28
Needle Holder Derf Curved 12.5cm # NHDC		28
Needle Holder Derf Straight 11.5cm T C # NH5032		29
Needle Holder Derf Straight 12.5cm # NHD		28
Needle Holder Halsey Straight 13.0cm T C # NH5037		29
Needle Holder Mayo Hegar Curved 16.0cm # NHMH		28
Needle Holder Mayo Hegar Straight 16.0cm # NHMH		28
Needle Holder Mayo Hegar Straight 16.0cm T C # NH5042		29
Needle Holder Mayo Hegar Straight 18.0cm # NHMH18		28
Needle Holder Olsen Hegar With Scissor 17.0cm # NHOH		28
Nerve-Canal Plier Curved # RFEC		70
Nerve-Canal Plier Straight # RFES		70

**O**

Ochsenbein Chisel # CO1		55,84
Ochsenbein Chisel # CO2		55,84
Ochsenbein Chisel # CO3	#6	84
Ochsenbein Chisel # CO4	#6	84
Offer Package Premium # OPP		124
Offer Package Standard # OPS		124
Oseotomes Instruments Plate Curved With Cassette # OSTMPUAWC6		122
Oseotomes Instruments Plate Straight With Cassette # OSTMPUWC5		122
Osteotome Bone Chisel Curved 4mm # OSS6519CS		55
Osteotome Bone Chisel Curved 7.5mm # OSS6521CS		55
Osteotome Bone Chisel Straight 4mm # OSS6518S		55
Osteotome Bone Chisel Straight 7.5mm # OSS6520S		55
Osteotome Bone Pusher Angulated 2.7 # OSTMPU27A		50
Osteotome Bone Pusher Angulated 3.2 # OSTMPU32A		50
Osteotome Bone Pusher Angulated 3.7 # OSTMPU37A		50
Osteotome Bone Pusher Angulated 4.2 # OSTMPU42A		50
Osteotome Bone Pusher Angulated 4.7 # OSTMPU47A		50
Osteotome Bone Pusher Angulated 5.4 # OSTMPU54A		50
Osteotome Bone Pusher Straight 2.7 # OSTMPU27		50
Osteotome Bone Pusher Straight 3.2 # OSTMPU32		50
Osteotome Bone Pusher Straight 3.7 # OSTMPU37		50
Osteotome Bone Pusher Straight 4.2 # OSTMPU42		50
Osteotome Bone Pusher Straight 4.7 # OSTMPU47		50
Osteotome Bone Pusher Straight 5.4 # OSTMPU54		50
Osteotome Tapered Concave Angulated 2.7 # OSTMSH27A		51
Osteotome Tapered Concave Angulated 3.2 # OSTMSH32A		51
Osteotome Tapered Concave Angulated 3.7 # OSTMSH37A		51
Osteotome Tapered Concave Angulated 4.2 # OSTMSH42A		51
Osteotome Tapered Concave Angulated 4.7 # OSTMSH47A		51
Osteotome Tapered Concave Angulated 5.0 # OSTMSH50A		51
Osteotome Tapered Convex Angulated 2.7 # OSTMSP27A		51
Osteotome Tapered Convex Angulated 3.2 # OSTMSP32A		51
Osteotome Tapered Convex Angulated 3.7 # OSTMSP37A		51
Osteotome Tapered Convex Angulated 4.2 # OSTMSP42A		51
Osteotome Tapered Convex Angulated 5.0 # OSTMSP50A		51

**P**

Platelet Rich Fibrin # PRF		54
Paquette Blade Handle # KPO		86
Perio Scaler in Pouch # PSPB	#3	125
Perio Surgical Instruments With Cassette # PSINC12		88, 119
Periodontal File Schluger # FS9/10S	#6	84
Periodontal File Sugarman File # FS1/2S	#6	84
Periodontal File Sugarman File # FS3/4S	#6	84
Periosteal Elevator # Beaf1	#3	40
Periosteal Elevator # P149	#6	39
Periosteal Elevator # P856	#6	40
Periosteal Elevator # P80	#6	39
Periosteal Elevator # PH2	#6	39
Periosteal Elevator # PHP3	#6	40
Periosteal Elevator Berque # PBQ1	#6	39
Periosteal Elevator Buser # PPBUSER	#6	39



ITEM NAME	HANDLE	PAGE #	ITEM NAME	HANDLE	PAGE #
Periosteal Elevator Freer # P16	#3	40	Plier Posterior Band Remover Long TC # 3000/50TC		108
Periosteal Elevator Freer # P884	#6	40	Plier Posterior Band Removing Short # 3000/49		102
Periosteal Elevator Glickman # P24G	#6	39,57	Plier Posterior Band Removing Short TC # 3000/49TC		108
Periosteal Elevator Goldman Fox # P14	#6	39	Plier Ribbon Arch Long # 3000/52		104
Periosteal Elevator Hirschfeld # P20	#6	39	Plier Ribbon Arch Long TC # 3000/52TC		107
Periosteal Elevator Howarth # PH		39,46	Plier Ribbon Arch Short # 3000/51		104
Periosteal Elevator Kramer Naviers # PKN	#3	40	Plier Ribbon Arch Short TC # 3000/51TC		107
Periosteal Elevator Kramer Naviers # PKN1	#6	39	Plier Separating Module # 3000/391		103
Periosteal Elevator Labanca # PPLABANCA2	#3,#6	40	Plier Separating Module Angular # 3000/392		103
Periosteal Elevator Molt # P883	#6	40	Plier Stand Acrylic # P5		115,128
Periosteal Elevator Molt-9 large # P9	#3	40	Plier Three Beak # 3000/57		104
Periosteal Elevator Molt-9 Small # P9S	#3	40	Plier Three Beak TC # 3000/57TC		108
Periosteal Elevator Pepila Elevator # PH2GM	#6	57	Plier Tweed O'brien # 3000/97		104
Periosteal Elevator Prichard # PPR3	#6	39	Plier Universal # 3000/55		106
Periosteal Elevator Seldin # P23		39	Plier V Bending # 3000/321		105
Periotome Flexi # PTF1		87	Plier Weingart Utility Slim TC # 3000/32TC		107
Periotome Flexi # PTF2		87	Plier Weingart Utility Short # 3000/31		102
Periotome Flexi # PTF3		87	Plier Weingart Utility Slim # 3000/32		102
Periotome Flexi # PTF4		87	Plier Weingart Utility Short TC # 3000/31TC		107
Periotome Flexi # PTF5		87	Plugger/Condenser # PLG1	#3,#4	65
Periotome Flexi # PTF6		87	Plugger/Condenser # PLG1L	#3,#4	65
Periotome Non Flexi # P1	#4	87	Plugger/Condenser # PLG2	#3,#4	65
Periotome Non Flexi # P2	#4	87	Plugger/Condenser # PLG2L	#3,#4	65
Periotome Non Flexi # P3	#4	87	Plugger/Condenser # PLG3	#3,#4	65
Periotome Non Flexi # P4	#4	87	Plugger/Condenser # PLG3L	#3,#4	65
Periotomes Set # PTS3		127	Plugger/Condenser # PLG4	#3,#4	65
PKT Carver Colour Coded in Pouch # PKTCCPS		126	Plugger/Condenser # PLG4L	#3,#4	65
PKT Carver Set # PKTSS	#3	127	PMT Set Premium # PMTP4	#4	127
Plaster Shears # PCS		99	PMT Set Standard # PMTS4	#3	127
Plastic Filling Instrument # PF2	#3,#4	65	PMT Set With Cassette # PMTWC3	#4	121
Plastic Filling Instrument # PF21	#3,#4	65	Prob University of North Carolina/William D/E # PCPUNC15/W6	#3,#6	83
Plastic Filling Instrument # PF11	#3	65	Probe CP12/Screening Probe D/E # PCP12/11.5B	#6	83
Plastic Filling Instrument # PF2	#3	65	Probe Explorer with Screening Probe # EXP23/11.5B	#6	4,83
Plastic Filling Instrument # PF3	#3	65	Probe Explorer with UNC15 # EXP23/UNC15	#6	4,83
Plastic Filling Instrument # PGFT1	#3	65	Probe Explorer with William # EXP23/QW6	#6	4,83
Plastic Filling Instrument # PGFT2	#3	65	Probe Goldman Fox # PGF	#6	82
Plastic Filling Instrument # PLGH3	#3,#4	65	Probe Goldman Fox/William D/E # PCPGF/W6	#6	83
Plier Adam's # 3000/33		106	Probe Marquis Colour Coded # PCP12	#6	82
Plier Adam's TC # 3000/33TC		108	Probe Nabers D/E # PQ2N	#3,#6,#7	83
Plier Anterior Band Remover # 3000/48		103	Probe Screening Probe # PCP11.5B	#3,#6	82
Plier Ball & Socket # 3000/4		103	Probe Thin William # POW6	#3,#4,#6	82
Plier Band Forming Left # 3000/41		106	Probe University of North Carolina # PCPUNC15	#3,#6	82
Plier Band Forming Right # 3000/42		106	Probe University of North Carolina/Screening Probe D/E # PCP15/11.5	#6	83
Plier Band Pinching Left # 3000/5		106	Probe William /Screening Probe D/E # PCPQOW/11.5B	#3,#6	83
Plier Band Pinching Right # 3000/6		106	Probe William/CP12 D/E # PCP12/QDW6	#6	83
Plier Band/Crown Crimping # 3000/225		103	Probe/Explorer D/E # EXD5	#3,#4	4
Plier Bird Beak Long # 3000/326		104	Probe/Explorer S/E # EXS11	#3	4
Plier Bird Beak Long TC Tip # 3000/326TC		107	Probe/Explorer S/E # EXS23	#3	4
Plier Bird Beak Short # 3000/325		104	Probe/Explorer S/E # EXS54	#3	4
Plier Bird Beak Short TC Tip # 3000/325TC		107	Probe/Explorer S/E # EXS55	#3	4
Plier Bracket Remover Curved # 3000/82		102	Probe/Explorer S/E # EXS59	#3	4
Plier Bracket Remover Straight # 3000/83		102	Probe/Explorer S/E # EXS6	#3	4
Plier Canine Contouring # 3000/310		103	Probe/Explorer S/E # EXS6L	#3	4
Plier Cinch Back # 3000/63		105	Probe/Explorer S/E # EXS6XL	#3,#4	4
Plier Crimpable Hook Placement # 3000/2		108	Probe/Explorer S/E # EXS8	#3	4
Plier De La Rosa Band Contouring # 3000/34		104	Proxitip Luxating Hybrids Bernard Left # EPTBER7		16
Plier De La Rosa Wire Contouring # 3000/35		104	Proxitip Luxating Hybrids Bernard Right # EPTBER6		16
Plier How Curved 2mm # 3000/72		102	Proxitip Luxating Hybrids Large Distal # EPTLGD		16
Plier How Curved 4mm # 3000/73		102	Proxitip Luxating Hybrids Large Mesial # EPTLGM		16
Plier How Curved TC 4mm # 3000/73TC		108	Proxitip Luxating Hybrids Small Curved # EPTSMC		16
Plier How Straight 2mm # 3000/84		102	Proxitip Luxating Hybrids Small Distal # EPTSMO		16
Plier How Straight 4mm # 3000/85		102	Proxitip Luxating Hybrids Small Mesial # EPTSMH		16
Plier How Straight TC 4mm # 3000/85TC		108	Proxitip Luxating Hybrids Small Straight # EPTSMS		16
Plier Jhonson Contouring # 3000/59		103	Proxitip Luxating Hybrids Straight Spade # EPTSSP		16
Plier Lingual Arch # 3000/58		105	Proxitip Luxating Hybrids Warwick James Left # EPTWJL		16
Plier Nance # 3000/80		105	Proxitip Luxating Hybrids Warwick James Right # EPTWJR		16
Plier Optical # 3000/53		106	Pterigoidous Chisel # PC		43
Plier Posterior Band Remover Long # 3000/50		102	Putty Cutter Left Handress # PCL		99

## INDEX

ITEM NAME	HANDLE	PAGE #	ITEM NAME	HANDLE	PAGE #
Putty Cutter Right Handress # PCR		99	Retractor Walsham Nasal Septum Forcep Left # WNSL		42
<b>R</b>			Retractor Walsham Nasal Septum Forcep Right # WNSR		42
Rachel Wrench Universal # RW		55	Retractor Weider Cheek & Tongue Depressor Large 15.0cm # CRW2		35
RCT Plugger Set # RCP56	#3	127	Retractor Weider Cheek & Tongue Depressor Small 14.0cm # CRW1		35
RCT Spreaders Set # RCS56	#3	127	Retractor Weider Cheek & Tongue Depressor With Suction Large 15.0cm # CRW25		36
Retractor Appuntito Blunt 16.5cm # TRAB2		38	Rhodes Back Action # C36/37	#6	84
Retractor Appuntito Sharp 16.5cm # TRAB1		38	Root Canal Plugger 1/3 .40mm/.45mm # RCP1/3	#3,#7	69
Retractor Asch Nasal Septum Forcep Large # ANS		42	Root Canal Plugger 5/7 .50mm/.75mm # RCP5/7	#3,#7	69
Retractor Austin Fork # TRA2		35	Root Canal Plugger 9/11 1.0mm/1.15mm # RCP9/11	#3,#7	69
Retractor Austin Half Moon # TRA3		35	Root Canal Plugger ISO size .30mm # RCP30	#3,#7	69
Retractor Austin Regular # TRA4		35	Root Canal Plugger ISO size .40mm # RCP40	#3,#7	69
Retractor Austin With Handle Fork # TRA1		35	Root Canal Plugger ISO size .50mm # RCP50	#3,#7	69
Retractor Austin With Suction Fork # TRA2S		36	Root Canal Plugger ISO size .60mm # RCP60	#3,#7	69
Retractor Austin With Suction Half Moon # TRA3S		36	Root Canal Plugger Luks1 .25mm # RCPL1	#3,#7	69
Retractor Barsky Pharyngeal Flap Elevator # BPF		43	Root Canal Plugger Luks2 .50mm # RCPL2	#3,#7	69
Retractor Bishop 13.5cm # CR8		36	Root Canal Plugger Luks3 .75mm # RCPL3	#3,#7	69
Retractor Bone Elevator Hohman Type # BEH		43	Root Canal Plugger Luks4 1.0mm # RCPL4	#3,#7	69
Retractor Bone Holding Forcep (Mandible) # BHM		43	Root Elevator Apex # E301		17
Retractor Bristow Elevator # BE		43	Root Elevator Apex # E302		17
Retractor Bruenings 19.0cm # TRB		36	Root Elevator Flohr Left # E2F (Standard)		19
Retractor Channel # CHR		43	Root Elevator Flohr Left # E2FP (Premium)		19
Retractor Chin # CR		43	Root Elevator Flohr Right # E3F (Standard)		19
Retractor Columbia 14.0cm # CRC		36	Root Elevator Flohr Right # E3FP (Premium)		19
Retractor Con Molla Flexible 2 16.5cm # TRCM2		38	Root Elevator Flohr Straight # E1F (Standard)		19
Retractor Con Molla Flexible 3 16.5cm # TRCM3		38	Root Elevator Flohr Straight # E1FP (Premium)		19
Retractor Con Molla Flexible 4 16.5cm # TRCM4		38	Root Elevator Bandley # E44C		18
Retractor Farabauf 12.0cm # CRF		35	Root Elevator Bein # EB1 (Standard)		18
Retractor Ferguson Gag Large # FMGA		42	Root Elevator Bein # EB1P (Premium)		18
Retractor Ferguson Gag Small # FMGC		42	Root Elevator Bein # EB2 (Standard)		18
Retractor Gillies 18.0cm # TRG1		38	Root Elevator Bein # EB2P (Premium)		18
Retractor Gillies 18.0cm # TRG2		38	Root Elevator Bein # EB3 (Standard)		18
Retractor Gillies 18.0cm # TRG3		38	Root Elevator Bein # EB3P (Premium)		18
Retractor Hartmann 15.0cm 72x24mm # CRH		35	Root Elevator Bein # EB4		18
Retractor Heister Mouth Gag 13.0cm # HMG		35	Root Elevator Bein # EB5		18
Retractor Hyton William Forcep Left # HWL		42	Root Elevator Bein # EB6		18
Retractor Hyton William Forcep Right # HWR		42	Root Elevator Cogswell # EA		18
Retractor Kilner C-Shape # KRC		43	Root Elevator Coupland # E1C		18
Retractor Kocher 40x11mm 21.0cm # CRK1		37	Root Elevator Coupland # E2C		18
Retractor Kocher 40x11mm 21.0cm # CRK3		37	Root Elevator Coupland # E3C		18
Retractor Kocher 55x11mm 21.0cm # CRK4		37	Root Elevator Coupland # EA1C (Standard)		18
Retractor Kocher 56x11mm 21.0cm # CRK2		37	Root Elevator Coupland # EA1CP (Premium)		18
Retractor Langenback 30x10 21.5cm # CRL1		37	Root Elevator Coupland # EA2C (Standard)		18
Retractor Langenback 40x10 21.5cm # CRL2		37	Root Elevator Coupland # EA2P (Premium)		18
Retractor Langenback 40x13 21.5cm # CRL3		37	Root Elevator Coupland # EA3C (Standard)		18
Retractor Langenback 50x10 21.5cm # CRL4		37	Root Elevator Coupland # EA3P (Premium)		18
Retractor Martin Simplex 10.5cm # CRMS		36	Root Elevator Cryer Left # E31 (Standard)		18
Retractor Middeldorpf 14x17mm 19.0cm # CRMD1		37	Root Elevator Cryer Left # E31P (Premium)		18
Retractor Middeldorpf 20x22mm 19.0cm # CRMD2		37	Root Elevator Cryer Left # EA30		18
Retractor Middeldorpf 26x30mm 19.0cm # CRMD3		37	Root Elevator Cryer Right # E32 (Standard)		18
Retractor Molt Mouth Gag 11.0cm Small # MGC		35	Root Elevator Cryer Right # E32P (Premium)		18
Retractor Molt Mouth Gag 14.0cm Large # MGA		35	Root Elevator Cryer Right # EA31		18
Retractor MSY 11.5cm # MSY		36	Root Elevator Heidbrink Straight # E1		17
Retractor Nasal Septum Osteotome With Guard U Shape # NSOU		43	Root Elevator Heidbrink Straight # E2		17
Retractor Oringer Large # CR03		35	Root Elevator Lindo Levain 2.5x5mm # LLS3		19
Retractor Oringer Medium # CR02		35	Root Elevator Lindo Levain 2.5x5mm # LLS6		19
Retractor Oringer Small # CR01		35	Root Elevator Lindo Levain 3x8mm # LLM2		19
Retractor Ragnell-Davis 8x4mm/15x5mm # TRRD		38	Root Elevator Lindo Levain 3x8mm # LLMCS		19
Retractor Ramus 35x11mm 21.5cm # CRR		37	Root Elevator Lindo Levain 4x8mm # LLL1		19
Retractor Row Maxillary Forcep Left # RMDL		42	Root Elevator Lindo Levain 4x8mm # LLLC4		19
Retractor Row Maxillary Forcep Right # RMDR		42	Root Elevator London Hospital # E20H (Standard)		19
Retractor Senn Miller Blunt 16.0cm # TRSM2		38	Root Elevator London Hospital # E20HP (Premium)		19
Retractor Senn Miller Sharp 16.0cm # TRSM1	38,45		Root Elevator Miller Apex Left # E71		17
Retractor Sternberg 16.0cm # CRCS		36	Root Elevator Miller Apex Left # E73		17
Retractor University of Minnesota 14.0mm # CRM		36	Root Elevator Miller Apex Right # E72		17
Retractor University of Minnesota With Suction 14.0mm # CRMS		36	Root Elevator Miller Apex Right # E74		17
Retractor Vestibulum 14.0cm # CRV		36	Root Elevator Seldin # E2S		17
			Root Elevator Seldin # E301W		17



ITEM NAME	HANDLE	PAGE #	ITEM NAME	HANDLE	PAGE #
Root Elevator Seldin # E304W		17	Rubber Dam Frame Stainless Steel Medium 5" # RDCF5		74
Root Elevator Seldin # E34		17	Rubber Dam Frame Stainless Steel Small 4" # RDCF4		74
Root Elevator Seldin # E34S		17	Rubber Dam Punch 16.0cm # RDP		73
Root Elevator Seldin # E3S		17	Rubber Dam Punch Ainsworth 17.0cm # RDPA		73
Root Elevator Seldin Left # E1L		17	Rubber Dam Sheet 5"x5" # RDLC		74
Root Elevator Seldin Left # E4L		17	Rubber Dam Sheet 6"x6" # RDLA		74
Root Elevator Seldin Right # E1R		17	Rubber Dam Template 6"x6" # RDT		74
Root Elevator Seldin Right # E4R		17			
Root Elevator with Cassette Premium # REPWC9		20, 119			
Root Elevator Premium In Pouch # REPP9		123			
Root Elevator Standard In Pouch # RESP9		123			
Root Elevator With Warwick James In Pouch # RESP10		123			
Root Fragment Steglitz # RF4S		70			
Root Fragment Steglitz # RFSTR		70			
Root Fragment Steglitz # RF90		70			
Root Tip Pick Heidbrink D/E # EHB13/14	#6	19	<b>S</b>		
Root Tip Pick Heidbrink Left # EHB2		19	Scaler Anterior # SCNEV11	#6	79
Root Tip Pick Heidbrink Right # EHB3		19	Scaler Anterior # SN137	#3,#4,#6	79
Root Tip Pick Heidbrink Straight # EHB1		19	Scaler Bilid # SB4/8	#3,#4	80
Root Tip Pick Howard D/E # EHB2/3	#6	19	Scaler Crane Kaplan # SCK6	#6	79
Rubber Dam Clamp Black Line 12A # RDCB12A		71	Scaler Curvine # SCM152	#3,#6	79
Rubber Dam Clamp Black Line 13A # RDCB13A		71	Scaler Hoe Orban # SO6/7	#3,#6	79
Rubber Dam Clamp Black Line 14A # RDCB14A		71	Scaler Hoe Orban # SO7/9	#3,#6	80
Rubber Dam Clamp Black Line 201 # RDCB201		71	Scaler Hoe Orban # SO8/9	#3,#6	79
Rubber Dam Clamp Black Line 202 # RDCB202		71	Scaler Hygienist # SH6/3	#6	79
Rubber Dam Clamp Black Line 205 # RDCB205		71	Scaler Manipal # SM5/6		80
Rubber Dam Clamp Black Line 210 # RDCB210		71	Scaler Manipal # SM137		80
Rubber Dam Clamp Black Line 211 # RDCB211		71	Scaler Manipal # SM2/3		80
Rubber Dam Clamp Black Line 2A # RDCB2A		71	Scaler Posterior # S129	#3,#6	80
Rubber Dam Clamp Black Line 3 # RDCB3		71	Scaler Posterior # S204S	#3,#6	79
Rubber Dam Clamp Black Line 4 # RDCB4		71	Scaler Posterior # SCNEW 2	#6	79
Rubber Dam Clamp Gold Line 12A # RDCG12A		71	Scaler Posterior # SN135	#3,#4,#6	79
Rubber Dam Clamp Gold Line 13A # RDCG13A		71	Scaler Remington # SRM3/4	#6	79
Rubber Dam Clamp Gold Line 14A # RDCG14A		71	Scaler Springfield # SSP2/3	#3,#6	79
Rubber Dam Clamp Gold Line 201 # RDCG201		71	Scaler Sub Gingival # SCU61		81
Rubber Dam Clamp Gold Line 202 # RDCG202		71	Scaler Sub Gingival # SEC1		81
Rubber Dam Clamp Gold Line 205 # RDCG205		71	Scaler Sub Gingival # SEC2		81
Rubber Dam Clamp Gold Line 210 # RDCG210		71	Scaler Sub Gingival # SF1		81
Rubber Dam Clamp Gold Line 211 # RDCG211		71	Scaler Sub Gingival # SF2		81
Rubber Dam Clamp Gold Line 2A # RDCG2A		71	Scaler Sub Gingival # SWV1		81
Rubber Dam Clamp Gold Line 3 # RDCG3		71	Scaler Sub Gingival # SWV2		81
Rubber Dam Clamp Gold Line 4 # RDCG4		71	Scaler Supra Gingival # SU15		81
Rubber Dam Clamp Matt Finish # RDC12A		72	Scaler Supra Gingival # SB2		81
Rubber Dam Clamp Matt Finish # RDC13A		72	Scaler Supra Gingival # SB3		81
Rubber Dam Clamp Matt Finish # RDC14A		72	Scaler Supra Gingival # SB6		81
Rubber Dam Clamp Matt Finish # RDC201		72	Scaler Supra Gingival # SC1		81
Rubber Dam Clamp Matt Finish # RDC202		72	Scaler Supra Gingival # SDP5		81
Rubber Dam Clamp Matt Finish # RDC205		72	Scaler Supra Gingival # SDP6		81
Rubber Dam Clamp Matt Finish # RDC210		72	Scaler Supra Gingival # SSP1		81
Rubber Dam Clamp Matt Finish # RDC211		72	Scaler Surface # MDY652	#3,#6	79
Rubber Dam Clamp Matt Finish # RDC2A		72	Scaler Towner-Jacquette # SU15/30	#3,#4,#6	79
Rubber Dam Clamp Matt Finish # RDC3		72	Scalpel Handle Bard Packer No.3 12.5cm # 10-130-03		34
Rubber Dam Clamp Matt Finish # RDC4		72	Scalpel Handle Bard Packer Style With Matrix Scale No.3 12.5cm # 10-130-03E		34
Rubber Dam Clamp Peedo # RDCA1		72	Scalpel Handle No.4 13.5cm # 10-100-04E		34
Rubber Dam Clamp Peedo # RDCA3		72	Scalpel Handle No.7 16.5cm # 10-100-07E		34
Rubber Dam Clamp Peedo # RDCA4		72	Scalpel Handle Round Angular 14.5cm # 10-130-5AE		34
Rubber Dam Clamp Peedo # RDCA5		72	Scalpel Handle Round Angular 14.5cm # 10-130-5AE #4	#4	34
Rubber Dam Clamp Peedo # RDCA6		72	Scalpel Handle Round Double Side Blade # 10-130-5DE		34
Rubber Dam Clamp Peedo # RDCA2		72	Scalpel Handle Round Orientable # 10-130-180		34
Rubber Dam Clamp Set Black Line # RDC08B11		72	Scalpel Handle Round Straight 14.5cm # 10-130-5EM		34
Rubber Dam Clamp Set Gold Line # RDC08G11		72	Scalpel Handle Round Straight 14.5cm # 10-130-5EM #4	#4	34
Rubber Dam Clamp Set Mat Finish # RDC0811		72	Scissor Castroviejo Large 14.0cm # Curved # S38		25
Rubber Dam Forcep Ivory 17.0cm # RDF		73	Scissor Castroviejo Large 14.0cm # Straight # S39		25
Rubber Dam Forcep Stoke 17.5cm # RDF4E		73	Scissor Castroviejo Small 12.0cm # Curved # S33		25
Rubber Dam Frame Plastic # RDAF6P		74	Scissor Castroviejo Small 12.0cm # Straight # S32		25
Rubber Dam Frame Stainless Steel Large 6" # RDAF6		74	Scissor Crown/Band Cutting 12.0cm # Curved # S0GC		23
			Scissor Crown/Band Cutting 12.0cm # Straight # S0GS		23
			Scissor Crown/Band Cutting Curved 12.0cm TC.Tip # S5039		26
			Scissor Crown/Band Cutting Straight 12.0cm TC.Tip # S5038		26
			Scissor Dean 16.5cm # S9		24
			Scissor Goldman Fox Curved 13.0cm # S16C		22
			Scissor Goldman Fox Curved 13.0cm T C # S5081		26

INDEX

ITEM NAME	HANDLE	PAGE #
Scissor Goldman Fox Double Curved 13.0cm # S10		22
Scissor Goldman Fox Double Curved 13.0cm T C # S5079		26
Scissor Goldman Fox Straight 13.0cm # S16		22
Scissor Goldman Fox Straight 13.0cm T C # S5080		26
Scissor Heath For Suture Cutting 15.5cm # S25		23
Scissor Iris Angular 11.5cm # S7		22
Scissor Iris Curved 11.5cm # S18		22
Scissor Iris Curved 11.5cm TC # S5083		26
Scissor Iris Micro Curved 11.5cm # S26C		22
Scissor Iris Micro Straight 11.5cm # S26		22
Scissor Iris Side Curved 11.5cm # S27		22
Scissor Iris Straight 11.5cm # S17		22
Scissor Iris Straight 11.5cm TC # S5082		26
Scissor Kelly Curved 16.0cm # S1		23
Scissor Kelly Straight 16.0cm # S2		23
Scissor Lestrage Saw Edge 12.0cm # S14		22
Scissor Lestrage Saw Edge 12.0cm T C # S5084		26
Scissor Lockin Saw Edge Curved Handle 16cm # S12		23
Scissor Lockin Saw Edge Straight Handle 16cm # S11		23
Scissor Marilyn Curved 10.0cm # S3548		24
Scissor Marilyn Curved 9.5cm # S3542		24
Scissor Marilyn Straight 10.0cm # S3547		24
Scissor Marilyn Straight 9.5cm # S3541		24
Scissor Mayo Curved 14.5cm # S3		23
Scissor Mayo Straight 14.5cm # S4		23
Scissor Metzenbaum 14.5cm Curved # S29		23
Scissor Metzenbaum 14.5cm Straight # S28		23
Scissor Metzenbaum Curved 12.0cm T C # S5055		26
Scissor Metzenbaum Straight 12.0cm T C # S5054		26
Scissor Micro Heath For Suture Cutting 10.0cm # S25S		23
Scissor Noyes Angular 11.0cm # S42		24
Scissor Noyes Curved 11.0cm # S31		25
Scissor Noyes Straight 11.0cm # S30		25
Scissor Quinby Curved 12.5cm # S8		22
Scissor Quinby Straight 12.5cm # S8S		22
Scissor Spencer Straight 13.0cm # S13		24
Scissor Spencer Straight 9.0cm # S13S		24
Scissor Spencer Angular 13.0cm # S13A		24
Scissor Spencer Angular 13.0cm T C # S3482		26
Scissor Spencer Angular 9.0cm # S13SA		24
Scissor Westcott Curved 11.0cm # S35		24
Scissor Spencer Angular 10.0cm T C # S3481		26
Scissor Spencer Straight 10.0cm T C # S3483		26
Scissor Spencer Straight 13.0cm T C # S3484		26
Sinus Lift Instruments in Pouch # IMPSLP9	#6	126
Sinus Lifting Instrument 1 # IMPSL1	#6	56
Sinus Lifting Instrument 2 # IMPSL2	#6	56
Sinus Lifting Instrument 3 # IMPSL3	#6	56
Sinus Lifting Instrument 4 # IMPSL4	#6	56
Sinus Lifting Instrument 5 # IMPSL5	#6	56
Sinus Lifting Instrument 6 # IMPSL6	#6	56
Sinus Lifting Instrument 7 # IMPSL7	#6	56
Sinus Lifting Instrument 8 # IMPSL8	#6	56
Sinus Lifting Instrument 9 # IMPSL9	#6	56
Spirit Lamp # ASL		128
Spreader 1 .25mm # RCSGP1	#3,#7	70
Spreader 11 .40mm # RCS011	#3,#7	70
Spreader 11S .35mm # RCS011S	#3,#7	70
Spreader 11T .30mm # RCS011T	#3,#7	70
Spreader 11TS .35mm # RCS011TS	#3,#7	70
Spreader 2 .25mm # RCSGP2	#3,#7	70
Spreader 3 .20mm # RCSGP3	#3,#7	70
Spreader 57 .20mm # RCSMA57	#3,#7	70
Sterilising Forcep Cheate 27.0cm # CF		29
Sterilising Forcep Rogge 20.0cm # RF20		29
Sterilising Forcep Rogge 30.0cm # RF30		29
Sub Gingival Scaler in Pouch # SGSP7		126

ITEM NAME	HANDLE	PAGE #
Suction Tube Coupland 1.5mm # ASPCT1		33
Suction Tube Coupland 2.0mm # ASPCT2		33
Suction Tube Coupland 2.5mm # ASPCT3		33
Suction Tube Coupland 3.0mm # ASPCT4		33
Suction Tube Coupland 5.0mm # ASPCT5		33
Suction Tube Frazier 1.0mm # ASPFR4		33
Suction Tube Frazier 1.5mm # ASPFR6		33
Suction Tube Frazier 2.0mm # ASPFR8		33
Suction Tube Frazier 2.5mm # ASPFR10		33
Suction Tube Frazier 3.0mm # ASPFR12		33
Suction Tube Frazier 3.5mm # ASPFR14		33
Supera Gingival Scaler in Pouch # SGSP6		125
Surgical Curette Lucas 84 # CL84	#3,#4,#6	86
Surgical Curette Lucas 85 # CL85	#3,#4,#6	86
Surgical Curette Lucas 86 # CL86	#3,#4,#6	86
Surgical Curette Lucas 87 # CL87	#3,#4,#6	86
Surgical Curette Lucas 88 # CL88	#3,#4,#6	86
Surgical Instruments in Pouch # SSP10		123
Surgical Instruments With Cassette # SWC10		119

T

Temporary Crown Holding Plier # CR4562		99
Tissue Forcep Adson 12.0cm # TP41		30
Tissue Forcep Adson 15.0cm # TP45		30
Tissue Forcep Adson 1x2 12.0cm # TP42		30
Tissue Forcep Adson 1x2 15.0cm # TP46		30
Tissue Forcep Adson Curved 1x2 12.0cm # TP44		30
Tissue Forcep Semken 1x2 Curved 12.5cm # TP34		30
Tissue Forcep Semken 1x2 Straight 12.5cm # TP33		30
Tissue Forcep Semken Curved 12.5cm # TP32		30
Tissue Forcep Semken Straight 12.5cm # TP31		30
Tissue Nipper # NIPS		86
Tissue Plier 11.5cm # TP36		30
Tissue Plier 11.5cm # TP37		30
Tissue Plier 12.5cm # TP35		30
Tissue Punch Curved 4mm/10.cm # TPC4		54
Tissue Punch Curved 5mm/10.cm # TPC5		54
Tissue Punch Curved 6mm/10.cm # TPC6		54
Tissue Punch Straight 4mm/10.cm # TPS4		54
Tissue Punch Straight 5mm/10.cm # TPS5		54
Tissue Punch Straight 6mm/10.cm # TPS6		54
Tooth Separator Elliot # TSE		74
Tooth Separator Ivory # TSI		74
Torquing Plier Female # 3000/19		105
Torquing Plier Male # 3000/18		105
Torquing Plier With Key # 3000/20 & 3000/21		105
Towel Clamp Backhaus 11.0cm # TC4		28
Towel Clamp Backhaus 13.0cm # TC5		28
Towel Clamp Backhaus 9.0cm # TC3		28
Trauma Kit # TRAUMA		122
Trophine Kit # IMPTS8		53
Tunnelling Knives 1 # TKN1	#6	80
Tunnelling Knives 2 # TKN2	#6	80
Tweezer Brush Holding # BHT5		6
Tweezer Curved # DP90		6
Tweezer GP Holding # EPL1S, EPL1		6
Tweezer London College # DPU17, DP2		6
Tweezer Marlam # DP31S, DP3		6
Tweezer Perry # DP1		6
Tweezer Self Locking # DP18L, DP18LS		6

V

Vista Tunneling Procedure # VISTA1	#6	88
Vista Tunneling Procedure # VISTA2	#6	88
Vista Tunneling Procedure # VISTA3	#6	88
Vista Tunneling Procedure # VISTA4	#6	88
Vista Tunneling Procedure # VISTA5	#6	88

ITEM NAME	HANDLE	PAGE #	PRODUCT CODE	HANDLE	PAGE #	PRODUCT CODE	HANDLE	PAGE #
Vista Tunneling Procedure # VISTA6	#6	88	673		113	3000/52TC		107
<b>W</b>			674		113	3000/53		106
Warwick James Elevator Left # EWJL		19	1006		112	3000/55		106
Warwick James Elevator Right # EWJR		19	1008		112	3000/57		104
Warwick James Elevator Straight # EWJS		19	1132		112	3000/57TC		108
Wax and Modelling Carver Gritman # WCG		90	1134		112	3000/58		105
Wax and Modelling Carver Hylin1 # WCH1		90	2811		114	3000/59		103
Wax and Modelling Carver Hylin2 # WCH2		90	2812		114	3000/6		106
Wax and Modelling Carver Hylin3 # WCH3		90	2815		115	3000/63		105
Wax and Modelling Carver Lecron # WCL1		90	2816		114	3000/67		110
Wax and Modelling Carver Lecron # WCL4		90	2817		114	3000/67H		110
Wax and Modelling Carver PKT # PKT1	#3	90	4870		113	3000/67S		110
Wax and Modelling Carver PKT # PKT2	#3	90	4871		114	3000/68		110
Wax and Modelling Carver PKT # PKT3	#3	90	4875		114	3000/7		110
Wax and Modelling Carver PKT # PKT4	#3	90	4982		113	3000/708		109
Wax and Modelling Carver PKT # PKT5	#3	90	10-100-04E		34	3000/72		102
Wax and Modelling Carver PKT Colour Coded # PKTCC1		90	10-100-07E		34	3000/73		102
Wax and Modelling Carver PKT Colour Coded # PKTCC2		90	10-130-03		34	3000/73TC		108
Wax and Modelling Carver PKT Colour Coded # PKTCC3		90	10-130-03E		34	3000/8		110
Wax and Modelling Carver PKT Colour Coded # PKTCC4		90	10-130-180		34	3000/80		105
Wax and Modelling Carver PKT Colour Coded # PKTCC5		90	10-130-5AE		34	3000/82		102
Wax and Modelling Carver Zahle1 # WCZ1		90	10-130-5AE #4	#4	34	3000/83		102
Wax and Modelling Carver Zahle2 # WCZ2		90	10-130-5DE		34	3000/84		102
Wax and Modelling Carver Zahle3 # WCZ3		90	10-130-5EM		34	3000/85		102
Wax Knife # DK5A		91	10-130-5EM #4	#4	34	3000/85TC		108
Wax Knife Gritman # DKG		91	1-015		57	3000/9		110
Wax Knife Spatula Fahrnenstock Large # WKFL		91	3000/100		110	3000/97		104
Wax Knife Spatula Fahrnenstock Small # WKFS		91	3000/101		110	3200/106		111
Wax Knife Spatula Lessmann Large # WKLL		91	3000/106		109	3200/109		111
Wax Knife Spatula Lessmann Small # WKLS		91	3000/107		109	3200/67B		111
Wax Spatula Curved # PS4		91	3000/109		109	3200/7		111
Wax Spatula Flexible # PS6		91	3000/110		109	3200/708		111
Wax Spatula Straight # PS3		91	3000/111		109	4876/1		114
Winter Cryer Laclose # EWL		20	3000/18		105	4876/2		114
Winter Cryer Left # EWC12L (Standard)		20	3000/19		105	672/1		113
Winter Cryer Left # EWC12LP (Premium)		20	3000/2		108	672/2		113
Winter Cryer Right # EWC12R (Standard)		20	3000/20 & 3000/21		105			
Winter Cryer Right # EWC12RP (Premium)		20	3000/210		106			
Wire Twister # WT		113	3000/220		108			
Wire Twister TC # 15.0cm # WTTC		113	3000/225		103	<b>A</b>		
<b>X</b>			3000/31		102	A5222		99
X Ray Film Hanger # AFH1		129	3000/310		103	AC504		68
X Ray Film Hanger # AFH10		129	3000/310C		107	ACS201		68
X Ray Film Hanger # AFH14		129	3000/31TC		107	ACS202		68
X Ray Film Hanger # AFH6		129	3000/32		102	ACS203		68
<b>Z</b>			3000/321		105	ACH		128
Zygomatic Bone Awl # ZBA		44	3000/325		104	AEB		128
			3000/325TC		107	AEBP		128
			3000/326		104	AEBR		128
			3000/326TC		107	AFH1		129
			3000/32TC		107	AFH10		129
			3000/33		106	AFH14		129
			3000/33TC		108	AFH6		129
			3000/34		104	AFX1		13
			3000/35		104	AFX13		13
			3000/391		103	AFX51A		13
			3000/392		103	AFX7		13
			3000/4		103	AFX74N		13
			3000/41		106	AFX79		13
			3000/42		106	ANH1		129
			3000/48		103	ANH2		129
			3000/49		102	ANS		42
			3000/49TC		108	ART1	#6	74
			3000/5		106	ART2	#6	74
			3000/50		102	ASL		128
			3000/50TC		108	APF2E		96
			3000/51		104	APF2C		96
			3000/51TC		107	ASPCT1		33
			3000/52		104	ASPCT2		33



INDEX

PRODUCT CODE	HANDLE	PAGE #	PRODUCT CODE	HANDLE	PAGE #	PRODUCT CODE	HANDLE	PAGE #	PRODUCT CODE	HANDLE	PAGE #
ASPCT3		33	CIB2	#6	60	CPR8		46	<b>D</b>		
ASPCT4		33	CIB3	#8	59	CPRD		45	DDK		73,127
ASPCT5		33	CIB3	#6	60	CPRFL		46	DDKP		73,127
ASPFR10		33	CIB4	#8	59	CPRFL		46	DE3		48
ASPFR12		33	CIB4	#6	60	CPRU		45	DE5		48
ASPFR14		33	CIB5	#8	59	CR		43	DMG		46
ASPFR4		33	CIB5	#6	60	CR1		98	DP1		6
ASPFR6		33	CIB6	#8	59	CR10		98	DP18L, DP18LS		6
ASPFR8		33	CIB6	#6	60	CR2		98	DP318, DP3		6
<b>B</b>			CIBWC6		121	CR3		98	DP90		6
BB18	#3,#4	66	CIG1	#6	60	CR4		98	DPU17, DP2		6
BB21	#3,#4	66	CIG2	#6	60	CR4562		99	DSS1.5SS		44
BB218	#3,#4	66	CIG3	#6	60	CR4566		99	DSS2.0SS		44
BB22/23	#3,#4	66	CIG4	#6	60	CR5		98	DSS2.5SS		44
BB27/29	#3,#4	66	CIG5	#6	60	CR6		98			
BBL3	#3	66	CIG6	#6	60	CR7 (for CRS)		98			
BC1331/1		53	CIGFT1	#6	62	CR8 (for CRP)		98	<b>E</b>		
BC1331/1F		53	CIGFT2	#6	62	CR9		98	E1		17
BC1331/2		53	CIGFT3	#6	62	CR8		36	E1C		18
BC1331/2F		53	CIGFT4	#6	62	CRC		36	E1F (Standard)		19
BC14614	#6	54	CIGFTMIN1	#6	62	CRCS		36	E1FP (Premium)		19
BC14615	#6	54	CIGFTMIN3	#6	62	CRF		35	E1L		17
BC15		53	CIGFTP6		125	CRH		35	E1R		17
BC35		53	CIGP6		125	CRK1		37	E2		17
BCC14616	#6	54	CIGWC6		121	CRK2		37	E20H (Standard)		19
BE		43	CP19		125	CRK3		37	E20HP (Premium)		19
Beal1	#3	40	CIP81	#8	59	CRK4		37	E2C		18
BEH		43	CIP82	#8	59	CRL1		37	E2F (Standard)		19
BF12	#3	33	CIP83	#8	59	CRL2		37	E2FP (Premium)		19
BF52	#3	33, 45	CIP84	#8	59	CRL3		37	E2S		17
BF52S	#3	33	CIP85	#8	59	CRL4		37	E301		17
BF64	#3	33	CIP86	#8	59	CRM		36	E301W		17
BFMS		43	CIPF5		125	CRMD1		37	E302		17
BHM		43	CIPWC5		121	CRMD2		37	E304W		17
BHP		44	CJWC19		120	CRMD3		37	E31 (Standard)		18
BHT5		6	CLC40C		53	CRMP		97	E31P (Premium)		18
BIC4882		55	CL84	#3,#4,#6	86	CRMP1		97	E32 (Standard)		18
BIC4883		55	CL85	#3,#4,#6	86	CRMP2		97	E32P (Premium)		18
BIC4887		55	CL86	#3,#4,#6	86	CRMP3		97	E34		17
BP6333		55	CL87	#3,#4,#6	86	CRMS		36	E34S		17
BP6373		55	CL88	#3,#4,#6	86	CRMS		36	E3C		18
BPF		43	CLC20C		53	CRD1		35	E3F (Standard)		19
BPP		44	CLC20L		53	CRD2		35	E3FP (Premium)		19
			CLC40L		53	CRD3		35	E3S		17
<b>C</b>			CLEFT		46,122	CRP		98	E44C		18
C36/37	#6	84	CLP1		96	CRR		37	E4L		17
CC7		33	CLP2		96	CRS		98	E4R		17
CC9		33	CO1		55,84	CRT		97	E71		17
CCP		46	CO2		55,84	CRT1		97	E72		17
CD 4/5	#3,#4	66	CO3	#6	84	CRT2		97	E73		17
CD89/92	#3	66	CO4	#6	84	CRT3		97	E74		17
CF		29	CP1/2	#3	68	CRT4		97	EA		18
CG1		33	CP10	#3	88	CRT5		97	EA1C (Standard)		18
CG2		33	CP24	#3,#4	68	CRT6		97	EA1CP (Premium)		18
CH876	#3,#4	86	CP3/4	#3	68	CRV		36	EA2C (Standard)		18
CH877	#3,#4	86	CP5/6	#3	68	CRW1		35	EA2P (Premium)		18
CH878	#3,#4	86	CP53/54	#3,#4,#6	68,74	CRW2		35	EA30		18
CH879	#3,#4	86	CP7	#3,#4	68	CRW2S		36	EA31		18
CHFRC		99	CP7/8	#3	68	CS3	#3	92	EA3C (Standard)		18
CHFRS		99	CP8	#3	68	CS4	#3,#4	92	EA3P (Premium)		18
CHR		43	CP9	#3	68	CS5	#3	92	EB1 (Standard)		18
CIAP5		125	CPDL	#3	45	CS6	#3	92	EB1P (Premium)		18
CIAWC5		121	CPDR		45	CVFR2/3	#3	66	EB2 (Standard)		18
CIB1	#6	59	CPE1		45	CVFR2/3H	#3,#4	66	EB2P (Premium)		18
CIB1	#6	60	CPE2		45	CVHL3S	#3,#4	66	EB3 (Standard)		18
CIB2	#8	59	CPHD		46	CVHL3W	#3	66	EB3P (Premium)		18
			CPKT		46	CVPCOA	#3	66	EB4		18
						CVWR1S	#3,#4	66	EBS		18





INDEX

PRODUCT CODE	HANDLE	PAGE #	PRODUCT CODE	HANDLE	PAGE #	PRODUCT CODE	HANDLE	PAGE #	PRODUCT CODE	HANDLE	PAGE #
IMSCA6E		117	<b>K</b>			MMRS	#6	5	OSTMPU27A		50
IMSCA6U		117	KB3/4	#6	85	MMT50SCS		5	OSTMPU32		50
IMSCA7E		117	KB5/6	#6	85	MMT50SCS6		5	OSTMPU32A		50
IMSCDEF		118	KCC		129	MMTDC/6		5	OSTMPU37		50
IMSCDL		118	KCEF		129	MMTDS5		5	OSTMPU37A		50
IMSCDRE		118	KCIMPT		129	MMTFS4		5	OSTMPU42		50
IMSCDS		118	KCK3	#6	85	MMTFS5		5	OSTMPU42A		50
ITRLDNPA		96	KCMNH		129	MMTM4		5	OSTMPU47		50
ITRLDNPL1		93	KCRE		129	MMTM5		5	OSTMPU47A		50
ITRLDNPL2		93	KCSYR		129	MMP3		5	OSTMPU54		50
ITRLDNPL3		93	KGF11	#6	85	MMP4		5	OSTMPU54A		50
ITRLDNPL4		93	KGF8	#6	85	MMP5		5	OSTMPUAWC6		122
ITRLDNPP6		126	KGF9	#6	85	MP3		128	OSTMPUWC5		122
ITRLDNPU1		93	KK15/16	#6	85	MR4972L		68	OSTMSH27A		51
ITRLDNPU2		93	KD1/2	#6	85	MR4972S		68	OSTMSH32A		51
ITRLDNPU3		93	KD12KPO3A	#6	85	MRDF		68	OSTMSH37A		51
ITRLDNPU4		93	KPQ		86	MRTU		68	OSTMSH42A		51
ITRLDPA		96	KRC		43	MS3		127	OSTMSH47A		51
ITRLDPH1		95	KT		128	MSW6		120	OSTMSH50A		51
ITRLDPH2		95				MSY		36	OSTMSP27A		51
ITRLDPH3		95	<b>L</b>			MT26H	#6	67	OSTMSP32A		51
ITRLDPH4		95	L2S		48	MT27H	#6	67	OSTMSP37A		51
ITRLDPHR1		95	L3C		48	MT28H	#3,#4,#6	67	OSTMSP42A		51
ITRLDPHR2		95	L3CA		48	MT29H	#3,#4,#6	67	OSTMSP50A		51
ITRLDPHR3		95	L3C		48	MTA4865		67			
ITRLDPHR4		95	L3C		48	MTA4866		67	<b>P</b>		
ITRLDPL1		93	L3S		48	MTA4869		67	P1	#4	87
ITRLDPL2		93	L5C		48				P13M	#3	34
ITRLDPL3		93	L5S		48	<b>N</b>			P14	#6	39
ITRLDPL4		93	LGWC7		120	NCGL		42	P149	#6	39
ITRLDPLA		96	LLL1		19	NGR		42	P16	#3	40
ITRLDPP6		126	LLC4		19	NH5020		29	P16H1.5SS		44
ITRLDPP1		95	LLM2		19	NH5021		29	P16H2.0SS		44
ITRLDPP2		95	LLMC5		19	NH5032		29	P16H2.5SS		44
ITRLDPP3		95	LLS3		19	NH5037		29	P2	#4	87
ITRLDPP4		95	LLSC6		19	NH5038		29	P20	#6	39
ITRLDPPU1		95	LWC7		120	NH5042		29	P23		39
ITRLDPPU2		95				NHA		28	P24G	#6	39.57
ITRLDPPU3		95	<b>M</b>			NHCVC		28	P2H1.5SS		44
ITRLDPPU4		95	MAL1		33.55	NHCVS		28	P2H2.0SS		44
ITRLDPU1		93	MAL2		33.55	NHCW		28	P2H2.5SS		44
ITRLDPU2		93	MCPR		45	NHD		28	P3	#4	87
ITRLDPU3		93	MDY652	#3,#6	79	NHDC		28	P4	#4	87
ITRLDPU4		93	MGA		35	NHMS024R		29.57	P4H1.5SS		44
ITRLDPUA		96	MGC		35	NHMC		112	P4H2.0SS		44
ITRLEDNPL1		94	MH1	#1	3	NHMF		112	P4H2.5SS		44
ITRLEDNPL2		94	MH2	#2	3	NHMH		28	P4HLL1.5SS		44
ITRLEDNPL3		94	MH3	#3	3	NHMH18		28	P4HLL2.0SS		44
ITRLEDNPL4		94	MH3	#3	3	NHMH		28	P4HLR1.5SS		44
ITRLEDNPP6		126	MH4	#4	3	NHMR		112	P4HLR2.0SS		44
ITRLEDNPU1		94	MHG1	#7	3	NHMS		112	P4H01.5SS		44
ITRLEDNPU2		94	MHG10	#7	3	NHDH		28	P4H02.0SS		44
ITRLEDNPU3		94	MHG2	#7	3	NPS		86	P856	#6	40
ITRLEDNPU4		94	MHG3	#7	3	NSOU		43	P883	#6	40
ITRLEDPL1		94	MHG4	#7	3				P884	#6	40
ITRLEDPL2		94	MHG5	#7	3	<b>O</b>			P8D	#6	39
ITRLEDPL3		94	MHG6	#7	3	OKSA		91	P9	#3	40
ITRLEDPL4		94	MHG7	#7	3	OKG		91	P9S	#3	40
ITRLEDPP6		126	MHG8	#7	3	OPP		124	P9D1	#6	39
ITRLEDPU1		94	MHG9	#7	3	OPS		124	PC		43
ITRLEDPU2		94	MHS		3	OSD1.5		44	PCL		99
ITRLEDPU3		94	ML		44	OSD2.0		44	PCP11.5B	#6	82
ITRLEDPU4		94	MM1	#6	5	OSD2.5		44	PCP12	#6	82
			MM2	#6	5	OSS6518S		55	PCP12/11.5B	#6	83
			MM3	#6	5	OSS6519C5		55	PCP12/QOW6	#6	83
			MMH50SCS	#4	5	OSS6520S		55	PCP1301	#6	54
			MMR	#6	5	OSS6521C5		55	PCP15/11.5	#6	83
						OSTMPU27		50	PCPGF/W6	#6	83

PRODUCT CODE	HANDLE	PAGE #	PRODUCT CODE	HANDLE	PAGE #	PRODUCT CODE	HANDLE	PAGE #	PRODUCT CODE	HANDLE	PAGE #
PCPD0W/11.5B	#3,#6	83	PS3		91	RDCF5		74	S29C		22
PCPUNC15	#3,#6	82	PS4		91	RDCG12A		71	S27		22
PCPUNC15/W6	#3,#6	83	PS6		91	RDCG13A		71	S28		23
PCR		99	PSWC12		119	RDCG14A		71	S29		23
PCS		99	PSP8		125	RDCG201		71	S3		23
PTF1	#6	87	PTAA		33	RDCG202		71	S30		25
PF2	#3,#4	65	PTS3		127	RDCG205		71	S31		25
PTF2	#6	87				RDCG210		71	S32		25
PF21	#3,#4	65	<b>R</b>			RDCG211		71	S33		25
PTF3	#6	87	R1A		21	RDCG2A		71	S3481		26
PTF4	#6	87	RBC		44	RDCG3		71	S3482		26
PTF5	#6	87	RBL		21	RDCG4		71	S3483		26
PTF6	#6	87	RBL45		21	RDC0811		72	S3484		26
PH1	#3	65	RBL90		21	RDC08B11		72	S35		24
PF2	#3	65	RBOE		21	RDC08G11		72	S3541		24
PF3	#3	65	RCP1/3	#3,#7	69	RDF		73	S3542		24
PFHS1	#1,#6	92	RCP30	#3,#7	69	RDF4E		73	S3547		24
PFHS2	#1,#6	92	RCP40	#3,#7	69	RDLA		74	S3548		24
PFHS3	#1,#6	92	RCP5/7	#3,#7	69	RDLG		74	S38		25
PGF	#6	82	RCP50	#3,#7	69	RDP		73	S39		25
PLGH3	#3,#4	65	RCP60	#3,#7	69	RDPA		73	S4		23
PGFT1	#3	65	RCP9/11	#3,#7	69	RDT		74	S42		24
PGFT2	#3	65	RCPL1	#3,#7	69	REPP9		123	S5038		26
PH		39,46	RCPL2	#3,#7	69	REPWC9	20,119		S5039		26
PH2	#6	39	RCPL3	#3,#7	69	RESP10		123	S5054		26
PH25M	#6	57	RCPL4	#3,#7	69	RESP9		123	S5055		26
PHP3	#6	40	RCPS6	#3	127	RF		21	S5079		26
PICH6	#3	67	RCSD11	#3,#7	70	RF20		29	S5080		26
PICH7	#3	67	RCSD11S	#3,#7	70	RF30		29	S5081		26
PICH8		67	RCSD11T	#3,#7	70	RF45		70	S5082		26
PKN	#3	40	RCSD11TS	#3,#7	70	RF90		70	S5083		26
PKN1	#6	39	RCSGP1	#3,#7	70	RREC		70	S5084		26
PKT1	#3	90	RCSGP2	#3,#7	70	RFES		70	S7		22
PKT2	#3	90	RCSGP3	#3,#7	70	RFSTR		70	S8		22
PKT3	#3,#4	66,90	RCSMA57	#3,#7	70	RMDL		42	S8S		22
PKT4	#3	90	RCSS6	#3	127	RMDR		42	S9		24
PKT5	#3	90	RDAF6		74	RMF		21	SAS1/2	#7	77
PKTCC1		90	RDAF6P		74	RMFM30		21	SAS1/2	#6	77
PKTCC2		90	RDC12A		72	RMFM45		21	SAS11/12	#7	77
PKTCC3		90	RDC13A		72	RMFM90		21	SAS11/12	#6	77
PKTCC4		90	RDC14A		72	RMFML		21	SAS13/14	#7	77
PKTCC5		90	RDC201		72	RW		55	SAS13/14	#6	77
PKTCCP5		126	RDC202		72				SAS15/16	#7	77
PKT55		127	RDC205		72				SAS15/16	#6	77
PLG1	#3,#4	65	RDC210		72	<b>S</b>			SAS3/4	#7	77
PLG1L	#3,#4	65	RDC211		72	S1		23	SAS3/4	#6	77
PLG2	#3,#4	65	RDC2A		72	S1.5		44	SAS5/6	#7	77
PLG2L	#3,#4	65	RDC3		72	S10		22	SAS5/6	#6	77
PLG3	#3,#4	65	RDC4		72	S11		23	SAS7/8	#7	77
PLG3L	#3,#4	65	RDCA1		72	S12		23	SAS7/8	#6	77
PLG4	#3,#4	65	RDCA2		72	S129	#3,#6	80	S82		81
PLG4L	#3,#4	65	RDCA3		72	S13		24	S83		81
PLGH3	#3,#4	85	RDCA4		72	S13A		24	S84/8	#3,#6	80
PMGF1		86	RDCA5		72	S13S		24	S86		81
PMGF2		86	RDCA6		72	S13SA		24	SC2R/2L	#3,#6,#7	78
PMP5	#6	54	RDC812A		71	S14		22	SC4R/4L	#3,#6,#7	78
PMTP4		127	RDC813A		71	S16		22	SCGC		23
PMTS4		127	RDC814A		71	S16C		22	SCGS		23
PMTWC3		121	RDC8201		71	S17		22	SG1		81
POW6	#3,#4,#6	82	RDC8202		71	S18		22	SCK6	#6	79
PPUSER	#6	39	RDC8205		71	S2		23	SCM152	#3,#6	79
PPLABANCA2	#3,#6	40	RDC8210		71	S2.0		44	SCNEV 2	#6	79
PPR3	#6	39	RDC8211		71	S2.5		44	SCNEV1	#6	79
PQ2N	#3,#6,#7	83	RDC82A		71	S2045	#3,#6	79	SCU61		81
PRF		54	RDC83		71	S25		23	SDFC		29
PS		115,128	RDC84		71	S25S		23	SDFP		29
			RDCF4		74	S26		22			



INDEX

PRODUCT CODE	HANDLE	PAGE #	PRODUCT CODE	HANDLE	PAGE #	PRODUCT CODE	HANDLE	PAGE #
SDFS		29	SSH1.5		44	TRA25		36
SDP5		81	SSH2.0		44	TRA3		35
SDP5/6K	#6	80	SSE		45	TRA3S		36
SDP6		81	SSH2.5		44	TRA4		35
SEC1		81	SSP1		81	TRAB1		38
SEC2		81	SSP2/3	#3,#6	79	TRAB2		38
SF1		81	SU15		81	TRAUMA	44, 122	
SF2		81	SU15/30	#3,#4,#6	79	TR8		36
SFX197.11		12	SWV1		81	TRCM2		38
SFX49.11		12	SWV2		81	TRCM3		38
SFX845.00		12	SYR1.8		32	TRCM4		38
SFX849.00		12	SYR2.2		32	TRG1		36
SFX897.00		12	SYRA1.8		32	TRG2		38
SFX944.00		12	SYRA2.2		32	TRG3		38
SFX944.01		12	SYRACC1.8R		32	TRRD		38
SFX945.00		12	SYRF1.8		32	TRSM1	38, 45	
SFX951.00		12	SYRI1		32	TRSM2		38
SFX959.01		12	SYRSA1.8		32	TSE		74
SFX974.00		12				TSI		74
SG1/2	#3,#6	76	<b>T</b>					
SG1/2	#7	77	T82	#3,#4	66	<b>V</b>		
SG11/12	#3,#6	76	TC3		28	VISTA1	#6	88
SG11/12	#7	77	TC4		28	VISTA2	#6	88
SG13/14	#3,#6	76	TC5		28	VISTA3	#6	88
SG13/14	#7	77	TFB		27	VISTA4	#6	88
SG15/16	#3,#6	76	TFS		27	VISTA5	#6	88
SG17/18	#3,#6	76	TKN1	#6	80	VISTA6	#6	88
SG3/4	#3,#6	76	TKN2	#6	80			
SG3/4	#7	77	TM		46	<b>W</b>		
SG5/6	#3,#6	76	TNCFM/L	#6	61	WA		68
SG5/6	#7	77	TNCFIS/M	#6	61	WCG		90
SG7/8	#3,#6	76	TNCGFT3	#6	61	WCH1		90
SG7/8	#7	77	TNCGFTM1	#6	61	WCH2		90
SG9/10	#3,#6	76	TNCPCL	#6	61	WCH3		90
SG9/10	#7	77	TNCPCM	#6	61	WCL1		90
SGSP7		126	TNCPCS	#6	61	WCL4		90
SGSP8		125	TNPLGSA	#6	61	WC21		90
SH6/7	#6	79	TNPLGH3	#6	61	WC22		90
SISP10		123	TNPLGOT	#6	61	WC23		90
SMC10		119	TP31		30	WKFL		91
SL1/2	#6	78	TP32		30	WKFS		91
SL3/4	#6	78	TP33		30	WKLL		91
SL5/6	#6	78	TP34		30	WKLS		91
SM11/12	#6	78	TP35		30	WNSL		42
SM11/12A	#6	78	TP36		30	WNSR		42
SM11/12S	#6	78	TP37		30	WT		113
SM13/14	#6	78	TP41		30	WTL		44
SM13/14S	#6	78	TP42		30	WTTG		113
SM137		80	TP44		30			
SM17/18	#6	78	TP45		30			
SM19/20	#6	78	TP46		30			
SM2/3		80	TPC4		54			
SM5/6		80	TPCS		54			
SN135	#3,#4,#6	79	TPC6		54			
SN137	#3,#4,#6	79	TPG1		31	<b>Z</b>		
S137		79	TPG2		31	ZBA		44
SO6/7	#3,#6	79	TPG3		31			
SO7/9	#3,#6	80	TPG4		31			
SO8/9	#3,#6	79	TPS4		54			
SP20		85	TPS5		54			
SP20SLEM		85	TPS6		54			
SP60616	#3	67	TPSLCOCMN		31			
SPM20		85	TPSLCOSM		31			
SPV		25, 57	TPSSTMBH		31, 57			
SQBL2	#6	80	TRA1		35			
SQMD2	#6	80	TRA2		35			
SRM3/4	#6	79						








## **ELEKTRO GENESIS INSTRUMENTS**

NO 1/24 Stebondale Street, london, United Kingdom

Authorize Distributor: Goldengate Holdings Limited, Dubai, U.A.E

 +44 7552016275

 +44 7552016275

 [info@elektrogenesis.com](mailto:info@elektrogenesis.com), [www.elektrogenesis.com](http://www.elektrogenesis.com)

**[www.elektrogenesis.com](http://www.elektrogenesis.com)**