


INFANT RADIANT WARMER



Infant radiant warmer with integrated bed and tray below for stocking essential items.

Display	<ul style="list-style-type: none"> -Large Seven segment LED display for patient temperature -Monochrome LED display for temperature setting -Large LED display for APGAR timer/ timer setting
Heater source	<ul style="list-style-type: none"> -Imported ceramic/quartz heater -Heater power – 600 Watts Max
Alarms	Audio Visual Alarms
Temperature settings	<ul style="list-style-type: none"> -High Temperature (+0.5° C of set temperature) -Low temperature (-0.5° C of set temperature) -Sensor failure -Power failure -Heater failure -Automatic cut- off -Servo Skin Mode -Skin Temperature sensor – Active -Accuracy +/- 0.1° C (31°C to 37°C) -Set temperature range : 34°C to 38°C



- Simple, elegant yet functional design. It fits easily in all NICUs.
- Easy to operate with all parameters cleared laid on an easy to view control panel.
- Very safe. Infant radiant range of warmers contain multiple safety features which include high temperature cutoffs and clear audible alarms.
- Powerful microprocessor works in the background calibrating the heater output efficiently based on the skin sensor inputs.
- Skin/Servo modes for heating provided for desired treatment setting.
- Easy baby monitoring and treating. The integrated bed can be tilted up to 150 and the side acrylic walls are easily collapsible. The heater head can be swivelled 90° either side for easy access.
- Halogen monitoring lights are integrated into the heater head.