

ELEKTRO PLASMA EXTRACTOR



Elektro Plasma Extractor

ELEKTRO Plasma Extractor is designed to Extract Blood component from centrifuged Bags. The spring loaded front panel applies pressure on the collected bag causing the liquid to contain in transfer bag. Its sturdy construction makes it durable and reliable, easy to use and portable. The unit is available in both versions automatic and manual.

Models	TPE-01	TPE-02
Sensor	-	Infra Red
Clamping Control	Manual	Motor Activated
Pressure Plate	Acrylic 10 mm Thick	Acrylic 10 mm Thick
Alarm	-	Audio Visual