



Cherish your life, Cherish your health!



ELEKTRO GENESIS ME LLC

Cherish Your Life, Cherish Your Health!

**ELECTRO
PROSTHETIC SOLUTIONS**

Who we are

Elektro Genesis ME LLC has emerged as a leading force in the medical equipment industry, rapidly growing from its modest origins in the Middle East to become a multinational powerhouse. As one of the largest distributors and suppliers of medical equipment across Africa, GCC, and the Middle East, we have established an unrivaled reputation for reliability and quality.

Our Reach:

Operating in multiple countries, including Ghana, Kenya, Uganda, Liberia, the United Arab Emirates, and Zambia, we have built a strong presence in key markets. Our extensive distribution network ensures that healthcare professionals in these regions have access to cutting-edge medical instruments and devices.

Our Commitment:

At Elektro Genesis, our unwavering mission is to provide customers with medical equipment built to the highest standards, without compromise. We take great pride in the longstanding relationships we have cultivated with our clients over the past 20 years. Their trust in our ability to deliver unmatched products and services sets us apart from the competition.





Our Mission :

Elektro Genesis is committed to delivering exceptional value to its customers worldwide by providing technologically advanced medical products of the highest quality at affordable prices. We aim to establish a comprehensive distribution network that serves clients across continents, including Asia, Europe, the Middle East, and the GCC. Through continuous innovation and a relentless focus on customer satisfaction, we strive to improve healthcare outcomes and enhance the well-being of individuals and communities globally.

Our Vision :

Our vision is to be a trusted leader in the global medical equipment industry, renowned for our unwavering dedication to excellence, integrity, and customer-centricity. We are driven to expand our product portfolio with state-of-the-art technologies that cater to the evolving needs of healthcare professionals worldwide. By cultivating strategic partnerships and leveraging our expertise, we aspire to shape the future of healthcare delivery, offering innovative solutions that positively transform patient care and outcomes. Together with our valued clients and partners, we envision a healthier and more accessible healthcare landscape spanning across continents.

Table of Contents

Foot and its accessories	7
Basic Accessories Section	12
Knee Joint	22
Hip Joint	27
Orthosis	29
Ankle Joint	33
BK/AK Foam Cover	35
Sock	38
Gel Cover	40
Gel Cover Lock	45
Other Accessories	48
Hand	58
Children Series	62
Protective Gear Series	67
Machine	74

Foot and its accessories

Sach Foot



Model Number	ELK 1D10
Size	15-29cm
Material	Polyurethane

Sach Foot Brown



Model Number	ELK 1D10-BR
Size	15-29cm
Material	Polyurethane

Double Axis Foot ---- No Split Toes



Model Number	ELK 1H41
Size	21-30cm
Material	Polyurethane

Double Axis Foot ---- Split Toes



Model Number	ELK 1H41-ST
Size	21-28cm
Material	Polyurethane

Single Axis Foot



Model Number	ELK 1H40
Size	20-30cm
Material	Polyurethane

Carbon Fiber Storage Energy Foot



Model Number	ELK 1D10-SE
Size	21-29cm
Material	Polyurethane

Carbon Fiber Chopart Foot



Model Number	ELK 1CF-001
Size	22-28cm
Material	Carbon Fiber

Carbon Fiber Syme Foot



Model Number	ELK 1CF-002
Size	22-28cm
Material	Carbon Fiber

Low Ankle Carbon Fiber Foot



Model Number	ELK 1CF-003
Size	22-28cm
Material	Carbon Fiber

High Ankle Carbon Fiber Foot



Model Number	ELK 1CF-004
Size	22-28cm
Material	Carbon Fiber

Low Ankle Carbon Fiber Foot



Model Number	ELK 1CF-005
Size	22-27cm
Material	Carbon Fiber

High Ankle Carbon Fiber Foot



Model Number	ELK 1CF-006
Size	22-27cm
Material	Carbon Fiber

High Ankle Carbon Fiber Storage foot



Model Number	ELK 1CF-007
Size	22-26cm
Material	Carbon Fiber

Carbon Fiber Foot Cover



Model Number	ELK 1FC-001
Size	22-28cm
Material	Polyurethane

Carbon Fiber Foot Cover



Model Number	ELK 1FC-002
Size	22-27cm
Material	Polyurethane

Cap For Sach Foot



Model Number	ELK CSF
Material	Polypropylene

Cap For Dynamic Foot



Model Number	ELK CDF
Material	Polypropylene
Size	S / L / M

Sach Foot Adaptor (not polished)



Model Number	ELK 2R8 = M10
Material	Stainless steel

Sach Foot Adaptor



Model Number	2R8-CNC-SS
Material	Stainless steel

Sach Foot Adaptor



Model Number	ELK 2R8-CNC-2-SS/2R8-CNC-2-Ti
Material	Stainless steel/Titanium

Double Axis Foot Adaptor



Model Number	ELK 2P10-SS / 2P10-Ti
Material	Stainless steel / Titanium

Single Axis Foot Adaptor



Model Number	ELK 2P11-SS / 2P11-Ti
Material	Stainless steel / Titanium

Basic Accessories Section

Integration Tube Adaptor



Model Number	ELK 2P2-SS / 2P2-Ti
Material	Aluminum tube + Stainless steel adaptor
	Aluminum tube + Titanium adaptor
Length	200 mm

Integration Tube Adaptor



Model Number	ELK 2P3-SS / 2P3-Ti
Material	Aluminum tube + Stainless steel adaptor
	Aluminum tube + Titanium adaptor
Length	400 mm

Carbon Fiber Integration Tube Adaptor



Model Number	ELK 2C1-220 / 2C1-420
Material	Carbon fiber tube + Stainless steel adaptor
Length	220mm / 420 mm

Carbon Fiber Integration Tube Adaptor



Model Number	ELK 2C2-220 / 2C2-420
Material	Carbon fiber tube + Titanium adaptor
Length	220mm / 420 mm

Aluminum Integration Tube Adaptor



Model Number	ELK 2A2-AL
Material	Aluminum tube + Aluminum adaptor
Length	20mm / 400 mm

Measuring tube



Model Number	ELK 4P10-M
Material	Aluminum

Lock Tube Adaptor



Model Number	ELK 4P21-SS / 4P21-Ti
Material	Stainless Steel / Titanium

Bent Tube Adaptor



Model Number	ELK 4P24-SS / 4P24-Ti
Material	Stainless Steel / Titanium

Pyramid Tube Adaptor



Model Number	ELK 4P25-SS / 4P25-Ti
Material	Stainless Steel / Titanium

Adjustable Short Tube Adaptor



Model Number	ELK 4P26-SS / 4P26-Ti / 4P26-A
Material	Stainless Steel / Titanium / Aluminum
Length	45 / 60 / 75 / 90 mm

Bi-directions Tube Adaptor



Model Number	ELK 4P27-SS / 4P27-Ti / 4P27-A
Material	Stainless Steel/ Titanium / Aluminum
Length	32 / 45 / 60 / 75 mm

Pyramid Bi-directions Tube Adaptor



Model Number	ELK 4P28-SS / 4P28-Ti / 4P28-A
Material	Stainless Steel/ Titanium / Aluminum
Length	32 / 45 / 60 / 75 mm

Lock Tube Adaptor Slidable



Model Number	ELK 4P29-SS / 4P29-Ti
Material	Stainless Steel/ Titanium

Double Adaptor Slidable



Model Number	ELK 4P30-SS / 4P30-Ti
Material	Stainless Steel/ Titanium

Adjustable Height of Tube Adaptor



Model Number	ELK 4P31-SS / 4P31-Ti
Material	Stainless Steel/ Titanium

Female Four Jaws



Model Number	ELK 4P43-SS / 4P43-Ti
Material	Stainless Steel/ Titanium

Male Four Jaws



Model Number	ELK 4P44-SS / 4P44-Ti
Material	Stainless Steel/ Titanium

Integration Male Three Jaws



Model Number	ELK 4P45-SS
Material	Stainless Steel

Female Three Jaws



Model Number	ELK 4P41-SS / 4P41-Ti
Material	Stainless Steel/ Titanium

Male Three Jaws



Model Number	ELK 4P42-SS / 4P42-Ti
---------------------	-----------------------

Material	Stainless Steel/ Titanium
-----------------	---------------------------

Round Adjustable Plate



Model Number	ELK 4P22-R-SS / 4P22-R-Ti
Material	Stainless Steel/ Titanium

Square Adjustable Plate



Model Number	ELK 4P22-SS / 4P22-Ti
Material	Stainless Steel/ Titanium

Square Adapter with Pyramid



Model Number	ELK 4P23-SS / 4P23-Ti
Material	Stainless Steel/ Titanium

Square Adapter with Pyramid



Model Number	ELK 4P23-1-SS / 4P23-1-Ti
---------------------	---------------------------

Material	Stainless Steel/ Titanium
-----------------	---------------------------

Male Four Jaws Plate



Model Number	ELK 4P63-M-SS / 4P63-M-Ti
Material	Stainless Steel/ Titanium

Female Four Jaws Plate



Model Number	ELK 4P63-F-SS / 4P63-F-Ti
Material	Stainless Steel/ Titanium

Rotatable Female Four Jaws Plate



Model Number	ELK 4P63-RF-SS / 4P63-RF-Ti
Material	Stainless Steel/ Titanium

Rotatable Male Four Jaws Plate



Model Number	ELK 4P63-RM-SS / 4P63-RM-Ti
Material	Stainless Steel/ Titanium

Elongated Rotatable Four Jaws Plate



Model Number	ELK 4P63-L
Material	Stainless Steel
Length	40 / 55 / 70 mm

Elongated Rotatable Female Three Jaws



Model Number	ELK 4P41-L
Material	Stainless Steel
Length	40 / 55 / 70 mm

Elongated Rotatable Square Plate



Model Number	ELK 4P65-L
Material	Stainless Steel
Length	40 / 55 / 70 mm

Adjustable Rotatable Square Plate



Model Number	ELK 4P65-SS / 4P65-Ti
Material	Stainless Steel / Titanium

Adjustable Rotatable Square Plate-Pyramid



Model Number	ELK 4P65-RM-SS / 4P65-RM-Ti
Material	Stainless Steel / Titanium

Rotator with Male Adaptor



Model Number	ELK 4P60
Material	Stainless Steel

Rotator with Male Adaptor



Model Number	ELK 4P60-1
Material	Stainless Steel

Sach Foot Adaptor Aluminum



Model Number	ELK 2R8-AL
Material	Aluminum

Male Four Jaws Aluminum



Model Number	4P63-M-AL
Material	Aluminum

Lock Tube Adaptor Aluminum



Model Number	ELK 4P21-AL
---------------------	-------------

Material	Aluminum
-----------------	----------

Square Adjustable Plate Aluminum



Model Number	ELK 4P22-AL-1
Material	Aluminum

Square Adjustable Plate Aluminum



Model Number	ELK 4P22-AL-2
Material	Aluminum

Thigh Connection Plate Aluminum



Model Number	ELK 4P19-BCP-1-AL
Material	Aluminum

Lap Clip Plate



Model Number	ELK 4P19-BCP
Material	Stainless steel

Crus Clip Plate



Model Number	ELK 4P19-CCP
---------------------	--------------

Material	Stainless steel
-----------------	-----------------

Thermoplastic Cavity Connection Plate



Model Number	ELK 4P19-BKTC
Material	Stainless steel

Wood Connection Seat



Model Number	ELK WCS-1
Material	Polyurethane

Knee Joint

Four Bar Knee Joint



Model Number	ELK 4P20-SS / 4P20-Ti
Material	Stainless steel/Titanium

Bearing Lock Knee Joint



Model Number	ELK 3P15-SS / 3P15-Ti
Material	Stainless steel/Titanium

Broken Knee Joint




Model Number	ELK 3P21-SS / 3P21-Ti
Material	Stainless steel/Titanium

Broken Knee Joint with lock



Model Number	ELK 3P21-L-SS / 3P21-L-Ti
Material	Stainless steel/Titanium

Single Axis Knee Joint with Adjustable Constant Friction



Model Number	ELK 3P14-SS
Material	Stainless steel
Weight	400g

Single Axis Knee Joint with Manual Lock



Model Number	ELK 3P17-SS / 3P17-Ti
Material	Stainless steel / Titanium

Single Axis Knee Joint without Lock



Model Number	ELK 3P17-NL-SS / 3P17-NL-Ti
Material	Stainless steel / Titanium

Lock Four Axis Knee Joint



Model Number	ELK 3P35-SS / 3P35-Ti
Material	Stainless steel / Titanium

Lock Broken Knee Joint



Model Number	ELK 3P35-LB
Material	Stainless steel

Modular Lightweight Single Axis Knee Joint-With Manual Lock



Model Number	ELK 3R40
Material	Aluminum

Four Bar Mechanical Knee Joint



Model Number	ELK 3M10A-1
Material	Aluminum alloy

Four Bar Mechanical Knee Joint



Model Number	ELK 3M10A-2
Material	Aluminum alloy

Four Bar Mechanical Knee Joint



Model Number	ELK 3M10A-3
Material	Aluminum alloy

Four Bar Pneumatic Knee Joint



Model Number	ELK 3P28A
Material	Aluminum alloy

Four Bar Pneumatic Knee Joint



Model Number	ELK 3P29
Material	Aluminum alloy

Four Bar Pneumatic Knee Joint Broken Adaptor



Model Number	ELK 3P29A
Material	Aluminum alloy

Five Bar Pneumatic Knee Joint



Model Number	ELK 3P31
Material	Aluminum alloy

Five Bar Pneumatic Knee Joint



Model Number	ELK 3P32
Material	Aluminum alloy

Hydraulic Knee Joint



Model Number	ELK 3P30
Material	Aluminum alloy

Hydraulic Knee Joint



Model Number	ELK 3P33
Material	Aluminum alloy

Hip Joint

Four Bar Hip Joint



Model Number	ELK 7E6
Material	Stainless steel / Titanium

Single Axis Hip Joint



Model Number	ELK 7E7
Material	Stainless steel

Four Bar Hip Joint



Model Number	ELK 7E5
Material	Stainless steel

Four Bar Hip Joint



Model Number	ELK 7E4
Material	Stainless steel

Aluminum Four Bar Hip Joint



Model Number	ELK 7E8
Material	Aluminum Alloy

Orthosis

Spring Lock




Model Number	Bar With	Material
ELK 17M20=17SS	17mm	Stainless steel
ELK 17M20=19SS	19mm	Stainless steel
ELK 17M20=17AL	17mm	Aluminum
ELK 17M20=19AL	19mm	Aluminum

Ring Lock



Model Number	Bar With	Material
ELK 17B38=17SS	17mm	Stainless steel
ELK 17B38=19SS	19mm	Stainless steel
ELK 17B38=17AL	17mm	Aluminum
ELK 17B38=19AL	19mm	Aluminum

Ring Lock -- Children



Model Number	Bar With	Material
ELK 17B38-OLSS	13mm	Stainless steel
ELK 17B38-OLAL	13mm	Aluminum

Ring Lock --Children (One-piece type)



Model Number	Bar With	Material
ELK 17B38-13SS	13mm	Stainless steel

Spring Lock



Model Number	Bar With	Material
ELK 17M21-17SS	17mm	Stainless steel
ELK 17M21-19SS	19mm	Stainless steel

Rear Lock



Model Number	Bar With	Material
ELK 17M22=17SS	17mm	Stainless steel
ELK 17M22=19SS	19mm	Stainless steel
ELK 17M22=17AL	17mm	Aluminum
ELK 17M22=19AL	19mm	Aluminum

Spring Lock



Model Number	Bar With	Material
ELK 17P23-17SS	17mm	Stainless steel
ELK 17P23-19SS	19mm	Stainless steel

Spring Lock



Model Number	Bar With	Material
ELK 17P36-19SS	19mm	Stainless steel

Drop Lock with Bead



Model Number	Bar With	Material
ELK 17B38-17BSS	17mm	Stainless steel
ELK 17B38-19BSS	19mm	Stainless steel

Children Hip Joint Hinge



Model Number	ELK HJH-1
Material	Stainless steel
Size	S / M / L

Traditional Calf Hinge



Model Number	ELK TCH-1
Material	Stainless steel

Elbow Hinge



Model Number	ELK ELH-1
Material	Stainless steel

Ankle Joint

Metal Ankle Hinge



Model Number	ELK 17AH1-S / 17AH1-L
Material	Stainless steel
Size	S / L

Bidirectional Herringbone Ankle Joint



Model Number	ELK 17AJ2-S / 17AJ2-L
Material	Stainless steel
Size	S / L

Adult Bidirectional Adjustment Ankle Joint



Model Number	ELK 17BA3-L
Material	Stainless steel

Child Bidirectional Adjustment Ankle Joint



Model Number	ELK 17BA4-S / 17BA4-L
Material	Stainless steel/Aluminum
Size	S (14mm) / L (17mm)

Fixed Ankle Joint



Model Number	Bar With	Side	Info
ELK 17B62-16SS	16mm	Left & Right	Bend & Straight
ELK 17B62-20SS	20mm	Left & Right	Bend & Straight

BK/AK Foam Cover

BK Cosmetic Foam Cover (Ordinary)



Model Number	ELK BKFC-6P17
Material	Foam / Sponge

BK Rubber and Plastic Cosmetic Foam Cover



Model Number	ELK BKFC-6P18RP
Material	Rubber and Plastic

BK Cosmetic Foam Cover (Strong)



Model Number	ELK BKFC-6P18
Material	Foam / Sponge

BK Cosmetic Foam Cover (Water Proof)



Model Number	ELK BKFC-6P19
Material	EVA

AK Cosmetic Foam Cover (Ordinary)



Model Number	ELK AKFC-6P20
Material	Foam / Sponge

AK Cosmetic Foam Cover (Strong)



Model Number	ELK AKFC-6P21
Material	Foam / Sponge

AK Cosmetic Foam Cover (Water Proof)



Model Number	ELK AKFC-6P22
Material	EVA

AK Cosmetic Foam Cover (Pre-shape)



Model Number	ELK AKFC-6P23
Material	Foam / Sponge

Prosthetic Leg EVA Liner



Model Number	ELK EVAL001-AK/EVAL001-BK
Material	EVA

Sock

Cotton Stump Sock



Model Number	ELK 3PSS-001
Material	Cotton
Length	30 / 40 / 50 / 70 cm

Nylon Stump Sock



Model Number	ELK 3PSS-002
Material	Nylon
Length	30 / 40 / 50 / 70 cm

Cosmetics Sock



Model Number	ELK 3PWGS-002
Material	Glass Fiber
Color	White/Black

Cosmetics Sock



Model Number	ELK 3PWWS-003
Material	Woolr
Size	30/36/40/45cm

Velvet Beauty Decoration BK/AK Sock



Model Number	ELK 3PCCS-BK / 3PCCS-AK
Material	Velvet
Length	45 / 60 cm
Color	Beige/Brown

Gel Sock



Model Number	ELK 3PGS-001
Material	Gel
Size	S/ L
Color	White/Black

Gel Cover

Universal Gel Cover (with lock/ no lock)



Model Number	ELK GCL-GPDT (with lock) / GCNL-GPFR
Material	Gel

Longitudinal Inelastic Gel Cover (with lock/ no lock)



Model Number	Viscous: GCL-AKDT(with lock) / GCNL-AKFR
	High Density: ELK GCL-AKDTHD (with lock) /GCNL-AKFRHD
Material	Gel

Antibacterial Simple Type Gel Cover (with lock/ no lock)



Model Number	ELK GCL-ELDT (with lock) / GCNL-ELFR
Material	Gel

Crus Special Gel Cover with lock



Model Number	ELK GCL-EZF
Material	Gel

High-Density Eco Gel Cover (with lock/ no lock)



Model Number	ELK GCL-ECDTHD (with lock) /GCNL-ECFRHD
Material	Gel

Gel Cover with Rings of Seal (with lock/ no lock)



Model Number	ELK GCL-SADT (with lock) /GCL-SBFR (with lock) /GCNL-SAFR
Material	Gel

Sweat-Proof Breathable Gel Cover with Lock



Model Number	ELK GCL-VSDT
Material	Gel

Superior Performance Gel Cover (with lock/ no lock)



Model Number	Viscous: GCL-SPDT (with lock) / GCNL-SPFR
	High Density: ELK GCL-SPDTHD (with lock) /GCNL-SPFRHD
Material	Gel

Thickness Gradient Gel Cover (with lock/ no lock)



Model Number	ELK GCL-APDT-TP (with lock) /GCNL-APFR-TP
Material	Gel

Gel Curved Sleeve



Model Number	ELK GCS-SFX
Material	Gel

Fabric Reinforced Suspension Sleeve



Model Number	ELK GSS-SFS
Material	Gel

LSL-LSC Silicone (with lock/ no lock)



Model Number	ELK SCL-LSL (with lock) /SCNL-LSC
Material	Gel

Gel Stump Sock



Model Number	ELK GS-ENCP
Material	凝胶

Silicone Pad



Model Number	ELK PDS
---------------------	---------

Material	Gel
Size	S / M / L / XL / XXL

Gel Cover Lock

Gel Cover Lock



Model Number	ELK SL001
Material	Aluminum
Color	Blue / Black/Red
Loud weight	100kg

Gel Cover Lock



Model Number	ELK SL002
Material	Aluminum
Color	Black
Loud weight	100kg

Gel Cover Lock



Model Number	ELK SL003
Material	Plastic+Stainless steel
Color	Black
Loud weight	100kg

Flat hanging Gel Cover Lock with Plate



Model Number	ELK P-496-T
Material	Aluminum
Color	Blue
Loud weight	100kg

Flat hanging Gel Cover Lock no plate



Model Number	ELK P-496-W
Material	Aluminum
Color	Blue
Loud weight	100kg

Flat hanging Gel Cover Lock with Pyramid



Model Number	ELK P498-P
Material	Aluminum
Color	Blue
Loud weight	100kg

Gel Cover Lock 601



Model Number	ELK GCL601
Material	Stainless steel
Color	Silver
Loud weight	100kg

Gel Cover Lock 602



Model Number	ELK GCL602-SS / GCL602-AL
Material	Stainless steel / Aluminum
Color	Silver
Loud weight	100kg

Gel Cover Lock



Model Number	ELK GCL001
Material	Aluminum
Color	Silver
Loud weight	100kg

Gel Cover Lock



Model Number	ELK GCL002
Material	Aluminum
Color	Silver
Loud weight	100kg

Gel Cover Lock



Model Number	ELK GCL003
Material	Aluminum
Color	Silver
Loud weight	100kg

Other Accessories

INK Pencil



Model Number	ELK INK-001
Color	red,blue,green, purple

Valve Seat with Tube-Manual



Model Number	ELK VM1
Material	Polyurethane
Kind	maunal valve with tube

Valve Seat with Tube-Automatic



Model Number	ELK INK-001
Material	Polyurethane
Kind	automatic valve with tube

Manual Valve Kit for Thermo Socket



Model Number	ELK VMK-3
Material	Polyurethane

Automatic Valve Kit for Thermo Socket



Model Number	ELK AVK-4
Material	Polyurethane

Screw Valve



Model Number	ELK SV-1
Material	Nylon

Seat Valve



Model Number	ELK VM3
Material	Polyurethane
Weight	20g

Plastic Inside Cup for AK Thermo Socket



Model Number	ELK PIC
Material	Polyurethane
Size	S/M/L

EVA-Sheet(Inner Plate for Prosthetic)



Model Number	ELK EVAS-001
Material	EVA
Size	1500*1000*mm
Thickness	3/4/5/6mm

EVA-Sheet(Inner Plate for Orthotic)



Model Number	ELK EVAS-OR-2
Material	EVA/ flower color
Size	1400*1000*mm
Thickness	2/3/4mm
we can custom the size , thickness, color and hardness	

EVA-Sheet for Shoe Pad



Model Number	ELK EVAS-SP	
Material	EVA	
Size	330*260mm	
Thickness	30mm	
we can custom the size , thickness, color and hardness		

PVA Sheeting



Model Number	ELK PVAS-001	
Material	PVA	
Size	1.1*0.00008m	
Length	50/100m	

PVA Bags




Model Number	ELK PVAB	
700*190*50mm	1000*190*50mm	
700*270*50mm	1000*260*50mm	
600*130*40mm	1000*500*50mm	
10PCS /Box		

PP Sheets



Model Number	ELK PPS-001	
Material	Polypropylene	
Size	300*300*12mm	1500*1000*3mm
	300*300*15mm	1500*1000*4mm
	400*400*12mm	1500*1000*5mm
	400*400*15mm	1500*1000*6mm

PE Sheets



Model Number	ELK PES-001
Material	Polyethylene
Size	2440*1220*3mm
	2440*1220*4mm
	2440*1220*5mm

Sleeve



Model Number	ELK PF-001S
Size	S,M,L,XL,XXL
Color	Blue , Green

Nylon square buckle for Orthotic



Model Number	ELK NBO-001
Material	Nylon
Size	2.5cm (100pcs)/box
	3.8cm (100pcs)/box
Eight colors: deep blue, pink, yellow, green, light blue, orange, white, black color; Various styles to choose from	

Set of cloth Velcro



Model Number	ELK CV-002
Material	Velcro
Size	2.5*28cm (10pcs)/box
	3.8*38cm (10pcs)/box
Nine colors: deep blue, pink, yellow, green, light blue, brown, orange, white, black color; Various styles to choose from	

Plate cloth orthopedic tape (Matte)



Model Number	ELK COT-003
Material	Velcro
Size	2.5cm*25m
	3.8cm*25m
Nine colors: deep blue, pink, yellow, green, light blue, brown, orange, white, black color; Various styles to choose from	

Wrench



Model Number	ELK MK-1
Material	steel
Size	4-6 / 7-9 / 8-10mm

copper rivet



Model Number	ELK CR-1
Material	steel
MOQ	100sets

Double Cap Rivet



Model Number	ELK DCR-1
Material	steel
MOQ	100sets

Single Cap Rivet



Model Number	ELK SCR-2
Material	Steel
MOQ	100sets

Plaster File Flat



Model Number	ELK PFF
Material	steel

Plaster File half Round



Model Number	ELK PFHR
Material	steel

Plaster File Round



Model Number	ELK PFR
Material	Nylon

Stockinette Cotton



Model Number	ELK COS
Material	Cotton
Size	6, 9,12,15,23,33cm

Stockinette Glass Silk



Model Number	ELK GSS
Material	Glass silk
Size	6,9,13,18,25cm

Stockinette Perlon



Model Number	ELK PES
Material	Terylene
Size	6,9,10,12,15,20,22,25cm

Stockinette Fiber Carbon



Model Number	ELK FCS
Material	Fiber carbon
Size	7,9 cm

Stockinette Nylon




Model Number	ELK NYS
Material	Nylon
Size	6cm skin color

Carbon fiber cloth



Model Number	ELK CFCL
Material	Carbon fiber

Print Paper



Model Number	ELK PPR
Material	Paper
A variety of patterns to choose from	

Elastic bandage



Model Number	ELK ELB
Size	10*450cm
	15*450cm
Pack	Type:roll

Plaster bandage



Model Number	ELK PLB
Size	15*460cm
Pack	72 rolls/ carton

Conical Grinding Roller



Model Number	ELK CGR-1

Sand Sleeve



Model Number	ELK SASL-1

Tools Grinding Head



Model Number	ELK TGH-1 / TGH-2

Tools Grinding Head



Model Number	ELK TGH-3

Silicone Foot Cover



Model Number	ELK SFC
Material	silicone
Color can be selected by color table.	

Customized Silicone Foot Cover (as the patient's actual sizes)



Model Number	ELK SFC-C
Material	Polyurethane
Color can be selected by color table.	

Cosmetic Gloves (with Bolt)



Model Number	ELK CGB
Material	silicone
Color can be selected by color table.	

Cosmetic Gloves (with Zipper and Filler)



Model Number	ELK CGFZ
Material	silicone
Color can be selected by color table.	

Customized Silicone Gloves (as the patient's actual sizes)



Model Number	ELK CSG-C
Material	silicone
Color can be selected by color table.	

Silicone finger set



Model Number	ELK SF
Material	silicone
Color can be selected by color table.	

Customized Silicone finger



Model Number	ELK SF-C
Material	silicone
Color can be selected by color table.	

Hand

Cable controlled mechanical hand for the below-elbow



Model Number	ELK FS-001S
Weight	215g
Length	140mm
Wrist circumference	144mm
Self-power source, using traction cable to control the opening and closing of artificial hand, wrist joint can be passively flexed and stretched and rotated.	

Shell elbow joint device, passive and self-lock



Model Number	ELK FS-003K
Weight	359g
Forearm length	280mm
upper arm length	143mm

Cable controlled mechanical hand for AE



Model Number	ELK FS-002S
Weight	430g
Forearm length	350mm
upper arm length	143mm

Fore Arm Cosmetic Hand System



Model Number	ELK FX-081
Weight	162g
Length	128mm
Grip	> 6N
Wrist circumference	144mm
Fingers can open passively. Wrist can be rotated passively 360° .	

Fore Arm Cosmetic Hand System with adaptor connector



Model Number	ELK FX-081L
Weight	175g
Length	153mm
Grip	> 6N
Wrist circumference	144mm
Fingers can open passively. Wrist can be rotated passively 360° .	

Arm skeleton for elbow disarticulation



Model Number	ELK FX-082GX
Weight	266g
Length	345mm
Forearm length	310mm
upper arm length	35mm
Fingers can open passively. Wrist can be rotated passively 360° . Elbow joint can be rotated passively.	

Cosmetic hand system for elbow disarticulation



Model Number	ELK FX-082G
Weight	424g
Length	470mm
Grip	> 6N
Forearm length	435mm
upper arm length	35mm
Fingers can open passively. Wrist can be rotated passively 360° . Elbow joint can be rotated passively.	

Above Elbow cosmetic hand system



Model Number	ELK FX-082B
Weight	515g
Length	750mm
Grip	> 6N
Forearm length	435mm
upper arm length	315mm
Upper arm circumference	143mm
Fingers can open passively. Wrist can be rotated passively 360° . Elbow joint can be rotated passively.	

Above Elbow arm skeleton



Model Number	ELK FX-082BX
Weight	375g
Length	620mm
Forearm length	310mm
upper arm length	315mm
Upper arm circumference	143mm
Fingers can open passively. Wrist can be rotated passively 360° . Elbow joint can be rotated passively.	

Cosmetic hand system for shoulder disarticulation



Model Number	ELK FX-083Y
Weight	565g
Length	780mm
Grip	>6N
Forearm length	435mm
upper arm length	345mm
Fingers can open passively. Wrist can be rotated passively 360° . Elbow joint and shoulder joint can be rotated passively.	

Arm skeleton for shoulder disarticulation



Model Number	ELK FX-083JG
Weight	408g
Length	655mm
Forearm length	310mm
upper arm length	345mm
Fingers can open passively. Wrist can be rotated passively 360° . Elbow joint and shoulder joint can be rotated passively	

Fore Arm Myoelectric Hand System with three fingers



Model Number	ELK FX--061S
Size	S / M

Myoelectric Hand System of Carbon Fiber with five fingers for Wrist Disarticulation



Model Number	ELK FX--071W
Size	S / M

Children Series

Aluminum Four Bar Knee Joint for Children



Model Number	ELK 3P65
Material	Aluminum alloy
Weight	380g
Loud weight	45kg
Bending angle	175°

Single Axis Knee Joint for Children No Lock



Model Number	ELK 3P66
Material	Aluminum alloy
Weight	175g
Loud weight	45kg
Bending angle	150°

Four Bar Knee Disarticulation for Children



Model Number	ELK 3P67
Material	Aluminum alloy
Weight	390g
Loud weight	45kg
Bending angle	150°

Four Bar Knee Joint for Children



Model Number	ELK 3P65SS
Material	Stainless Steel
Weight	330g
Loud weight	45kg
Bending angle	120°

Children Four Bar Hip Joint



Model Number	ELK C7E6-A
Material	Aluminum alloy
Weight	480g
Loud weight	45kg
Bending angle	175°

Children Four Bar Hip Joint



Model Number	ELK C7E6-S
Material	Stainless steel
Weight	550g
Loud weight	45kg
Bending angle	175°

Foot Adaptor for Children



Model Number	ELK 2R8=M8
Material	Stainless steel
Loud weight	45kg

Four Jaws for Children



Model Number	ELK 3S63
Material	Stainless steel
Loud weight	45kg

Lock Tube Adaptor for Children



Model Number	ELK 3S21
Material	Stainless steel
Loud weight	45kg

Socket Square Plate for Children



Model Number	ELK 3S22
Material	Stainless steel
Loud weight	45kg

Square Adapter with Pyramid for Children



Model Number	ELK 3S23
Material	Stainless steel
Loud weight	45kg

Pyramid Tube Adaptor for Children



Model Number	ELK 3S25
Material	Stainless steel
Loud weight	45kg

Integrated Tube for Children



Model Number	ELK 3S2. 3S4
Material	Aluminum alloy
3S2	Length 160mm
3S4	Length 310mm
Loud weight	45kg

Prosthetic Sach Foot



Model Number	ELK 1D10(yellow)
Weight	140g-350g
Size	15-21cm
Loud weight	50-75kg

Short Tube Adaptor for Children



Model Number	ELK 3S26
Material	Stainless steel
Length	38/46/54/62mm

Bi-directions Tube Adaptor for Children



Model Number	ELK 3S27
Material	Stainless steel
Length	22/30/38mm

Double Pyramid Tube for Children



Model Number	ELK 3S28
Material	Stainless steel
Length	22/30/38mm

Lock Tube Adaptor for Children



Model Number	ELK 3S24
Material	Stainless steel
Weight	90g

Children Plastic Inside Cup for AK Thermo Socket



Model Number	ELK SPIC
Material	Polyurethane

Rotatable Female Four Jaws for children



Model Number	ELK 3S63-F
Material	Stainless steel
Weight	85g

AK Cosmetic Foam Cover For Children(Pre-shape)



Model Number	ELK AKFC-6P23-36
Material	Polyurethane
Weight	295g
Height	70cm

Protective Gear Series

Head and cervical orthosis



Model Number	ELK POR-HC-1
Specification	Free size
Kinds	Female/Male

Head and cervical orthosis



Model Number	ELK POR-HC-2
Specification	Free size
Kinds	Female/Male

Soft Waist belt



Model Number	ELK POR-WS-1
Specification	XL/L/M
Kinds	Female/Male

Gray Waist belt



Model Number	ELK POR-WG-2
Specification	L/M/S
Kinds	Female/Male

Shoulder abduction



Model Number	ELK POR-SA-1
Specification	L/S
Kinds	Female/Male

Shoulder abduction



Model Number	ELK POR-SA-2
Specification	L/S
Kinds	Female/Male

Inflatable Shoulder protector



Model Number	ELK POR-ISP-1
Specification	L/M
Kinds	Female/Male

Forearm Support



Model Number	ELK POR-FAS-1
Specification	Free size, Left/Right
Kinds	Female/Male

Hand Separator



Model Number	ELK POR-HS-1
Specification	L/S Left/Right
Kinds	Female/Male

Hand Exerciser



Model Number	ELK POR-HE-2
Specification	L/S Left/Right
Kinds	Female/Male

Children Hip Orthosis



Model Number	ELK POR-CHO-1
Specification	Free size
Kinds	Boy/Girl

Children Hip Orthosis



Model Number	ELK POR-CHO-2
Specification	Free size
Kinds	Boy/Girl

shoulder abduction orthopedic immobilizer



Model Number	ELK POR-SAO-1
Specification	Free size
Kinds	Female/Male

Adjustable Knee orthosis



Model Number	ELK POR-AKO-1
Specification	L/M
Kinds	Female/Male

Knee Immobilizer with a straight position



Model Number	ELK POR-KISP-1
Specification	XL/L/M/S
Kinds	Female/Male

Adjustable Knee orthosis



Model Number	ELK POR-AKO-1
Specification	L / M
Kinds	Female/Male

Adjustable Knee orthosis



Model Number	ELK POR-AKO-2
Specification	L / M
Kinds	Female/Male

Anti-rotation foot suport



Model Number	ELK POR-ARFS-1
Specification	L/M/S
Kinds	Female/Male

Hallux Valgus Orthosis



Model Number	POR-HVO-1
Specification	L/M
Kinds	Female/Male

Hallux Valgus Orthosis



Model Number	ELK POR-HVO-2
Specification	L/M
Kinds	Female/Male

Ankle Protect



Model Number	ELK POR-AP-1
Specification	L/M/S
Kinds	Female/Male

Child Denis Shoes



Model Number	ELK POR-CDS-1
Specification	Shoes No. 8-16
Kinds	Boy / Girl

Child Denis Shoes



Model Number	ELK POR-CDS-2
Specification	Shoes No. 8-16
Kinds	Boy / Girl

Abdominal Corset



Model Number	ELK POR-AC-1
Specification	L/M/S
Kinds	Female/Male

Ankle foot orthopedic boots



Model Number	ELK POR-AOB-1
Specification	L/M/S
Kinds	Female/Male

Ankle foot orthopedic boots



Model Number	ELK POR-AOB-2
Specification	L/M/S
Kinds	Female/Male

Ankle foot orthopedic boots



Model Number	ELK POR-AOB-3
Specification	L/M/S
Kinds	Female/Male

Single belt Diabetic shoes



Model Number	ELK POR-SBDS-1
Specification	No. 39-43
Kinds	Female/Male

Yellow Diabetic Shoes



Model Number	ELK POR-YDS-1
Specification	No. 36-44
Kinds	Female/Male

Double Belt Diabetic Shoes



Model Number	ELK POR-DBDS-1
Specification	NO.39-43
Kinds	Female/Male

leather Diabetic shoes



Model Number	ELK POR-AP-1
Specification	NO.39-43
Kinds	Female/Male

Machine

Stepless speed regulating column lifting grinding machine



Voltage	220 V / 50HZ
Power	1000W
Speed	0-3000
Angle	0-180°
Size	Length560mm * wight380mm * height1330mm
Up and Down	400-1100mm
Characteristic	Manual or intelligent setting speed regulation Manual lifting , angle adjustment Double protection, low power consumption and low noise.

Oil-free silent double machine vacuum pump station



Power Supply	220 V / 50HZ
Power	550W*2=1100W
Maximum flow	single pump200L/min
extreme pressure	-91KPa
noise	<=60dB
Use environment	5-40℃
Volume of pumping station	600mm*475mm*1150mm
Air extraction channe	8mm*2

Duster



Voltage	220 V / 50HZ
power	550W
suction	1000Pa
Rated air flow	1150 cubic metre/hour
Interface diameter	100mm
Net weight	24kg

surface grinding machine



Voltage	220 V / 50HZ
power	750W
Sand belt	150*1219mm
motor speed	1400 Turn/branch

distant infrared-ray oven



Voltage	380 V / 50HZ
power	6000W
Dimensions	1600mm*1100mm*1220mm
Working dimensions	1250mm*1000mm(PTFE coating)
Heating tube	1000W*6
temperature range	0-300°C
working temperature	160-220°C
Fan	50W
floodlight	45W*2



Contact Us

Surgimed Liberia Limited



Street 19 & 20 St. Monrovia Liberi



+231 188 198 0951



info.lr@surgimedgrp.com

Vigor Medical Zambia Ltd.



Shop No. 10 5106 Luanshya Road
Lusaka Zambia



+260 763 947 565



info.zm@vigormedical.org

Surgimed Ghana Ltd.



20A, Old Fadama Road, Accra
Ghana



+233 241 584 310



info.gh@surgimedgrp.com

Surgimed Uganda Ltd.



Shop No.1, Plot No.2, Bombo
Road, Kampala Uganda



+256 704 010 522, +256 707 400 400



info.ug@surgimedgrp.com

Vigor Medical Kenya Ltd.



A-1 Premier Industrial Park,
Babadogo, Nairobi Kenya,



+254 784 314 341, +254 785 925 733



info.ke@vigormedical.org

Golden Gate Holdings



1438 Tamani Arts Office
Building, Business Bay, Dubai
United Arab Emirates, P.O. 691403



+971 567-797-033

