

ELK-35/50/75/100/120 LT VERTICAL PRESSURE STEAM STERILIZER



ELK35/50/75/100/120 LT Vertical Pressure Steam Sterilizer

Main Technical Data

Valid volume	23L
Working pressure	0.22MPa
Working temperature	134°C
Heat average	≤ ±1°C
Timer	0 – 60 min
Pressure-temperature controller	109 – 134°C, 0.04 – 0.22Mpa
Power	2.5kW/220v, 50Hz
Overall dimension	480mm x 460mm x 850 mm (D x W x H)
Transport dimension	570 mm x 550 mm x 970 mm (L x W x H)
Weight	gross: 56kg/76/96/116/136 (Size Based)