

MOBILE EXAMINATION LIGHT



Model No. ELK-849EM

For use in consultation rooms, during delivery or in other medical services.

Goose neck mobile in LED examination lamp

Diameter of Lamp head: 89mm

Color Temperature: 4500-5500K lux

Color Rendering Index (CRI): 90 Ra

Control System: Digital Pressing Control on The Lamp head

Light Field Diameter: 12cm at D=50cm 18cm at D=100cm

Extension of Arm: 760cm Length Goose Neck (Q7), Articulated Spring

Light Source: 7 LED

LED Average Life: 50,000 Hours

Power Consumption: 10W

Input Voltage: 110V-265V 50/60Hz

Standard: Mobile Type with Detachable Metal Base

Adjustable height: up to 130cm

Castors: 5 castors

Optional Color: White

Light Spot: Variable Only by Changing the Distance of Light

Lamp head: Integrated Handle and Digital Control

Transport Package: Standard Export Carton Package

Specification: 79*52*14cm (L*W*H)/ctn