

HEMOGLOBIN METER



Model No. ELK-HL23M Hemoglobin Meter

Use:

An apparatus for measuring the amount of hemoglobin in a fluid by comparing it with a solution of known strength and of normal colour

Features

- World's leading technology
- Siphon method to ensure quantitative sample
- Fast hemoglobin and HCT test
- Accurate and reliable result
- Suitable for various workplaces
- Bluetooth optional, Alarm Reminder
- Samples: Capillary Whole blood, venous or arterial
- Measuring range 0 to 25.6 g / dL (0 to 256 g / L, 0 to 15.9 mmol / L)
- Results In 5 seconds
- Sample volume ~ 10 µL
- Batteries Included
- Storage Temperature Analyzer: 0-50 ° C (32-122 ° F)
- Microcuvettes: unopened vial 15-30 ° C (59-86 ° F);
- Power Supply AC Adapter or 4 AA-size Batteries
- Printer Interface
- Quality Control "Auto - test integrated, optional control solutions

Gluko Doctor



*For Accurate And Reliable Results
Trust Your
Gluko Doctor*