

# ***NEONATAL BED (CRIBETTE)***



# Neonatal Bed (Cribette)

## SPECIFICATIONS:

- ICU model with 4 sides openings for 360-degree patient access
- High resistance to corrosion (tropical environment)
- Top Deck for sleeping with mattress included.
- Bottom Deck for accessories storage.
- Bed overall dimensions: 105.5 cm(L)\* 62.5 cm(W)\* 110 cm(H)
- Easy access ICU Safety Gate System
- With 4 Casters (two braking; ball bearing design)
- Solid pan sleeping system (non-helical design)
- Configured with two IV poles in stainless
- **Material:** Steel framework, top class ABS & PP plastic, non-toxic, odourless and environmentally friendly.
- Rev. trendelenburg: 0-12°
- Net weight: 32KG
- Safe load capacity: 10KG

## FUNCTIONS:

- Rev. trendelenburg positions operated by gas spring lever helps improve baby spit-up.
- 360° transparent side rails helps parents and caregivers observe the baby easily and create a safe area for the baby. Flip to open the side rails, which is convenient for caregivers. Additional small side rails inside protects the baby and keep safer.

**\*Delivered with 1x mattress with cover and all necessary accessories for it's good operation**

**\*On site installations and maintenance services can be provided.**

