

IV STAND



Model No. ELK-IV

Use:

A mobile, height adjustable stand for fluid administration.

Products Name	IV Stand
Dimensions	Ø480*H1200-2500mm
Brand	ELEKTRO
Model	JDASY152W2
Hook	4
Castor	5
Castor with break	2



Features

- Hook and pole material type: medical grade stainless steel with four hooks on top
- Load per hook; 2Kg
- Base material type: Plastic
- Five maintenance free Swivel casters with black rubber wheel; two with cross brake
- Telescopic adjustable height
- Pole sleeve welded securely to the base for reinforced stability
- Tube diameter 25mm
- 16mm stainless steel tube AISI 304

