

WARD SCREEN



Model No. ELK-226DWS Ward Screen

Use:

Foldable screen to position between hospital beds or to separate a space to create privacy for patients

Features

- Structure is made from non-magnetic stainless-steel material
- All Folding Panels are 60 wide. PVC Flame Retardant & Washable Curtains for Each Section.
- Dimensions (w x d x h), (2.40 x 0.05 x 1.70)

Specifications

- Total Wheels: 8
- Antistatic wheels

Technical Specification:

- Mobile Four section bed screen
- Comprising a metal tube frame mounted on 8 castors.
- Casters positioned in a broad stance for stability of the frame.
- Frame suspends a curtain material for privacy.
- Each side of the frame has a hinged section that can be angled for privacy. Curtains of white plastic material Dimensions (wxdxh), (2.40 x 0.05 x 1.70)
- Frame stainless steel
- Fireproof and washable plastic fabric

System Configuration Accessories, Spares

- 3x spare of castors
- All standard accessories included