

ELK 9000I VENTILATOR



ELK9000i Ventilator

- Super clear 18.5-inch color touch screen, German and English operation and display.
- Five wave four loops display, making doctors' understanding patients in a more intuitive way.
- Detachable suction and expiratory valve for easy cleaning and disinfection.
- Duration for the backup battery for 6 hours, avoiding power interruption during operation.
- With 4G and WiFi modules for more convenient use.



Parameter setting range	Monitoring parameters
Tidal volume: 2ml~4000ml	Tidal volume: 0-6000mL
Frequency: 1~100 times/minute	Frequency: 0~200 times/minute
Inspiration and Expiration ratio: 10:1~1:10	Inspiration and Expiration ratio: 150:1~1:150
Psupp: 0-80cm H ₂ O	Oxygen concentration: 18%~100%
Pinsp: 5-80cm H ₂ O	Leak% (Leakage rate: 0%~100%)
PEEP: OFF,1-50cmH ₂ O	PEEP: 0cmH ₂ O
Oxygen concentration: 21-100vol%	Minute ventilation: 0-100L/min

Invasive ventilation: VCV, PCV, VSIMV, PSIMV, CPAP/PSV, V+SIMV, BPAP, ASV, APRV, VS, PRVC

Non-invasive ventilation: PCV, SIMV-P, CPAP/PSV, APRV, BPAP, nCPAP, S/T