



## **ELEKTRO GENESIS**

Cherish your life, Cherish Your Health!



## **PATIENT CARE AND MEDICAL FURNITURE**



# ***PATIENT CARE AND MEDICAL FURNITURE***



## Patient Care and Medical Furniture

### CONTENTS:

01	ELK-EB Hospital Bed	3 - 10
02	ELK-EC Hospital Bed	11 - 22
03	ELK-OB Hospital Bed	23 - 24
04	ELK-CB Hospital Bed	25 - 28
05	Hospital Bed Parts	29 - 30
06	ELK-ET Hospital Bed	31 - 32
07	ELK-FS Stretcher	33 - 37
08	ABS Bedside Cabinet	38
09	Overbed Table	39
10	ABS Medical Trolley	40 - 48
11	Stainless Steel Products	49 - 58
12	Surgical Table	59 - 65
13	Electric dialysis chair	66 - 67
14	Medical Cart Series	68 - 73
15	Steam Sterilizer	74 - 75

## Patient Care and Medical Furniture

### ELK-EB-001-1 Electric Bed

- Detachable ABS head & foot board
- Tuck away big ABS railing
- Perforated powder coated steel bed platform
- Bed frame & surface made of cold-rolled steel tubes
- Patient & Nurse control panel built-in both railing
- Central locking 6" castors
- Angle degree indicator
- CPR quick release (manual & electrical)



### Feature display



### Technical Parameters

External Length	2150mm
External Width	1050mm
Bed Height Max.	800mm
Bed Height Min.	370mm
Mattress Platform	1950×900mm
Safe Working Load	250Kg

### Customized Products For You

#### Casters

No casters	<input type="checkbox"/>
Silent casters	<input type="checkbox"/>
Cover casters	<input type="checkbox"/>
CC casters Middle	<input type="checkbox"/>
CC casters Side	<input type="checkbox"/>

#### Bed slats

ABS Bed slats	<input type="checkbox"/>
Strips Bed slats	<input type="checkbox"/>
Mold Bed slats	<input type="checkbox"/>

#### Railing

Aluminum Alloy	<input type="checkbox"/>
4C <input type="checkbox"/> 5C <input type="checkbox"/> 6C <input type="checkbox"/>	
Stainless Steel	<input type="checkbox"/>
4C <input type="checkbox"/> 5C <input type="checkbox"/> 6C <input type="checkbox"/>	
ABS Railing	<input type="checkbox"/>
4 Tablets	<input type="checkbox"/>
2 Big	<input type="checkbox"/>
New Molds	<input type="checkbox"/>
Color	<input type="checkbox"/>

#### More features

Electric CPR	<input type="checkbox"/>
Manual CPR	<input type="checkbox"/>
Spare batteries	<input type="checkbox"/>
Angular display	<input type="checkbox"/>
Weighing System	<input type="checkbox"/>
X-ray Plates	<input type="checkbox"/>
Nurse controller	<input type="checkbox"/>
Railing controller	<input type="checkbox"/>
Soft connections	<input type="checkbox"/>
Electric with manual	<input type="checkbox"/>

## Patient Care and Medical Furniture

### ELK-EC-001 Electric Chair Hospital Bed



- Detachable ABS head & foot board
- Tuck away big ABS railing
- Perforated powder coated steel bed platform
- Bed frame & surface made of cold-rolled steel tubes
- Patient & Nurse control panel built-in both railing
- Central locking 6" castors
- Angle degree indicator
- CPR quick release (manual & electrical)
- One button chair position
- Bed platform can do X-ray



## Patient Care and Medical Furniture

### ELK-EC-001-2 Electric Chair Hospital Bed



- Detachable ABS head & foot board with nurse control
- Powder coated steel bed frame
- Bed framework made from mild steel ERW rectangular tubes.
- ABS siderails with build-in controller
- Central control castor with 150mm
- Drainage hook & S.S.I.V. pole
- X-ray, foot control: (optional)
- Electronic automatic weighing system (optional).

#### 5 Functions:

- Back-rest section: 80°+5"
- Leg section: 40+5"
- Height adjustment 50-75cm
- Trendelenburg 0-12°
- Anti-trendelenburg 0-12°



## Patient Care and Medical Furniture

### ELK-EC-002 Electric Chair Hospital Bed



- Detachable ABS head & foot board
- Collapsed type aluminum alloy 5-stuff railing
- Perforated powder coated steel bed platform
- Central locking 5"castors
- Chair position

### ELK-EC-003 Electric Chair Hospital Bed



- Detachable ABS head & foot board
- Tuck away ABS railing
- Perforated powder coated steel bed platform
- Control panel built-in railing
- Nurse control panel built-in foot board
- Central locking 5"castors

## Patient Care and Medical Furniture

### ELK-EC-004 Electric Chair Hospital Bed



- Detachable ABS head & foot board
- Collapsed type aluminum alloy 6-stuff railing
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake



### ELK-EC-005 Electric Chair Hospital Bed



- Detachable ABS head & foot board
- Tuck away ABS railing
- Perforated powder coated steel bed platform
- Platform with multi-board
- Control panel built-in railing
- Central locking 5" castors



### ELK-EC-005 Electric Chair Hospital Bed



- Detachable ABS head & foot board
- Tuck away ABS railing
- Wire mesh powder coated steel bed platform
- Central locking 5" castors
- Turn-over position

### ELK-EC-006 Electric Chair Hospital Bed



- Detachable ABS head & foot board
- Tuck away big ABS railing
- Powder coated steel bed platform
- Powerful actuator system.
- Central locking 5" castors
- Nurse control below foot board



## Patient Care and Medical Furniture

### ELK-EC-007 Electric Chair Hospital Bed



- Detachable ABS head & foot board
- Tuck away ABS railing
- ABS bed platform
- Powerful actuator system.
- Central locking 5" castors



### ELK-EC-008-2 Electric Chair Hospital Bed



- Detachable ABS head & foot board
- Tuck away ABS railing
- Perforated powder coated steel bed platform with soft connection
- Powerful actuator system.
- Central locking 5" castors
- Electric with two function adjustment by crank



### ELK-EC-009 Electric Chair Hospital Bed



- Detachable ABS head & foot board
- Collapsed type aluminum alloy 6-stuff railing
- Perforated powder coated steel bed platform
- Powerful actuator system.
- Central locking 5" castors



### ELK-EC-010 Electric Chair Hospital Bed



- Painted wooden head & foot board & side board
- Collapsed type aluminum alloy 6-stuff railing
- ABS bed platform
- Central locking 5" castors



## Patient Care and Medical Furniture

### ELK-CB-001 Four Crank Five Function

#### Feature display



#### Exquisite

Simple structure, exquisite workmanship

#### Simple

Crank bed structure is relatively simple

#### Practical

Simple and practical medical products

Crank bed is a relatively simple hospital bed. Its structure is simple, but the workmanship is not simple. It is usually composed of head & foot board (optional), bed frame (optional), bed surface (optional), guardrail (optional), crank system (optional), casters (optional) and other major components.

Crank system is the simplest structure, but can meet the basic needs of patient. It can also be a complex structure to deal with changes in various special situations.



- Detachable ABS head & foot board
- Tuck away 2-piece big ABS railing with dining table
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake

#### Technical Parameters

External Length	2080mm
External Width	970mm
Bed Height Max.	800mm
Bed Height Min.	370mm
Mattress Platform	1950*830mm
Safe Working Load	250Kg


### ELK-CB-001-2 Four Crank Five Function



- Detachable ABS head & foot board
- Collapsed type aluminum alloy 6-stuff railing
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake



### ELK-CB-001-3 Four Crank Five Function



- Detachable ABS head & foot board
- Tuck away 4-piece ABS railing
- Perforated powder coated steel bed platform
- Central locking 5" castors



### ELK-CB-002 Four Crank Five Function



- Detachable ABS head & foot board
- Collapsed type aluminum alloy 6-stuff railing
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake



### ELK-CB-002-1 Three Crank Hospital Bed



- Detachable ABS head & foot board
- Tuck away ABS railing
- Perforated powder coated steel bed platform with soft connection
- Central locking 5" castors



### ELK-CB-002-2 Three Crank Hospital Bed



- Detachable ABS head & foot board
- Collapsed type aluminum alloy 5-stuff railing
- Perforated powder coated steel bed platform
- Noiseless 5" castors & 4 feet
- Oxygen bottle holder



## Patient Care and Medical Furniture

### ELK-CB-002-3 Three Crank Hospital Bed



- Detachable ABS head & foot board
- Collapsed type aluminum alloy 5-stuff railing with dining table
- Perforated powder coated steel bed platform
- Noiseless 5" castors & 4 feet



### ELK-CB-002-4 Three Crank Hospital Bed



- Detachable ABS head & foot board
- Collapsed type aluminum alloy 6-stuff railing
- Perforated powder coated steel bed platform
- Noiseless 5" castors & 4 feet



## Patient Care and Medical Furniture

### ELK-CB-002-5 Three Crank Hospital Bed



- Detachable ABS head & foot board
- Tuck away 2-piece big ABS railing
- Perforated powder coated steel bed platform with soft connection
- Central locking 5" castors




### ELK-CB-002-6 Three Crank Hospital Bed



- Detachable stainless steel head & foot board
- Collapsed type aluminum alloy 5-stuff railing
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake



### ELK-CB-003 Two Crank Hospital Bed




- Detachable ABS head & foot board
- Collapsed type aluminum alloy 4-stuff railing with dining table
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake




## Patient Care and Medical Furniture

### ELK-CB-003-1 Two Crank Hospital Bed




- Detachable stainless steel head & foot board
- Collapsed type stainless steel 6-stuff railing
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake

### ELK-CB-003-2 Two Crank Hospital Bed




- Detachable aluminum alloy head & foot board
- Collapsed type aluminum alloy 4-stuff railing with dining table
- Perforated powder coated steel bed platform
- Central locking 5" castors

### ELK-CB-003-3 Two Crank Hospital Bed



- Detachable ABS head & foot board
- Collapsed type aluminum alloy 6-stuff railing with dining table
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake

### ELK-CB-003-4 Two Crank Hospital Bed



- Detachable aluminum alloy head & foot board
- Collapsed type aluminum alloy 4-stuff railing with dining table
- Perforated powder coated steel bed platform
- Central locking 5" castors




### ELK-CB-003-5 Two Crank Hospital Bed (Stainless Steel)



- Detachable stainless steel head & foot board
- Collapsed type aluminum alloy 4-stuff railing with dining table
- Painted wooden batten bed platform, stainless steel bed frame
- Noiseless 5" castors with brake



### ELK-CB-003-6 Two Crank Hospital Bed



- Detachable powder coated steel and wooden head & foot board
- Collapsed type aluminum alloy 6-stuff railing
- Strip powder coated steel bed platform
- Noiseless 5" castors with brake



### ELK-CB-003-7 Three Crank Hospital Bed (Stainless Steel)



- Detachable Stainless head & foot board
- Stainless steel side railing
- Stainless steel bed platform
- Bed framework made from stainless steel
- Noiseless castor dia. 125mm with brake
- With drainage bag holders & telescopic SS 2pcs hook IV
- Three functions by hand crank 4 inch protective with roller bumper in every corner

Backresting angle:  $80^{\circ} \pm 5^{\circ}$

Knoeresting angle:  $4025^{\circ}$

Height adjustment: 500-750mm

- Dimension: 2000\*800\*500-750mm



### ELK-CB-004 One Crank Hospital Bed



- Detachable aluminum alloy head & foot board
- Collapsed type aluminum alloy 4-stuff railing with dining table
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake



Wire Mesh

### ELK-CB-004-1 One Crank Hospital Bed



- Detachable ABS head & foot board
- Perforated powder coated steel bed platform



### ELK-CB-004-2 One Crank Hospital Bed



- Detachable aluminum alloy head & foot board
- Collapsed type aluminum alloy 4-stuff railing with dining table
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake



### ELK-CB-004-3 One Crank Hospital Bed



- Detachable ABS head & foot board
- Collapsed type aluminum alloy 6-stuff railing
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake



## Patient Care and Medical Furniture

### ELK-CB-004-4 One Crank Hospital Bed



- Stainless steel head & foot board
- Insert stainless steel guardrail
- Perforated powder coated steel bed platform



### ELK-CB-004-5 One Crank Hospital Bed



- Stainless steel head & foot board
- Wire mesh powder coated steel bed platform



### ELK-CB-005 Flat Hospital Bed

- Detachable stainless steel head & foot board
- Collapsed type stainless steel 5-stuff railing
- Perforated powder coated steel bed platform



### ELK-CB-005-1 Flat Hospital Bed

- Powder coated steel head & foot board
- Perforated powder coated steel bed platform



## Patient Care and Medical Furniture

### ELK-CB-005-2 Flat Hospital Bed

- Detachable ABS head & foot board
- Bar powder coated steel bed platform



### ELK-CB-005-3 Flat Hospital Bed

- Stainless steel / powder coated steel head & foot board
- Bar powder coated steel bed platform




### ELK-CB-005-4 Flat Hospital Bed

- Stainless steel / powder coated steel head & foot board
- Bar powder coated steel bed platform



## Patient Care and Medical Furniture


### ELK-OB-003-1 Orthopaedic Traction Bed



- Detachable aluminum alloy head & foot board
- Collapsed type aluminum alloy 4-stuff railing with flexible dining table
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake
- Aluminum alloy portal frame and direction frame



### ELK-OB-003-2 Orthopaedic Traction Bed




- Collapsed type aluminum alloy 6-stuff railing
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake
- Stainless steel head & foot board & portal frame & direction frame



## Patient Care and Medical Furniture


### ELK-OB-004-1 Orthopaedic Traction Bed



- Detachable ABS head & foot board
- Collapsed type aluminum alloy 4-stuff railing with dining table
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake
- Stainless steel double arm lifting pole



### ELK-OB-004-2 Orthopaedic Traction Bed



- Detachable ABS head & foot board
- Collapsed type aluminum alloy 4-stuff railing with dining table
- Perforated powder coated steel bed platform
- Noiseless 5" castors with brake
- Stainless steel double arm lifting pole




### ELK-CB-004-C1 Single Crank Children Bed



- Powder coated steel head & foot boards with laminated panels.
- Full length powder coated side rail, drop type
- Covered castors dia75mm with brake
- Dimension: 1800\*800\*580-1100mm
- With mattress (optional)
- With backrest lifting
- Noiseless 4" castors with brake



### ELK-CB-004-C2 Height Adjustable Children Bed



- Detachable aluminum alloy head & foot board
- Collapsed type aluminum alloy 4-stuff railing
- Perforated powder coated steel bed platform
- Noiseless 5 castors with brake



**ELK-CB-005-C1 Baby Crib**



Size: 890\*510\*780-1000mm

**ELK-CB-005-C2 Baby Crib**



Size: 875\*525\*650\*900mm

**ELK-CB-005-C3 Baby Crib**



Size: 875\*525\*650-900mm

**ELK-CB-005-C4 Baby Crib**



Size: 1600\*700\*550mm

**ELK-CB-005-C5 Baby Crib**



Size: 1600\*640\*1040mm

**ELK-CB-005-C6 Baby Crib**



Size: 1100\*640\*1040mm

## Patient Care and Medical Furniture

### ELK-CB-005-C1 Baby Crib



#### SPECIFICATION:

- Baby basin 875x 525 mm
- Tilt angle 0-15
- Hi-Low 650-900mm

### ELK-CB-005-C1A Baby Crib



#### SPECIFICATION:

- Baby basin 800 x 420 mm.
- Tilt angle 0-15
- Hi-Low 855-1045mm
- Hand rail angle 0-90

## Patient Care and Medical Furniture

### ELK-CB-004-C3 Baby Crib



- Steel structure. Plastic parts adopt high quality ABS & PP
- Aluminum alloy side rail ABS head & foot board; Spray iron bed board
- Overall size 1880mm 900mmx1210mm
- Bed board dimension 1770mm\*720mm
- Height from ground to bed board 585mm
- Castor: 5" double side castor
- Net Weight: 10KGS
- Load capacity: 200KGS
- Back-rest rising
- Tilt angle 0-75° Foot-rest rising
- Tilt angle 0-40°
- Side rail With three gears adjustable
- First gear Height from side rail to bed board) 80mm
- Second gear (Height from side rail to bed board) 310mm
- Third gear (Height from side rail to bed board) 500mm

## Patient Care and Medical Furniture

**ELK-FB-001**



Flat Bed mattress

**ELK-FB-003**



Double Crank Bed Mattress

**ELK-FB-005**



Air Mattress

**ELK-FB-002**



Single Crank Bed Mattress

**ELK-FB-004**



Orthopaedic Bed Mattress

**ELK-FB-006**



Detachable Air Mattress



Army Green



Blue



Grey

## Patient Care and Medical Furniture

### Head and Foot Board



### Soft Connection



### Guardrail



## Patient Care and Medical Furniture

### ELK-ET-004-01 Electric Exam Table

- Main technical parameters
- Size: 1800mm x 580mm
- Overall lifting adjustment: 500mm-700mm
- Back panel folding angle: 0°-60°



### ELK-ET-004-02 Electric Exam Table

#### Technique Data

- One electric motor for table Hi-Low function
- Height adjustable: 47cm - 92cm
- Angle of headrest: -20° to +75°
- Table length and width: 195cm x 66cm
- Max loading weight: 225KGS
- Certificate: ISO13485/CE/FDA/93 42 EEC
- Wheels: 75mm castors for safety and mobility
- Packing size: 190cm x 72cm x 54cm



### ELK-ET-004-03 Electric Exam Table

#### Description

- Two electric motors for table Hi-Low function
- Height adjustable: 64.5cm - 106cm
- Angle of backrest: 0° to +80°
- Angle of seat rest: 0° to +20°
- Table length and width: 185cm x 60cm
- Max loading weight: 225KGS
- Certificate: ISO13485/CE/FDA/93 42 EEC
- Wheels: 75mm castors for safety and mobility
- Packing size: 140cm x 74cm x 69cm



## Patient Care and Medical Furniture

### ELK-CB-004-Z1



**Exam Table (Stainless Steel)**

Size: 1900\*600\*680mm

### ELK-CB-004-Z2



**Exam Table (Stainless Steel)**

Size: 1900\*600\*680mm

### ELK-CB-004-Z3



**Exam Table (Stainless Steel)**

Size: 1900\*600\*680mm

### ELK-CB-004-Z4



**Parturition Bed**

Size: 1870\*600\*680mm

### ELK-CB-005-Z1



**Massage Table**

Size: 1900\*600\*680mm

### ELK-CB-005-Z2



**Exam Table (Powder Coated Steel)**

Size: 1900\*600\*680mm

### ELK-CB-005-Z3



**Exam Table (Powder Coated Steel)**

Size: 1900\*600\*680mm

### ELK-CB-005-Z4



**Exam Table (Stainless Steel)**

Size: 1900\*600\*680mm

## Patient Care and Medical Furniture

### ELK-FS-001 Foldable Stretcher

Full size (L\*W\*H): 215\*54\*12cm

Folded size (L\*W\*H): 58\*21\*14cm



### ELK-FS-002 Foldable Stretcher

Full size (L\*W\*H): 208\*55\*13cm

Folded size (L\*W\*H): 104\*17\*9cm



### ELK-FS-003 Foldable Stretcher

Full size (L\*W\*H): 185\*50\*18cm

Folded size (L\*W\*H): 96\*54\*16cm



## Patient Care and Medical Furniture

### ELK-FS-004 Spine Board Stretcher



Full size (L\*W\*H): 185\*47\*6.5cm

### ELK-FS-005 Foldable Stretcher



Full size(L\*W\*H): 220\*55\*15cm

Folded size(L\*W\*H): 95\*20\*9cm

### ELK-FS-006 Scoop Stretcher



Full size(L\*W\*H): 214\*42\*7cm

Folded size(L\*W\*H): 167\*42\*7cm

### ELK-FS-007 Stair Stretcher



Full size (L\*W\*H): 69\*51\*90cm

Folded size (L\*W\*H): 18\*51\*90cm

### ELK-FS-008 Stair Stretcher



Full size (L\*W\*H): 69\*51\*90cm

Folded size (L\*W\*H): 18\*51\*90cm

### ELK-FS-009 Stair Stretcher



Full size (L\*W\*H): 69\*51\*90cm

Folded size (L\*W\*H): 18\*51\*90cm

## Patient Care and Medical Furniture

### ELK-FS-010 Ambulance Stretcher

- Load Bearing: 159kg
- Product Size (L\*W\*H):1900\*550\*500mm
- Max. angle of the back: 70°
- N.W.: 25.5KG



### ELK-FS-011 Ambulance Stretcher

- Load Bearing: 159kg
- Product Size (L\*W\*H):1900\*550\*500mm
- Max. angle of the back: 70°
- N.W.: 29KG



### ELK-FS-012 Ambulance Stretcher

- Product Size(L\*W\*H)(high position):  
1900\*550\*920mm
- Product Size(L\*W\*H)(low position):  
1900\*550\*260mm



## Patient Care and Medical Furniture

### ELK-FS-013 Ambulance Stretcher

- Product Size(L\*W\*H)(high position):  
1900\*550\*920mm
- Product Size(L\*W\*H)(low position):  
1900\*550\*250mm



### ELK-FS-014 Ambulance Stretcher

- Product Size(L\*W\*H)(high position):  
1900\*550\*920mm
- Product Size(L\*W\*H)(low position):  
1900\*550\*260mm



### ELK-FS-015 Stretcher Trolley for Ambulance

- Load Bearing: 159kg
- Product Size (L\*W\*H):1900\*550\*860mm
- Max. angle of the back: 85°
- N.W.: 45KG



## Patient Care and Medical Furniture

### ELK-FS-016 Stretcher Trolley for Ambulance

- Load Bearing: 159kg
- Product Size (L\*W\*H)high position): 1900\*550\*880mm
- Product Size (L\*W\*H)low position): 1900\*550\*540mm
- Max, angle of the back: 70°
- N.W.: 35KG



### ELK-FS-017 Stretcher Trolley for Ambulance

- Load Bearing: 159kg
- Product Size (L\*W\*H)high position): 1900\*550\*880mm
- Product Size (L\*W\*H)low position): 1900\*550\*540mm
- Max, angle of the back: 70°
- N.W.: 35KG



### ELK-FS-018 Stretcher Trolley for Ambulance

- Load Bearing: 159kg
- Product Size (L\*W\*H)high position): 1900\*550\*880mm
- Product Size (L\*W\*H)low position): 1900\*550\*540mm
- Max, angle of the back: 70°
- N.W.: 35KG



## Patient Care and Medical Furniture

### ELK-BC-001

ABS Bedside Cabinet

Size 480x480x760mm



### ELK-BC-002

ABS Bedside Cabinet

Size 480x480x870mm



### ELK-BC-003

ABS Bedside Cabinet

Size 480x480x760mm



### ELK-BC-004

Stainless Steel Cabinet

Size 480x40x780mm



### ELK-BC-005

Metal Cabinet

Size 480x400x790mm



### ELK-BC-006

Metal Cabinet

Size 470x440x780mm



### ELK-BC-007

Stainless Steel Cabinet

Size 480x480x760mm



### ELK-BC-008

Metal Cabinet

Size 480x480x760mm



### ELK-BC-009

Metal Cabinet

Size 480x480x870mm



## Patient Care and Medical Furniture

### ELK-CKT-032 Overbed Table



### ELK-CKT-033 Overbed Table

Size: 800\*410\*730-1130



ABS Table Top with Cup Socket  
Aluminum alloy columns  
Steel base with surface treatment  
Air Spring for Height Adjustment

### ELK-CKT-034 Overbed Table

Size: 800\*410-730-1130mm



### ELK-CKT-035 Overbed Table

Size: 800\*410-730-1130mm



### ELK-CKT-036 Overbed Table

Size: 800\*410-730-1130mm



Wooden table top with edging  
Stainless steel frame or Iron frame  
Air Spring for Height Adjustment

### ELK-CKT-037 Overbed Table

Size: 800\*410\*730-1130mm



### ELK-CKT-038 Overbed Table

Size: 800-410-730-1130mm



### ELK-CKT-039 Electric Patient Lift



**ELK-MTS-001-a1**  
**ABS Medicine Trolley**

Size: 730x480x900mm



**ELK-MTS-001-b1**  
**ABS Medicine Trolley**

Size: 730x480x900mm



**ELK-MTS-001-b2**  
**ABS Medicine Trolley**

Size: 730x480x900mm



**ELK-MTS-003**  
**Treatment Trolley**



**ELK-MTS-002-a1**  
**ABS Medicine Trolley**

Size: 539x366x900mm



**ELK-MTS-002-a2**  
**ABS Medicine Trolley**

Size: 539x366x900mm



**ELK-MTS-002-b1**  
**ABS Medicine Trolley**

Size: 539x366x900mm



**ELK-MTS-004**  
**Treatment Trolley**



**ELK-MTS-005**  
**Treatment Trolley**



**ELK-MTS-006**  
**Treatment Trolley**



**ELK-MTS-007**  
**Treatment Trolley**



**ELK-MTS-008**  
**Treatment Trolley**



**ELK-MTS-009**  
**Treatment Trolley**



**ELK-MTS-010**  
**Treatment Trolley**



**ELK-MTS-011**  
**Treatment Trolley**



**ELK-MTS-013**  
**Treatment Trolley**



## Patient Care and Medical Furniture

### ELK-MTS-012 Emergency Trolley

Size:750x475x930mm



#### Technical Configuration:

Defibrillator shelf	1PC	IV pole	1PC	Central Lock	1PC
Sliding Side Shelf	1PC	CPR Board	1PC	File Box	1PC
Sharp Container	1PC	Oxygen Bottle Holder	1PC		
Dust Baskets	2PCS	Power Outlet	1PC		

## Patient Care and Medical Furniture

### ELK-MTS-014 Emergency Trolley

Size 750x475x930mm



### ELK-MTS-015 Emergency Trolley

Size 750x475x930mm



### ELK-MTS-016 Emergency Trolley

Size: 625x475x920mm



### ELK-MTS-017 Emergency Trolley

Size 750x475x930mm



### ELK-MTS-018 Medicine Trolley



### ELK-MTS-019 Medicine Trolley



### ELK-MTS-020 Medicine Trolley



### ELK-MTS-021 Medicine Trolley



## Patient Care and Medical Furniture

### **ELK-MTS-022** **Anesthesia Trolley**

Size 750x475x930mm



### **ELK-MTS-024** **Anesthesia Trolley**

Size 750x475x930mm



### **ELK-MTS-026** **Anesthesia Trolley**

Size: 625x475x930mm



### **ELK-MTS-027** **Anesthesia Trolley**

Size: 625x475x930mm



### **ELK-MTS-023** **Anesthesia Trolley**

Size: 625x475x930mm



### **ELK-MTS-025-1** **Record Holder Trolley**



### **ELK-MTS-025-2** **Record Holder Trolley**



### **ELK-MTS-028** **Disease Record**



## Patient Care and Medical Furniture

### **ELK-ITC-001 Instrument Trolley**

Size:670×450×900mm



### **ELK-ITC-002 Instrument Trolley**

Size:670×450×900mm



### **ELK-ITC-003 Instrument Trolley**

Size:670×450×900mm



### **ELK-ITC-004 Fan Instrument Trolley**

Size:1400×460×900mm



### **ELK-ITC-005 Instrument Trolley**

Size:670×450×900mm



### **ELK-ITC-006 Instrument Trolley**

Size:670×450×900mm



### **ELK-ITC-007 Instrument Trolley**

Size:670×450×900mm



### **ELK-ITC-008 Instrument Trolley**

Size:670×450×900mm



**ELK-ITC-009**  
**Emergency Trolley**

Size: 700x400x850mm



**ELK-ITC-011**  
**Emergency Trolley**

Size: 700x400x850mm



**ELK-ITC-013**  
**Instrument Trolley**

Size: 700x400x850mm



**ELK-ITC-015**  
**Emergency Trolley**

Size: 700x400x850mm



**ELK-ITC-010**  
**Anaesthetic Trolley**

Size: 600x400x850mm



**ELK-ITC-012**  
**Emergency Trolley**

Size: 600x400x850mm



**ELK-ITC-014**  
**Anaesthetic Trolley**

Size: 600x400x850mm



**ELK-ITC-016**  
**Instrument Trolley**

Size: 670x450x900mm



## Patient Care and Medical Furniture

### ELK-ITC-017 Medical Record Holder

Size: 340x400x900mm



### ELK-ITC-018 Medical Record Holder

Size: 630x430x900mm



### ELK-ITC-019 Infusion Treatment Trolley

Size: 670x450x1500mm



### ELK-ITC-020 Infusion Treatment Trolley

Size: 670x450x1500mm



### ELK-ITC-021 Mayo Trolley

Size: 450x350x1050mm



### ELK-ITC-022 Mayo Trolley

Size: 600x420x1030mm



### ELK-ITC-023 Soak Forceps Bucket Trolley

Size: dia. 300mm



### ELK-ITC-024 Medicine-change Bowl

Size: dia. 160mm



## Patient Care and Medical Furniture

### ELK-ITC-025 Movable Waste Bucket

Size: 14L



### ELK-ITC-026 Dressing Trolley

Size: 450x450x860mm



### ELK-ITC-027 Dressing Trolley

Size: 800x600x840mm



### ELK-ITC-028 Dressing Trolley

Size: 800x600x840mm  
Size: 600x840mm



### ELK-ITC-029 Dressing Trolley

Size: 900x450x1050mm



### ELK-ITC-030 Dressing Trolley

Size: 900x450x1050mm



### ELK-ITC-031 Dressing Trolley

Size: 670x550x890mm



### ELK-ITC-032 Dressing Trolley

Size: 700x400x850mm



## Patient Care and Medical Furniture

### ELK-ITC-032 Dressing Trolley

Size: 750x500x950mm



### ELK-ITC-033 Dressing Trolley

Size: 450x450x860mm



### ELK-ITC-034 Dressing Trolley

Size: 1280x840x980mm



### ELK-ITC-036 Treatment Trolley

Size: 900x550x950mm



### ELK-ITC-037 Oxygen Cylinders Trolley

Size: 550x380x1330mm



### ELK-ITC-038 Oxygen Cylinders Trolley

Size: 550x380x1330mm



### ELK-ITC-035 Dressing Trolley

Size: 1350x600x800mm



### ELK-ITC-039 Medical Diagnostic Lights

Power: 60W



## Patient Care and Medical Furniture

### **ELK-ITC-040** **Stainless Steel Stand**

Size : 450\*450\*1400-2200




### **ELK-ITC-042** **Stainless Steel Stand**

Size : 450\*450\*1400-2200



### **ELK-ITC-044** **Ward Screen (S.S.)**

Size: 2000x1870mm



### **ELK-ITC-046** **Trash Can (S.S.)**

Size: dia. 220x570mm



### **ELK-ITC-041** **Stainless Steel Stand**

Size : 450\*450\*1400-2200



### **ELK-ITC-043** **Ceiling I.V. Pole**

Size : 450\*450\*1400-2200



### **ELK-ITC-045** **Ward Screen (S.S.)**

Size: 1500x1870mm



### **ELK-ITC-047** **Trash Can (S.S.)**

Size: 12L



## Patient Care and Medical Furniture

**ELK-ITL-001**  
**One Step Stool**

Size: 460\*410\*250mm



**ELK-ITL-002**  
**Two Step Stool**

Size: 460\*230/490\*280mm



**ELK-ITL-003**  
**Two Step Stool**

Size: 460\*410\*430mm



**ELK-ITL-004**  
**Two Step Stool**

Size: 405\*230/490\*275mm



**ELK-ITL-005**  
**Two Step Stool**

Size: 460\*230/490\*280mm



**ELK-ITL-006**  
**Luxury Transfusion Chair**

Size: 710\*800\*1020mm



**ELK-ITL-009**  
**Luxury Transfusion Chair**

Size: 710\*800\*1020mm



**ELK-ITL-010**  
**Luxury Transfusion Chair**

Size: 710\*800\*1020mm



**ELK-ITL-011**  
**Luxury Transfusion Chair**



**ELK-ITL-012**  
**Luxury Transfusion Chair**



**ELK-ITL-013**  
**Luxury Transfusion Chair**



**ELK-ITL-014**  
**Nursing Folding Chair**



## Patient Care and Medical Furniture

### **ELK-ITL-015** **Nursing Folding Chair**

Open Size: 600\*600\*920mm

Close Size: 1850\*600\*420mm



**ELK-ITL-016  
Dental Chair**

Size: 350x450-560mm



**ELK-ITL-017  
Nurse Chair**

Size: 400x500-660mm



**ELK-ITL-018  
Nurse Chair**

Size: 450x540-750mm



**ELK-ITL-019  
Nurse Chair**

Size: 400x500-660mm



**ELK-ITL-020  
Nurse Chair**

Size: 410x500-870mm



**ELK-ITL-021  
Nurse Chair**

Size: 400x500-660mm



**ELK-ITL-022  
Lab Chair**

Size: 350x435-575mm



**ELK-ITL-023  
Lab Chair**

Size: 330x230x450-590mm



**ELK-ITL-024  
Doctor Chair**

Size: 410x500-870mm



**ELK-ITL-025  
Doctor Chair**

Size: 410x500-870mm



**ELK-ITL-026  
Doctor Chair**

Size: 400x500-660mm



**ELK-ITL-027  
Training Chair**

Size: 560x650x750mm



**ELK-ITL-028  
Doctor Chair**

Size: 410x500x870mm



**ELK-ITL-029  
Doctor Chair**

Size: 400x500-660mm



**ELK-ITL-030  
Doctor Chair**

Size: 380x280x475mm



**ELK-ITL-031  
Doctor Chair**

Size: 380x280x475mm



## Patient Care and Medical Furniture

**ELK-ITL-038**  
Wheelchair

Size: 670\*1050\*510/900



**ELK-ITL-039**  
Commode Chair

Size: 560\*66/470-920mm



**ELK-ITL-040**  
Commode Chair

Size: 710 63/460-1130mm



**ELK-ITL-041**  
Treat-waiting Chair

Size: 1800x740x800mm



**ELK-ITL-042**  
Treat-waiting Chair

Size: 1780x740x800mm



**ELK-ITL-043**  
Treat-waiting Chair

Size: 1800x740x800mm



**ELK-ITL-044**  
Treat-waiting Chair

Size: 1800x740x800mm



**ELK-ITL-045**  
Treat-waiting Chair

Size: 1800x710x800mm



**ELK-ITL-046**  
Treat-waiting Chair

Size: 1800x710x800mm



### ELK-ITP-001 Medical Square Tray (S.S.)

Size(cm): 23×17 2.8  
26×20×2.9  
40×30×3.3  
31×24×3.1  
45×33×3.4



### ELK-ITP-002 Bed Pan (S.S.)

Size: 1.2L



### ELK-ITP-003 Bed Pan (S.S.)

Size: 36×29×9cm



### ELK-ITP-004 Cream Cans (S.S.)

Size: 12cm/14cm/16cm



### ELK-ITP-005 Sterilized Bottles (S.S.)

Size: φ50×53mm  
φ50×83mm  
φ50×93mm  
φ50×103mm



### ELK-ITP-006 Measuring Cup (S.S.)

Size: 500ml, 1000ml



### ELK-ITP-007 Medical Tray (S.S.)

Size:  
30×20×2.8cm 35×25×2.8cm  
40×30×2.8cm 45×30×2.8cm  
45×35×2.8cm 50×35×2.8cm  
50×40×2.8cm 60×40×2.8cm



### ELK-ITP-008 Medical Tray with Lid (S.S.)

Size: 8", 9", 12"



### ELK-ITP-009 Deep Kidney Dish (S.S.)

Size:  
225×125×40mm  
210×110×35mm  
200×105×30mm



### ELK-ITP-010 Shallow Kidney Dish (S.S.)

Size:  
204×130×25mm  
190 116×25mm  
178×104×25mm



### ELK-ITP-011 Treatment Tray (S.S.)

Size:  
25×16×5cm  
30×20×5cm  
31×24×3.1cm  
40×30×3.3cm



### ELK-ITP-012 Test Tube Holder (S.S.)

Size:  
(φ18mm Hole)  
30 holes; 40 holes; 60 holes

Size:  
(φ18mm Hole)  
30 holes; 40 holes; 60 holes



## Patient Care and Medical Furniture

### ELK-ITM-001 Hickey Cabinet

Material: Stainless steel or powder coated steel  
Size: 900\*400\*1750mm



### ELK-ITM-002 Medical Cabinet

Material: stainless steel or powder coated steel  
Size: 900\*400\*1750mm



### ELK-ITM-003 Medical Cabinet

Material: Stainless steel or powder coated steel  
Size: 1800\*500/300\*1750mm



### ELK-ITM-004 Hickey Cabinet

Material: Stainless steel or powder coated steel  
Size: 900\*400\*1750mm



### ELK-ITM-005 Medical Cabinet

Material: Stainless steel or powder coated steel  
Size: 900\*400\*1750mm



### ELK-ITM-006 Gremfree Cabinet

Material: Stainless steel or powder coated steel  
Size: 900\*400\*1750mm



## Patient Care and Medical Furniture

### **ELK-ITM-007** **Six-gateway Cabinet**

Material: stainless steel or powder coated steel  
Size: 900\*400\*1750mm



### **ELK-ITM-007-1** **Four-gateway Cabinet**

Material: Powder coated steel  
Size: 900\*450\*2100mm



### **ELK-ITM-008** **Eighteenth-gateway Shoes Cabinet**

Material: Powder coated steel  
Size: 900\*450\*1750mm



### **ELK-ITM-009** **Filling Cabinet**

Material: Powder coated steel  
Size: 850\*400\*1850mm



### **ELK-ITM-010** **Filling Cabinet**

Material: Powder coated steel  
Size: 850\*400\*1850mm



### **ELK-ITM-011** **Filling Cabinet**

Material: Powder coated steel  
Size: 850\*400\*1850mm



## Patient Care and Medical Furniture

### ELK-ITM-012 Single-sided Medicine Shelf

Material: Powder coated steel  
Size: 800\*460\*2000mm



### ELK-ITM-013 Two-sided Medicine Shelf

Powder coated steel  
Size: 800\*910\*2000mm



### ELK-ITM-014 Single-sided Medicine Shelf

Powder coated steel  
Size: 800\*910\*2000mm



### ELK-ITM-015 Two-sided Medicine Shelf

Material: Powder coated steel  
Size: 1200\*460\*2040mm



### ELK-ITM-016 Two-sided Medicine Shelf

Powder coated steel  
Size: 800\*800\*2000mm



### ELK-ITM-017 Adjustable Component Medicine Shelf

Powder coated steel  
Size: 1030\*500\*2000mm



### ELK-ITM-018 Medicine Shelf

Material: Powder coated steel  
Size: 1000\*500/300\*2000mm



### ELK-ITM-019 Stainless Steel Medicine Shelf

Powder coated steel  
Size: 1200\*450\*2000mm



### ELK-ITM-020 Medicine Shelf

Powder coated steel  
Size: 1200\*450\*2000mm



## Patient Care and Medical Furniture

### ELK-ITM-021 Rotary Medicine Shelf

Material: Powder coated steel  
Size: 1200\*450\*2000mm



### ELK-ITM-022 Stainless Steel Work Table

Powder coated steel  
Size: 1800\*700\*850mm



### ELK-ITM-023 Stainless Steel Work Table

Powder coated steel  
Size: 1800\*700\*850mm



### ELK-ITM-024 Stainless Steel Work Table

Material: Powder coated steel  
Size: 1800\*800/300\*1600mm



### ELK-ITM-025 Chinese Medicine Operation Desk

Powder coated steel  
Size: 1600\*550/300\*1700/850mm



### ELK-ITM-026 Chinese Medicine Operation Desk

Powder coated steel  
Size: 1600\*550\*2000mm



### ELK-ITM-027 Stainless Steel Work Shelf

Material: Powder coated steel  
Size: 1500\*450\*2000mm



### ELK-ITM-028 Stainless Steel Work Shelf

Powder coated steel  
Size: 1500\*450\*2000mm



### ELK-ITM-029 Work Shelf

Powder coated steel  
Size: 2000\*480\*2000mm



### **ELK-CD-001** **Ordinary Parturition Bed**

#### **Main Technical Parameters**

- Back section raised from the horizontal:  $\geq 55^\circ$
- Seat section raised from the horizontal:  $\geq 20^\circ$
- Leg section lowered from the horizontal:  $\geq 85^\circ$
- Size: 1800\*500\*800mm



### **ELK-CD-002** **Ordinary Parturition Bed**

#### **Main Technical Parameters**

- Back section raised from the horizontal:  $\geq 55^\circ$
- Leg section lowered from the horizontal:  $\geq 85^\circ$
- Size: 1820\*500\*800mm



### **ELK-CD-003** **Ordinary Parturition Bed**

#### **Main Technical Parameters**

- Back section raised from the horizontal:  $\geq 55^\circ$
- Leg section lowered from the horizontal:  $\geq 85^\circ$
- Size: 1800\*500\*800mm



### **ELK-CD-004** **Ordinary Parturition Bed**

#### **Main Technical Parameters**

- Back section raised from the horizontal:  $\geq 55^\circ$
- Seat section raised from the horizontal:  $\geq 20^\circ$
- Leg section lowered from the horizontal:  $\geq 85^\circ$
- Size: 1820\*500\*800mm



## Patient Care and Medical Furniture

### ELK-CD-007A Operation Table

- The YKD007A Operation Table electric operating table is used for the operation of various parts of the human body in the operating room of the hospital.
- The table can be moved back and forth, and all of them are made of high-strength composite plates.
- X-ray diagnosis can be conducted.
- The leg plate can be turned, spread out and disassembled, and is easily adjusted.
- It is very convenient for urology surgery.
- The operating table adopts 24V DC voltage, which is safe and reliable, and can be used with a C-arm.



## Patient Care and Medical Furniture

### ELK-CD-007B Operation Table

#### Accessories:

- Main bed 1 set
- Head section 1 set
- Body support 2 set
- Shoulder support 2 set
- Anaesthetic Screen 1 set
- Knee crutch 2 set
- Armrest 2 set
- Cushion 1 set



### ELK-CD-008 Operation Table

#### Accessories:

- Anaesthetic Screen 1 set
- Shoulder Rest 2 set
- Arm support 2 set
- Body Support, Large 1 set
- Body Support, Small 1 set
- Knee Crutch 2 set
- Cushion 1 set
- Body Restraint Strap set (Optional Components)



This table are compactly designed and sturdily constructed, and can be raised or lowered hydraulically by means of an oil pump. They are suitable for surgical procedures of the head, neck, thorax, abdomen, perineum and extremities as well as for gynaecological, obstetrical, urological and oto-rhino-laryngological operations. Various positions can be quickly and safely adjusted.


The controlling devices, hand wheels and crank handles are mounted at the head-end of the table. All positions of the table are controlled entirely from the head-end. The table is double dimensions. Chassis Inside & outside cover are stainless steel.

## Patient Care and Medical Furniture

### ELK-CD-009

#### Accessories:

- Anaesthetic Screen 1 set
- Shoulder Rest 2 set:
- Arm support 2 set
- Body Support, Large 1 set
- Body Support, Small 1 set
- Knee Crutch 2 set
- Cushion 1 set
- Body Restraint Strap 1 set (Optional Components)



### ELK-CD-010

#### Accessories:


- Anaesthetic Screen 1 set
- Shoulder Rest 2 set
- Arm support 2 set
- Body Support, Large 1 set
- Body Support, Small 1 set
- Knee Crutch 2 set
- Cushion 1 set
- Body Restraint Strap 1 set (Optional Components)



## Patient Care and Medical Furniture

### ELK-CD-005A

- Stand for anaesthetic and for hanging intravenous drip bottle 1 set
- Shoulder rest and mattress 2 set
- Arm rest and mattress 2 set
- Handle 2 set
- Leg support, foot board and mattress 2 set
- Mattress on the table top 1 set
- The contaminative basin can prevent the amniotic fluid leaks, it is easy to assemble and disassemble.



## Patient Care and Medical Furniture

### ELK-CD-005B

- Stand for anaesthetic and for hanging intravenous drip bottle 1 set
- Shoulder rest and mattress 2 set
- Arm rest and mattress 2 set
- Handle 2 set
- Leg support, foot board and mattress 2 set
- Mattress on the table top 1 set
- The contaminative basin can prevent the amniotic fluid leaks, it is easy to assemble and disassemble



### ELK-CD-008

#### Accessories:

- Electric operation bed 1 set
- Shoulder support 2 set
- Anaesthetic Screen 1 set
- Handle 2 set
- Knee crutch 2 set
- Arm rest 2 set
- Cushion 1 set
- Contaminative basin 1 set

#### Range of application:

This product is applicable to antenatal care, maternal surgical childbirth, dystocia and obstetrics and gynecology and other evacuation operation.



## Patient Care and Medical Furniture

### ELK-CD-011

#### Multi-functional Electric Obstetric Bed

- External Length 2050mm
- External Width 980mm
- Table Height Min. 580mm
- Table Height Max. 950mm

#### Functions:

- Back-rest Adjustment: 0°-70°
- Trendelenburg 0°-12°
- Anti-trendelenburg 0°-5°
- Height Adjustment: 580-950mm

#### Standard Accessories:

- Pull Handle: 2pcs
- Leg Support 2pcs
- Foot Support 2pcs
- Hand Remote Controller 1pc
- Stainless steel Basin 1pc
- Instruction Manual



## Patient Care and Medical Furniture

### ELK-CL-047

#### Electric dialysis chair

##### Main technical parameters

- Block size face: 1850mm \*600mm
- Adjustment of chair height: 550 mm to 670 mm
- To break back point of view: 0°-70°
- Legs folded on board to angle: 0°-35°
- Handrail point of view before: 0°-30°
- Handrails up and down adjustment: 100mm
- Handrails plate size: 450 x 160mm



### ELK-CL-048

#### Electric dialysis chair

##### Main technical parameters

- Block size face: 1800mm \*580mm
- Block surface height: 550 mm
- To break back point of view: -10° - 70°
- Legs folded on board to angle: 0° to 85°
- Handrail point of view before: 0° to 30°
- Handrails up and down adjustment: 100mm
- Handrails plate size: 450 x 160mm



## Patient Care and Medical Furniture

### ELK-CL-049 Electric dialysis chair

#### Main technical parameters

- Block size face: 1800mm \*580mm
- Block surface height: 550 mm
- To break back point of view: -10°-70°
- Legs folded on board to angle: 0°-85°
- Handrail point of view before: 0°-30°
- Handrails up and down adjustment: 100mm
- Handrails plate size: 450 x 160mm



### ELK-CL-050 Dialysis chair

#### Main technical parameters

- Block size face: 1800mm \*580mm
- Block surface height: 550 mm
- To break back point of view: 5° - 70°
- Legs folded on board to angle: 0° to -90°
- Handrails plate size: 450 x 160mm



## Patient Care and Medical Furniture

### **ELK-CJ-001** **Stainless Steel Patient Trolley**

- Material: Stainless steel
- Stretcher and trolley can be separated
- Size: 1960\*600\*750mm



### **ELK-CJ-002** **Stainless Steel Patient Trolley**

- Material: powder coated steel
- Back rest lifting by gas spring
- Size: 1960\*600\*750mm



### **ELK-CJ-003** **Stainless Steel Patient Trolley**

- Material: stainless steel
- Stretcher and trolley can be separated
- Size: 1980\*600\*750mm



## Patient Care and Medical Furnitur

### **ELK-CJ-004** **Patient Trolley**

- Material: powder coated steel
- Back rest lifting by gas spring
- Size: 1960\* 700\*620-800mm



### **ELK-CJ-005** **Stainless Steel Patient Trolley**

- Three-crank back-rest lifting leg-rest lifting: height adjustment
- Stainless steel frame, fifth-aluminium alloy guardrail
- 5 inch silent caster, with brake
- Accessories: I.V.pole, oxygen cylinder holder, debris basket, mattress
- Size: 1930\*750\*620-880mm



## Patient Care and Medical Furniture

### ELK-CJ-007 Hydraulic Transfer Stretcher

Size: 1930\*640\*540-840mm

Crank System	1 Set
Directional Wheel	1pc
Central Locking Break Pedal	2pcs
Back Rest Gas Spring	1pc
IV Pole Previsions	2pcs
Transfer Mattress	1pc
ABS Side Rail	1set
IV Pole	1set
Directional Pedal	2pcs
Oxygen Cylinder Holder	1pc
6" Castors with Central Locking	4pcs



### ELK-CJ-008 Hydraulic Transfer Stretcher

Size: 2100-850-630-930mm

Hydraulic Pump	2pcs
74mm Wall Bumper Wheel	4pcs
Crank System	1 set
IV Pole Previsions	4pcs
Back Rest Gas Spring	2pcs
Mattress	1pc
8" Castors, Central Locking	4 pcs
IV Pole	1set
Emergency Push Handle	1set
Central Locking Break Pedal	2sets
Directional Wheel	1pc
Directional Wheel Release Pedal	2pcs
Folding Side Rail	1set
X-ray Translucent Bed Platform	1set




## Patient Care and Medical Furniture

### ELK-CJ-010 Hydraulic Transfer Stretcher

Size: 1930\*640\*630-930mm


Hydraulic Pump	2pcs
Directional wheel	1pc
Directional Pedal	2pcs
Back Rest Gas Spring	1pc
IV Pole Previsions	2pcs
Transfer Mattress	1pc
ABS Side Rail	1set
IV Pole	1set
8" Castors, Central Locking	4pcs
Central Locking Break Pedal	2 set



### ELK-CJ-011 Manual Transfer Stretcher

Size: 1930\*640\*630-930mm

Crank System	1 set
Directional Wheel	1pc
Central Locking Break Pedal	2pcs
IV Pole Previsions	2pcs
Transfer Mattress	1pc
Aluminum Alloy Side Rail	1 set
IV Pole	1set
Directional Pedal	2pcs
Oxygen Cylinder Holder	1pc
6" Castors with Central Locking	4 pcs



## Patient Care and Medical Furniture

### ELK-CJ-012

#### Hydraulic Transfer Stretcher

Size: 2100-850-630-930mm

Hydraulic Pump	2pcs
Wall Bumper Wheel	4pcs
Crank System	1 set
IV Pole Previsions	4pcs
Back Rest Gas Spring	2pcs
Mattress (optional)	1pc
6" Castors, Central Locking	4 pcs
IV Pole	1set
Emergency Push Handle	1set
Central Locking Break Pedal	2sets
Directional Wheel	1pc
Directional Wheel Release Pedal	2pcs
Folding Side Rail	1set
X-ray Translucent Bed Platform	1set




### ELK-CJ-013

#### Hydraulic Transfer Stretcher

Size: 2100\*850\*630-930mm

Hydraulic Pump	2pcs
Directional wheel	1pc
Directional Pedal	2pcs
Back Rest Gas Spring	1pc
IV Pole Previsions	2pcs
Transfer Mattress	1pc
Stainless steel Side Rail	1set
IV Pole	1set
8" Castors, Central Locking	4pcs
Central Locking Break Pedal	2 set



## Patient Care and Medical Furniture

### ELK-CJ-014

#### Patient Transportation Trolley

Size: 1920\*640\*535-860mm




Crank System	1 set
Directional Wheel	1pc
Oxygen Cylinder Holder	1pc
Back Rest Air Spring	1pc
Central Locking Break Pedal	4pcs
Plastic Side Rail	1 set
IV Pole Previsions	2pcs
Transfer Mattress	1pc
6" Castors with Central Locking	4pcs
Aluminum Alloy Side Rail	1set
IV Pole	1 set
Safe Lock for Side Rail	2 pcs



### ELK-CJ-015

#### Connecting Transfer Stretcher For Operation Room

Size: 3900\*640\*665-915mm



Crank System	2pcs
Directional Wheel Release Pedal	4pcs
Aluminum Alloy Side Rail	1set
Directional Wheel	2pcs
6" Castors with Central Locking	8pcs
IV Pole Previsions	2pcs
IV Pole	1set
Central Locking Break Pedal	4pcs



## Patient Care and Medical Furniture

### Vertical Steam Sterilizer Series



Model	Volume	Consumed power	Chamber Volume	Basket dimension	Outer dimension
DGL-35B	35L	3.5kW	φ350 x 400	φ330 x 320	498x498x940
DGL-50B	50L	3.5kW	φ350 x 525	φ330 x 460	498x498x1110
DGL-75B	75L	4.5kW	φ400 x 625	φ380 x 560	570x570x1185
DGL-100B	100L	4.5kW	φ450 x 650	φ420 x 540	588x588x1255
DGLS-35B	35L	2.5kW	φ300 x 500	φ280 x 400	400x436x840
DGLS-50B	50L	3.5kW	φ350 x 525	φ330 x 460	480x525x920
DGLS-75B	75L	4.5kW	φ400 x 625	φ380 x 560	522x567x1020
DGLS-100B	100L	4.5kW	φ450 x 650	φ420 x 540	555x599x1080
DGLS-120B	120L	4.5kW	φ450 x 760	φ420 x 660	555x599x1200
DGLS-150B	150L	6.0kW	φ500 x 780	φ470 x 720	605x649x1210
DGL-35G	35L	3.5kW	φ350 x 400	φ330 x 320	498x498x940
DGL-50G	50L	3.5kW	φ350 x 525	φ330 x 460	498x498x1110
DGL-75G	75L	4.5kW	φ400 x 625	φ380 x 560	570x570x1185
DGL-100G	100L	4.5kW	φ450 x 650	φ420 x 540	588x588x1255
DGL-35X	35L	3.5kW	φ350 x 400	φ330 x 320	698x498x940
DGL-50X	50L	3.5kW	φ350 x 525	φ330 x 460	698x498x1110
DGL-75X	75L	4.5kW	φ400 x 625	φ380 x 560	770x570x1185
DGL-100X	100L	4.5kW	φ450 x 650	φ420 x 540	788x588x1255
DGL-35M	35L	3.5kW	φ350 x 400	φ330 x 320	698x498x940
DGL-50M	50L	3.5kW	φ350 x 525	φ330 x 460	698x498x1110
DGL-75M	75L	4.5kW	φ400 x 625	φ380 x 560	770x570x1185
DGL-100M	100L	4.5kW	φ450 x 650	φ420 x 540	788x588x1255

## Patient Care and Medical Furniture

### Table Top Steam Sterilizer Series



Model	Volume	Consumed power	Chamber Volume	Outer dimension
DGT-8A	8L	1500W	φ170 x 360	445x580x400
DGT-12A	12L	1800W	φ200 x 385	445x580x400
DGT-18A	18L	2200W	φ250 x 370	480x590x450
DGT-23A	23L	2200W	φ250 x 470	480x710x450
DGT-8B	8L	1500W	φ170 x 360	445x580x400
DGT-12B	12L	1800W	φ200 x 385	445x580x400
DGT-18B	18L	2200W	φ250 x 370	480x590x450
DGT-23B	23L	2200W	φ250 x 470	480x710x450
DGT-30B	30L	3500W	φ325 x 370	710x655x575
DGT-45B	45L	3500W	φ325 x 550	890x655x575

### Portable Pressure Steam Sterilizer

Model	Volume	Consumed power	Chamber Volume
DGS-280A	18L	2200W	φ280x295
DGS-280A+	24L	2200W	φ280x390
DGS-280B	18L	2200W	φ280x295
DGS-280B+	24L	2200W	φ280x390
DGS-280C	18L	2200W	φ280x295
DGS-280C+	24L	2200W	φ280x390





 **Elektro**  
GENESIS

Cherish your life, Cherish your health!

